


The Invisible Akin  
I N V I S I B L E O N X - R A Y

PEEK COMPRESSION STAPLE

# OS2<sup>®</sup>-vp

Low Profile Design



Straight 90°

Angled 26°

9mm Compression Staples  
Compression on Release Design  
Teeth Ensure Stability

## PEEK COMPRESSION STAPLE


# OS2<sup>®</sup>-vp

Low Profile Design / 90° Straight & 26° Angled

- Designed for Akin type osteotomies with compression on release
- PEEK radiolucency and biocompatibility
- No nickel or metal sensitivity concerns
- 26° Angled option designed to match the metaphysis of the phalanx

This staple is **invisible on x-ray** (radiolucent). It is manufactured from high strength, medical-grade polyetheretherketone (PEEK). PEEK is an inert polymer that exhibits excellent mechanical and thermal properties. These properties make PEEK the most suitable material for an invisible Akin staple. Further benefits of PEEK are as follows:


- Radiolucent for easy, artifact-free monitoring and assessment of the healing site.
- Bone-like modulus



- A** COMPRESSION - The compression on release design legs of the Staple converge 20°
- B** PRIMARY STABILITY - The compression Staple's teeth design ensure primary stability on insertion

### STAPLE SETUP


After opening the sterile pack, release the Loader from the Holder/Drill Guide. Free the K-wires from the Impactor/Tamp. The OS2<sup>®</sup>-vp Staple is nested within the Holder.



### ORDERING INFO

CAT NO.	SIZE [mm]	STYLE
T40 SP190	OS2 <sup>®</sup> -vp - Staple Kit	90° Interaxis 9mm
T40 SP126	OS2 <sup>®</sup> -vp - Staple Kit	26° Interaxis 9mm

### TECHNIQUE FOR AKIN OSTEOTOMY



#### Bone Prep

**1** Position the Drill-guide end of the combo Inserter on the bone bridging the osteotomy. Use the 1.6 x 70mm K-wire to drill the first hole. Leave the K-wire embedded and drill the second hole with the 1.6 x 100mm K-wire. Remove the Drill Guide / Inserter and both wires and prepare to insert the Staple.

#### Load and Insert the Staple

**2** Using the Inserter end of the Drill Guide combo, carefully **load**

**the Staple** by sliding it into the Holder well at 90°. A firm quarter turn clockwise should engage the Staple, spreading the legs and loading the compression Staple. **3** Insert the Staple into the prepared holes. **4** When the Inserter / Staple combo is flush to the bone, release the Staple with a quarter turn counterclockwise. When released from the Inserter, the Staple legs converge to provide compression across the fusion site. **5** Use the Impactor to seat the Staple flush to the bone as desired.

All content contained herein is furnished for informational purposes only. In2Bones does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgment and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and precautions are listed in the implant package insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure. Availability of these products might vary from a given country or region to another as a result of specific local regulatory approval or clearance requirements for sale in such country or region. **CAUTION:** Federal law (USA) restricts this device to sale and use by, or on the order of a physician.



**Corporate Headquarters**  
In2Bones Global, Inc. • Memphis, TN • USA  
844. 602. 6637 • Info@i2b-USA.com

**International Office**  
In2Bones SAS • Lyon • France  
+33 (0)4 72 29 26 26

In2Bones.com



OS2, the In2Bones name and logo are trademarks of In2Bones or its affiliates. In2Bones USA Memphis TN, 38119 USA / In2Bones SAS, 69130 Ecully, France - © 2019 In2Bones USA, Memphis, TN • Patent Pending • All rights reserved • OS2VP - ST - D, 0420