

THE REVO® / MINI-REVO® SHOULDER FIXATION SYSTEM

SURGICAL TECHNIQUE

THE SUTURE FRIENDLY, EASY TO USE, EASY TO RE-THREAD, PRE-LOADED

BIG EYED REVO

RETRIEVABLE ANCHOR




CONMED
LINVATEC

REVO®




The following techniques are described by Stephen J. Snyder, M.D., Van Nuys, CA.


Arthroscopic repair of the rotator cuff is now a well accepted technique for the experienced arthroscopic surgeon. This has been made possible by technological advancements in the field, specifically the specialized rotator cuff **Big-Eyed Revo** Anchor and the new suturing devices that are currently




The rotator cuff tear is carefully evaluated with the arthroscope on both the articular and bursal sides. The frayed edges of the cuff are debrided. The best view of the rotator cuff is usually "the 50 yard line" view with the arthroscope in a lateral subacromial portal which is located at the center point of the rotator cuff tear.



A crescent suture hook with a **Shuttle Relay™** or a **Blitz®** suture passer is used to repair any side-to-side tears in the rotator cuff tendon.



The **Revo** bone punch is inserted through a small puncture wound adjacent to the lateral border of the acromion to create pilot holes for the anchors. The punch is directed to enter the bone in a medial direction below the subchondral bone at approximately a 45-degree angle. The pilot holes are angled away from the center of the trough in a fan-like pattern.




The 4mm **Big Eye Revo** anchor, loaded with either one or two strands of #2 braided suture, is inserted directly through the percutaneous puncture wound (**no cannula is needed to insert this anchor**).

ROTATOR CUFF REPAIR SURGICAL TECHNIQUE


Steps in repair technique

available. The surgeon must have an excellent understanding of the technique and must practice the suture passing and knot tying steps before attempting the operation. The following illustrations highlight the important steps in a typical rotator cuff repair. ConMed Linvatec will be happy to provide


more comprehensive videotape instructions. You may also use an "Alex The Shoulder Professor"™* shoulder model to practice these techniques prior to surgery. Information can be obtained by calling your local ConMed Linvatec representative or Customer Service at (800) 925-4255.




3
After passing the curved suture hook across the tear, a #2 nonabsorbable braided suture (or other chosen suture) is carried across the tear and the suture tails tied together.



4
Light decortication or a superficial trough is created at the anatomical neck of the humerus, adjacent to the articular cartilage, using a high speed burr or shaver.





7
The **Revo**® anchor is inserted into the bone until the seating ring on the driver is just below the level of the cortical bone. The vertical seating mark is aligned perpendicular to the cuff edge.





8
The anchor security is tested by pulling on the suture strand. The eyelet of the anchor is best aligned with the opening perpendicular to the cuff to avoid having a twist (in the sutures after passing thru the cuff).

REVO




9
The arthroscope can be positioned in the lateral, anterior, or posterior portals, whichever affords the best view of the torn cuff and suturing materials. The anterior anchor is usually inserted first.

10
One limb of the suture is carried out through a posterior operating cannula (**ConMed Linvatec** 6mm operating cannula) using a crochet hook.




12a
ALTERNATIVE METHOD: With the scope viewing from the posterior portal, a modified Caspari suture punch can be inserted through an 8.4 mm cannula in the lateral portal to pass a **Shuttle Relay**[™] suture passer from bottom to top through the cuff. The suture passer is carried out the anterior cannula.


12b
MODIFIED CASPARI SUTURE PUNCH CONTINUED: The eyelet of the **Shuttle Relay** suture passer is loaded with the suture outside the lateral cannula and carried through the cuff from bottom to top by pulling on the opposite end. The surgeon may choose to use a simple suture by passing one of the strands through the tendon, or a mattress suture by passing both suture strands through the tendon. For additional security, two suture strands can be used in the **Big Eye Revo** anchor and one strand of each suture passed through the tendon at a position approximately 45 degrees away from the anchor to create a very strong, double-suture anchor point.




The **Spectrum**[®] suture hook is passed through the torn rotator cuff from top to bottom. If a crescent type hook is used, it may be inserted through the cannula. If a 90-degree hook is used, a larger cannula may be required to accommodate the larger size. The **Shuttle Relay**[™] suture passer is passed through the hook and retrieved with a grasping clamp through the alternative (posterior) cannula.



The suture strand is loaded into the eyelet of the shuttle relay suture passer outside the posterior cannula and carried through the cuff from the articular to the bursal side by retracting the opposite end of the suture passer out the anterior cannula.



At the completion of the repair, the torn end of the rotator cuff is tightly compressed to the bone to promote strong rotator cuff tendon healing.




This illustration of the final repair shows a mattress suture anterior, a suture in the middle and a double suture posterior. Notice the final side-to-side repair.


mini-REVO®




Anterior dislocation of the shoulder is a very common orthopedic injury, mostly seen in young athletes, that often includes a labral avulsion associated with a capsuloligamentous elongation, or stretch injury. Arthroscopic repair of this injury can now be easily performed with minimal morbidity on an outpatient basis using **Mini-Revo** suture anchors, strong nonabsorbable sutures, and secure arthroscopic knots. The surgeon must decide in each individual case how much capsular plication surgery is necessary with the Bankart repair. Using the described capsular tuck procedure, appropriate ligamentous tension can be restored while repairing the labrum to its anatomical position on the edge of the articular cartilage. Although each case requires individual assessment, the steps described here represent the basic repair procedure. The interested arthroscopic surgeon should review the videotape and rehearse the steps preoperatively using "Alex the Shoulder Professor"™.




With the arthroscope inserted in the anterior superior portal, the Bankart lesion is evaluated and debrided with the shaver and the liberator elevator is used to completely mobilize the capsule and labrum from the neck of the glenoid.




A **Mini-Revo** suture anchor loaded with a #2 braided nonabsorbable suture is inserted into the inferior-most hole and seated until the black seating ring is just below the surface and the vertical alignment mark is facing the articular surface so that the eyelet opens towards the labrum.




2 A 4.5 mm spherical Linvatec bur lightly decorticates the anterior glenoid beginning at the edge of the articular cartilage.



3 A **Mini-Revo**®* punch or drill bit is used to create three pilot holes at the 4:30, 4:00 and 3:30 positions on the anterior inferior glenoid. These holes are located **on** the articular surface just above the corner edge of the anterior inferior glenoid. (N.B. In hard bone a **Mini-Revo Tap** may be used after creating the pilot hole with a punch or drill to facilitate insertion of the **Mini-Revo** anchor.)




5 The lower limb of the suture is retrieved with a crochet hook out through the posterior cannula.




6 A curved suture hook is passed in a "pinch-tuck-stitch fashion" through the capsule and **under** the torn labrum inferior to the level of the anchor. A **Shuttle Relay**™ suture passer is passed through the needle and retrieved out through the posterior cannula with a grasping clamp.

mini-REVO[®]




7
The suture is loaded in the eyelet of the suture passer, which is then carried down below the labrum and through the pinch-tuck-stitch of the capsule.




8
By tightening the two sutures, the capsule is folded up and on to the labrum.


Using The BLITZ[®] Suture Retriever




1
Advance the curved tip of the **Blitz**[®] Suture Retriever through the soft tissue.




9 An arthroscopic Tennessee Slider or **Revo**® knot is tied through the midglenoid cannula, using the knot pusher.



10 Additional anchors are inserted and sutures passed as needed to completely restore tension in the capsule and reconstruct the anterior labral buttress.



2 Deploy the wire loop of the **Blitz**® suture retriever. Then, insert a crochet hook through the wire loop, and retract the suture strand.



3 With the wire loop almost fully retracted, pull the **Blitz** suture retriever back through the soft tissue and out the cannula, carrying the suture strand through the tissue. The suture can now be tied and secured.

Arthroscopic REVO KNOT

The arthroscopic Revo® knot is an extremely important knot for all surgeons performing advanced shoulder reconstruction procedures. This knot can be used in any and all situations, whether or not the suture material slides freely through the tissue and anchor eyelet. If a complex stitch such as a figure-of-eight or double-pass mattress stitch is used, this knot is preferable to any sliding knot. In addition, when capsular plication is performed, it is important not to use a sliding knot because of the possible trauma to the labrum as the suture is pulled through.




❶ Both suture tails are the same length and the loop-handled knot pusher is threaded onto the suture which has been passed through the soft tissue. This original "post" is positioned on the left side. The knot pusher is passed down the original post suture to ensure that there are no twists or soft tissue obstructions.



❷ An **underhand** half-hitch is placed around the original post and advanced into position on the edge of the soft tissue.

Tennessee SLIDER KNOT

If the suture slides easily through the soft tissue and the anchor eyelet, a Tennessee slider knot will be a useful alternative to the "Revo knot." The Tennessee slider cannot be used if the sutures do not easily slide through the soft tissues. If there is any doubt about the freedom of suture passage, then the Revo knot should be used.



❶ The knot pusher is threaded on the first post strand (left side) and a clamp is placed to secure it. The knot pusher is passed into the joint to ensure that there are no twists or obstructing soft tissue. The suture is then arranged so that the original post suture is short, with only 5cm of suture material outside the cannula.



❷ A "buntline hitch" sliding knot is tied with the loop strand around the original post.



3 Tension is held on the post suture while a second **underhand** half-hitch is worked down the post suture to reinforce the first hitch.



4 An **overhand** half-hitch is next placed on the same initial post and worked down into position on the other two throws.



5 The knot pusher and clamp are changed to the opposite suture and after checking for twists and soft tissue, an **underhanded** throw is advanced down onto the knot stack.




6 The knot pusher is advanced to "past point" to lock the half-hitch securely.



7 A fifth **overhand** half-hitch is placed over the second post and worked down into position on the knot stack.



8 Sometimes a sixth half-hitch can be used as the surgeon prefers, and the suture tails are cut with microscissors.



9 By pulling firmly on the post, the suture slides through the tissues and anchor eyelet, and the knot is carried down to the soft tissues. The hitch is locked in position by pulling on the non-post strand.




10 An additional **underhand** half-hitch knot is then placed on the original post and carried down into position on the soft tissue to add additional friction to the sliding knot.



11 The knot pusher and clamp are changed to the opposite post. After clearing the suture, an **underhand** half-hitch is worked into the joint.



12 A final **overhand** half-hitch is placed on the same post and worked down into position on the knot stack. Past-pointing is used to permanently lock the half-hitches in place.



13 An additional **underhand** hitch can be used for further security if the physician chooses. The suture tails are cut short.

ORDERING INFORMATION

THE REVO®/ MINI-REVO® SHOULDER FIXATION SYSTEM

Implants

Revo Suture Anchor, 4mm O.D. (1.8mm I.D.) x 12mm (pre-loaded with #2 braided, polyester suture)	C6101A
Mini-Revo Suture Anchor, 2.7mm O.D. (1.5mm I.D.) x 8.5mm (pre-loaded with #2 braided, polyester suture)	C6109A

Revo Instrument Set

Revo Driver	C6102B
Suture Threader	C6113
Revo Drill Bit, 2mm Diameter	C6106A
Revo Bone Punch, 2mm Diameter	C6107
Revo Remover Guide	C6126
Revo Remover	C6127
Instrument Tray with Lid	C6108

Mini-Revo Instrument Set

Mini-Revo Driver	C6110A
Suture Threader	C6113
Mini-Revo Bone Punch, 1.5mm Diameter	C6114
Mini-Revo Drill Bit, 1.5mm Diameter	C6115A
Mini-Revo Drill Guide	C6116
Mini-Revo Remover Guide	C6124
Mini-Revo Remover	C6125
Mini-Revo Tap	C6128
Instrument Tray with Lid	C6117

Suture Passing Instruments

Slotted Jaw, Suture Punch, 4mm needle	18.1008
---------------------------------------	---------

Spectrum® Instrument Set

Suture Hook Handle	27.00011
Suture Hook, Straight	97.10015
Suture Hook, 45° Left Curve	97.14115
Suture Hook, 45° Right Curve	97.14215
*Suture Hook, 90° Left Curve	97.19115
*Suture Hook, 90° Right Curve	97.19215
Crescent Suture Hook, Small Curve, 3 x 15mm	C8740
Crescent Suture Hook, Medium Curve, 4 x 20mm	C8741
Crescent Suture Hook, Large Curve, 6x25mm	C8742
Shuttle-Relay™ Suture Passer (10 per box)	C6004
Blitz® Suture Retriever, Straight (6 per box)	C6111
Blitz Suture Retriever, 45° Left (6 per box)	C6211
Blitz Suture Retriever, 45° Right (6 per box)	C6311
Hawkeye® Suture Needle (6 per box)	C6001

Accessories

Loop Handle Knot Pusher	C6112
Crochet Hook	C6105
Microscissors, 2.75mm Diameter, Straight	2.10011
Suture Retrieval Forceps, 3.4mm Diameter	16.1018
Grasping Forceps, 3.4mm Diameter, Straight with ratchet	11.1001
Liberator™ Knife	25.50014
Rasp Liberator Knife	25.50015

Surgical Procedure Videos

The Revo Rotator Cuff Fixation System	VT65
The Mini-Revo Labral Repair System	VT73
Using Mini-Revo Bankart, Inferior Capsular Shift, Mini-Open, Superior Labrum Repairs	VT89



Revo Instrument Set



Mini-Revo Instrument Set

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY

‡Build to order, not CE Marked "Alex The Shoulder Professor" is a trademark of Pacific Research Laboratories.



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
 www.linvatec.com