

Clinical Summary

Early Experience in Collateral Ligament Repair of the Hand with the Use of a Novel Bio-Composite Scaffold: Technique and Outcomes

Authors: Devin C Potter, Zachary Hanson, M.D., Gary M Lourie, M.D.

Journal: Journal of Orthopaedic Experience & Innovation



Aim:

To evaluate the early clinical and functional outcomes of a small cohort using a novel bio-composite scaffold (BioBrace®) in the surgical repair of hand collateral ligament injuries, particularly in high-level athletes.

Diagnosis & Clinical History

Condition: Acute or chronic collateral ligament injuries of the hand, especially the thumb UCL and RCL.

Population: 13 patients (mean age 21), 85% of whom were high-level athletes.

Symptoms: Pain, instability, reduced pinch strength, and limited range of motion.

Surgical Technique

- **Technique:**
 - » Surgical repair or reconstruction using BioBrace®, a bio-inductive, bioresorbable scaffold
 - » Anchored with all-suture anchors to minimize bone damage
 - » Post-op care included splinting and early mobilization

Results

- **QuickDASH:** Improved from 81.5 pre-op to 4.1 at 6 weeks.
- **VAS Pain:** Decreased from 4.3 to 0.6.
- **ROM:** Non-limiting by 4–6 weeks.
- **Pinch Strength:** >85% of contralateral side.
- **Return to Sport:** 82% returned by 8 weeks; all by 12 weeks.
- **Complications:** No infections, osteolysis, or re-injuries. Two cases of mild erythema resolved with splint adjustment.

Key Takeaways

- BioBrace® offers a promising alternative to non-absorbable suture augmentation by:
 - » Reducing stress shielding
 - » Avoiding complications from large polymer anchors
 - » Supporting early mobilization and return to sport
- Early outcomes show excellent functional recovery, minimal complications, and high patient satisfaction.