



SEE-THROUGH HUB WITH DYNAMIC, TRANSVERSE COMPRESSION

Lapidus Plates

CoLink® **View** Plating System



Low Profile, Anatomic Design, Type II Anodized



Mechanical Compression Designed to Stimulate the Fusion Process



A GLOBAL EXTREMITY COMPANY

Lapidus Plates

CoLink® **View** Plating System




Dynamic, Transverse Compression

Use of the Transverse, Cross-joint Screw provides mechanical compression across the fusion site and to stimulate the fusion process.



Lapidus Std. and +1 Plates

Surgical Technique



PROVISIONAL PLACEMENT & TRIAL PLATE EVALUATION

After the appropriate incisions, reduce the fragments and temporarily fix using cross-joining fixation pins. Position plate trial to confirm placement.

NOTE: Choose the correct plate. (Std, 1mm Step).







PLATE PREP AND POSITIONING

Once the correct Plate has been determined, open the sterile package to retrieve the Sterile Plate. If necessary bend the Plate to the required shape using the Plate Benders provided within the Instrument Set. Do not bend the Plate across any Screw holes. Plates should only be bent in one direction. Never re-bend Plates. Temporarily fix in place with the Olive Wires in the Wire Slots. Position the distal Olive Wire in the most proximal placement in the Wire Slot.






SCREW PREP

Begin Screw placements with the most distal hole and follow the suggested sequence at right. NOTE: All Plate Screw holes can accommodate both locking and non-locking screws (3.0mm and 3.5mm diameter). The Transverse Hole can only accommodate a 3.5mm Transverse Lag Screw.



- 3.0/3.5mm Locking/
Non-locking option ○
- 3.5mm Transverse
ONLY ○

SCREW SIZE & SEQUENCE





TRANSVERSE SCREW PREP

Once the most proximal screws are in place, the Transverse Drill Guide can be used to prepare the Transverse screw hole. For optimum lag screw orientation, ensure the Drill Guide is seated completely in the Transverse hole. Correct placement is with Handle vertical (90°) to plate with Screw trajectory at ~40° to plate.



MEASURE TRANSVERSE SCREW LENGTH

Use the Laser Markings on the provided Reamer with the provided Drill Guide for correct Transverse Screw length. OPTIONAL: The provided Depth Gauge can be used, however the depth reading must be adjusted -2mm to compensate for the space between Depth Gauge on Plate surface and actual bone surface. NOTE: Also, a Screw length subtraction of 1-2mm may be required if the joint is not completely reduced when the measurement is performed.



TRANSVERSE SCREW INSERTION / COMPRESSION

The lag screw should be tightened in a clock-wise motion. Once the joint is compressed, the remaining proximal screws are inserted and all temporary fixation wires are removed.

SEE-THROUGH HUB WITH DYNAMIC, TRANSVERSE COMPRESSION

Lapidus Plates

CoLink® View Plating System



CoLink® View Lapidus Plates Options:

Packaged with PEEK Insert:

CATALOG NO	DESCRIPTION
P40 ST141 ...	Lapidus Plate, Std., Right
P40 ST241 ...	Lapidus Plate, Std., Left
P40 ST151 ...	Lapidus Plate, +1 mm, Right
P40 ST251 ...	Lapidus Plate, +1 mm, Left

Also available



CoLink® View MTP Plates:

Packaged with PEEK Insert:

P40 ST131 ...	MTP Plate, 5-Hole, Right
P40 ST231 ...	MTP Plate, 5-Hole, Left
P40 ST136 ...	MTP Plate, 6-Hole, Right
P40 ST236 ...	MTP Plate, 6-Hole, Left




Titanium Compression Slot

CoLink® View, Titanium Insert

P40 ST614 ... MTP / Lapidus Compression Slot


P04 S0001.....CoLink® Disposable Sterile Instruments

For 3.0/3.5 Screws; Drills, Olive Wires, Guide Pins





CoLink® Plate Screws Non-Locking

CATALOG NO	DIA x LENGTH, STYLE
V30 ST208 ...	3.0 x 8mm, Non-Locking
V30 ST210 ...	3.0 x 10mm, Non-Locking
V30 ST212 ...	3.0 x 12mm, Non-Locking
V30 ST214 ...	3.0 x 14mm, Non-Locking
V30 ST216 ...	3.0 x 16mm, Non-Locking
V30 ST218 ...	3.0 x 18mm, Non-Locking
V30 ST220 ...	3.0 x 20mm, Non-Locking
V30 ST222 ...	3.0 x 22mm, Non-Locking
V30 ST224 ...	3.0 x 24mm, Non-Locking
V30 ST226 ...	3.0 x 26mm, Non-Locking
V30 ST228 ...	3.0 x 28mm, Non-Locking
V30 ST230 ...	3.0 x 30mm, Non-Locking



CoLink® Plate Screws Locking

CATALOG NO	DIA x LENGTH, STYLE
V30 ST308 ...	3.0 x 8mm, Locking
V30 ST310 ...	3.0 x 10mm, Locking
V30 ST312 ...	3.0 x 12mm, Locking
V30 ST314 ...	3.0 x 14mm, Locking
V30 ST316 ...	3.0 x 16mm, Locking
V30 ST318 ...	3.0 x 18mm, Locking
V30 ST320 ...	3.0 x 20mm, Locking
V30 ST322 ...	3.0 x 22mm, Locking
V30 ST324 ...	3.0 x 24mm, Locking
V30 ST326 ...	3.0 x 26mm, Locking
V30 ST328 ...	3.0 x 28mm, Locking
V30 ST330 ...	3.0 x 30mm, Locking



V35 ST208 ...	3.5 x 8mm, Non-Locking
V35 ST210 ...	3.5 x 10mm, Non-Locking
V35 ST212 ...	3.5 x 12mm, Non-Locking
V35 ST214 ...	3.5 x 14mm, Non-Locking
V35 ST216 ...	3.5 x 16mm, Non-Locking
V35 ST218 ...	3.5 x 18mm, Non-Locking
V35 ST220 ...	3.5 x 20mm, Non-Locking
V35 ST222 ...	3.5 x 22mm, Non-Locking
V35 ST224 ...	3.5 x 24mm, Non-Locking
V35 ST226 ...	3.5 x 26mm, Non-Locking
V35 ST228 ...	3.5 x 28mm, Non-Locking
V35 ST230 ...	3.5 x 30mm, Non-Locking
V35 ST232 ...	3.5 x 32mm, Non-Locking
V35 ST234 ...	3.5 x 34mm, Non-Locking
V35 ST236 ...	3.5 x 36mm, Non-Locking
V35 ST238 ...	3.5 x 38mm, Non-Locking
V35 ST240 ...	3.5 x 40mm, Non-Locking



V35 ST308 ...	3.5 x 8mm, Locking
V35 ST310 ...	3.5 x 10mm, Locking
V35 ST312 ...	3.5 x 12mm, Locking
V35 ST314 ...	3.5 x 14mm, Locking
V35 ST316 ...	3.5 x 16mm, Locking
V35 ST318 ...	3.5 x 18mm, Locking
V35 ST320 ...	3.5 x 20mm, Locking
V35 ST322 ...	3.5 x 22mm, Locking
V35 ST324 ...	3.5 x 24mm, Locking
V35 ST326 ...	3.5 x 26mm, Locking
V35 ST328 ...	3.5 x 28mm, Locking
V35 ST330 ...	3.5 x 30mm, Locking
V35 ST332 ...	3.5 x 32mm, Locking
V35 ST334 ...	3.5 x 34mm, Locking
V35 ST336 ...	3.5 x 36mm, Locking
V35 ST338 ...	3.5 x 38mm, Locking
V35 ST340 ...	3.5 x 40mm, Locking




CoLink® View 3.5mm Transverse Screws

CATALOG NO	DIA x LENGTH, STYLE
V35 ST418	3.5 x 18mm, Transverse
V35 ST420	3.5 x 20mm, Transverse
V35 ST422	3.5 x 22mm, Transverse
V35 ST424	3.5 x 24mm, Transverse
V35 ST426	3.5 x 26mm, Transverse
V35 ST428	3.5 x 28mm, Transverse

CoLink® View 3.5mm Transverse Screws

CATALOG NO	DIA x LENGTH, STYLE
V35 ST430	3.5 x 30mm, Transverse
V35 ST432	3.5 x 32mm, Transverse
V35 ST434	3.5 x 34mm, Transverse
V35 ST436	3.5 x 36mm, Transverse
V35 ST438	3.5 x 38mm, Transverse
V35 ST440	3.5 x 40mm, Transverse



INSTRUMENT TRAY

- 1 Driver Handle
- 2 Color Coded Non-locking Drill Guides
- 3 Transverse Drill Guide
- 4 Color Coded Locking Drill Guides
- 5 Insert Tool
- 6 Depth Gauge
- 7 Plate Benders
- 8 Cup & Cone Reamers
- 9 T8 Driver
- 10 Plate Trials

STERILE, OR READY

CoLink® Implants (Plates and Screws) are individually packaged in sterile boxes or tubes for OR efficiency and economy.

The color-coded package labels identify the implant type, style, configuration and length



All content contained herein is furnished for informational purposes only. In2Bones does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgment and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and precautions are listed in the implant package insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure. Availability of these products might vary from a given country or region to another as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.



Corporate Headquarters

In2Bones Global, Inc. • Memphis, TN • USA
844. 602. 6637 • Info@i2b-USA.com

International Office

In2Bones SAS • Lyon • France
+33 (0)4 72 29 26 26

FOLLOW US....

