## Hercules® XE Extended Eyelet Knotless Anchor for BioBrace®

Extended eyelet specifically designed for the 5mm x 250mm BioBrace® Reinforced Implant (sold separately) SURGICAL TECHNIQUE



# Hercules® XE

Extended Eyelet Knotless Anchor for BioBrace®

### **Product Information**

### **Product Description**

The Hercules<sup>®</sup> XE Knotless Anchor System provides a radiolucent PEEK implant for use in a range of soft tissue fixation applications. The extended eyelet (XE) allows for use with BioBrace<sup>®</sup> 5 x 250mm, which is a biocomposite soft tissue scaffold intended for use in a range of surgical procedures to reinforce soft tissue where weakness exists.


The BioBrace® Reinforced Implant is sold separately.

**Important:** BioBrace<sup>®</sup> is not intended to be a standalone construct, but an augment to other tissue/ligament repair procedures.

### Indications

The Hercules<sup>®</sup> Knotless XE Anchors with BioBrace<sup>®</sup> are intended for use in the following applications for:

- **Shoulder:** Rotator Cuff, Bankart and SLAP lesion repair, Biceps tenodesis, Acromio-Clavicular separation and Deltoid Repair, Capsular shift and Capsulolabral reconstruction
- **Foot and Ankle:** Lateral and medial ankle instability, Achilles tendon and Metatarsal ligament repair, Hallux Valgus and Midfoot reconstruction
- **Knee:** Medial collateral and Lateral collateral ligament repair, Patellar tendon, and Posterior oblique ligament repair
- **Hand/Wrist:** Scapholunate ligament, Radial collateral ligament and Ulnar collateral ligament reconstruction
- **Elbow:** Biceps tendon reattachment, Tennis elbow repair, Ulnar and Radial collateral ligament reconstruction.



### **Contraindications**

General contraindications for the use of these implants for joint reconstruction include:


- · Significant bone demineralization
- · Inadequate neurovascular status
- · Poor blood supply, which may impair the healing process
- · Inadequate skin or musculotendinous system
- · Physiologically or psychologically unsuitable patient
- · Surgical procedures other than those indicated
- · Possibility for conservative treatment
- Bone, musculature, tendons, or adjacent soft tissue compromised by disease, infection, or prior implantation, which cannot provide adequate support or fixation for the prosthesis.
- · Known allergies
- · Diabetes
- · Active infection to materials
- · Patients with high level of activity
- Pathological conditions of the soft tissue which would compromise secure fixation

## Hercules<sup>®</sup> XE Extended Eyelet Knotless Anchor for BioBrace<sup>®</sup> Surgical Technique | ATFL Repair



### Site Preparation

Expose and prepare the site following standard procedures (**Figure 1**), using available equipment and instrumentation, per surgeon discretion and expertise. Implantation of the Hercules® XE Knotless Anchors with BioBrace® should be performed as an augment to tissue/ligament repair procedures. Identify good quality bone for effective Anchor placement and determine appropriate Anchor size and corresponding Instrument Kit for the indicated area.



#### **Insert Guide Wires**

Insert the first Guide Wire ( $\emptyset$ .062mm x 9" length) perpendicular to the bone surface into the talus (**A**) at the insertion point of the anterior talofibular ligament (ATFL). Then insert the second Guide Wire into the fibula (**B**) at the origin of the ATFL (**Figure 2**).



### **Create Pilot Holes**


Starting with the talus, place the Cannulated Drill and Drill Guide over the Guide Wire (**Figure 3**). Advance the Drill until the drill stop reaches the Guide, then remove the Drill and Guide Wire.

Repeat for the fibula.

**Note:** The Guide Wires may be left in the site for reference until ready to insert the Anchors, so that the pilot holes are not mislocated.

Where softer bone quality is present, the Drill contained in the Hercules<sup>®</sup> 4.5mm Fully Threaded Anchor Instrument Kit (A02 S0004) can be used as the primary drill for the 3.5mm Hercules<sup>®</sup> XE Knotless Anchor.

## Hercules<sup>®</sup> XE Extended Eyelet Knotless Anchor for BioBrace<sup>®</sup> Surgical Technique | ATFL Repair



### Thread BioBrace<sup>®</sup> into the Anchor

Using the Suture Threader attached to the Driver assembly, pass the BioBrace<sup>®</sup> Reinforced Implant through the extended eyelet of the Hercules<sup>®</sup> Knotless XE Anchor and pull the Suture Threader through the eyelet until one end of BioBrace<sup>®</sup> is made available opposite the point of entry (**Figure 4**).

#### Alternative Method:

If needed, cut the end(s) of BioBrace<sup>®</sup> to a shallow taper. With the Anchor loaded onto the driver, thread the tapered end of BioBrace<sup>®</sup> manually through the eyelet of the Anchor and grab with fingertips or forceps on the opposite side of the eyelet and pull through (**Figure 5**).





### **Talus Anchor Insertion**

Ensure that only a portion of the BioBrace<sup>®</sup> material is pulled through the eyelet of the initial Anchor for the Talus, allowing the majority of the material to remain available for placement in the fibula. Using the Driver assembly, insert the Knotless Anchor into the talus through the prepared hole and advance the Anchor with a mallet until fully seated and flush to the bone surface. Pull the Driver away from the Anchor to disengage (**Figure 6**).

Note: Do not apply clockwise quarter turn for use with BioBrace®.




### **Fibula Anchor Insertion**

Using the longer, loose tail, thread a second Anchor with the BioBrace<sup>®</sup> and manually tension the BioBrace<sup>®</sup>, via surgeon-preferred method, prior to final seating of the second Anchor (**Figure 7**). Trim away any excess BioBrace<sup>®</sup> (**Figure 8**).

### **Close the Site**

Close the site via surgeon-preferred method.



## Hercules<sup>®</sup> XE Extended Eyelet Knotless Anchor for BioBrace®

#### Hercules® XE Knotless Anchor Implant Kits

ON	CATALOG NO.
ıles® XE Knotless Anchor, 3.5mm - Sterile	A40 SP135
ıles® XE Knotless Anchor, 4.5mm - Sterile	A40 SP145
ıles <sup>®</sup> Knotless Anchor, 5.5mm - Sterile	A40 SP055
	امم ules® XE Knotless Anchor, 3.5mm - Sterile ules® XE Knotless Anchor, 4.5mm - Sterile ules® Knotless Anchor, 5.5mm - Sterile

### Hercules<sup>®</sup> XE Knotless Anchor Instrument Kits

DESCRIPTION	CATALOG NO.
Hercules® XE Knotless Instrument Kit, Size 3.5mm	A04 S0135
Hercules <sup>®</sup> XE Knotless Instrument Kit, Size 4.5mm, 4.7mm Drill	A04 S0147
Hercules <sup>®</sup> Fully Threaded Anchor Instrument Kit, 4.5mm - Sterile	A02 S0004
Hercules <sup>®</sup> Knotless Anchor Instrument Kit, 5.5mm - Sterile	A04 S0055

### Guide Wire Kit

BioBrace<sup>®</sup> FA 40mm x 60mm

DESCRIPTION CATALOG NO. Guide Wires (2), S. Trocar, Smooth, .062" x 9", Sterile, Single-Use P06 S2301

5mm x 250mm





23mm x 30mm



BioBrace® FA Reinforced Implants (Sold Separately)	
DESCRIPTION	CATALOG NO.
BioBrace <sup>®</sup> FA 5mm x 250mm	BB5X250-FA
BioBrace <sup>®</sup> FA 23mm x 30mm	BB23X30-FA



In2Bones USA

844 602 6637

In2Bones SAS 6000 Poplar Ave., Suite 115 28, chemin du Petit Bois 69130 Ecully Memphis, TN 38119 France www.conmed.com/in2bones +33 (0)4 72 29 26 26

CONMED Corporation 11311 Concept Blvd. Largo, FL 33773 USA

™ and ® denote Trademarks and Registered Trademarks of CONMED or its affiliates. ©2024 CONMED Corporation. All Rights Reserved.

Some products not available in all markets. This document provides information regarding medical devices and instruments in surgical procedures. It is not medical advice and healthcare professionals should use their own judgment before using to treat a particular patient. Healthcare professionals should be trained in the use of such devices before surgery and should always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any product. BB40X60-FA