

Shoulder Restoration System™

PopLok® Knotless Suture Anchor

Simple, Secure, Versatile – all-PEEK knotless anchor system for rotator cuff and instability repairs









Arthroscopic Repairs using the PopLok® Knotless Suture Anchor

The all-PEEK PopLok® knotless suture anchor provides easy and reliable fixation for the repair of both labral and rotator cuff pathology. The characteristic (audible) "POP" indicates that the suture is securely locked between the two sliding components, independent of the bone-anchor interface, and deployable wings ensure repair security. The unique ability to modify the amount of tension applied to the soft tissues once the anchor is inserted lessens the risk of compromising the vascular supply, a key feature for healing. Furthermore, the knotless feature coupled with a simple insertion technique can lead to improved surgical efficiency.

Technique Contributions by John Randle, MD, Newmarket, Ontario, Canada See animation and surgical videos at srs.linvatec.com

INSTABILITY SURGICAL TECHNIQUES

Arthroscopic Bankart Repair using the 2.8 or 3.3mm PopLok® Knotless Suture Anchor









Pass Suture

A free strand of Hi-Fi® suture is passed through the capsule and full thickness of labrum using the desired $Spectrum^{\text{\tiny{\$}}}\ hook\ or\ Spectrum^{\text{\tiny{\$}}}\ MVP^{^{\text{\tiny{TM}}}}$ suture passers. Both ends of the Hi-Fi® suture are loaded through the eyelet of the PopLok® Anchor. A hemostat is then attached to the sutures and the anchor is allowed to hang outside the cannula.

Drill

The drill guide is placed just onto articular surface of glenoid. The drill is advanced until the proximal depth line is just below the surface of the subchondral bone.


Seat, Tension, POP, and Go!

Orient the suture eyelet on the anchor to face the labrum or repair site. The lever on the anchor handle stop makes contact with the drill is co-linear with the suture eyelet guide ensuring that the distal laser to aid in achieving this desired orientation. "Pretension" the suture to estimate the tension needed for the final repair. The PopLok® is then fully seated into the pilot hole using a mallet and the individual suture strands are then "micro-tensioned" as required. After disengaging the red safety lever, the anchor is deployed by squeezing the trigger, making an audible "POP". These steps are repeated for subsequent anchors.



INSTABILITY SURGICAL TECHNIQUES

Arthroscopic **SLAP Repair** using the 2.8 or 3.3mm PopLok® Knotless Suture Anchor (Left Shoulder)



Pass Sutures

Using a 45° Spectrum® II Hook and Super Shuttle® Relay suture passers, one end of a free #2 Hi-Fi® Suture is passed from the superior aspect of the labrum toward the articular surface just anterior to the biceps tendon and retrieved out an accessory anterior portal. The second end of the suture is then "shuttled" in a similar manner just posterior to the biceps tendon creating a "loop" just inferior to the biceps anchor. A looped grasper is then used to pull the two ends of the suture back through the loop, which is then tightened. The two ends of the suture are then loaded into the eyelet of the PopLok® anchor.



Drill

Place the drill guide just beyond the articular surface of glenoid at the midpoint of the biceps "anchor" on the superior glenoid tubercle. The drill is advanced until the proximal depth stop makes contact with the drill guide ensuring that the distal laser line is just below the surface of the subchondral bone.




Seat, Tension, POP, and Go!


Orient the suture eyelet on the anchor to face the labrum. The lever on the anchor handle is co-linear with the suture eyelet to aid in achieving this desired orientation. "Pretension" the suture to estimate the tension needed for the final repair. The PopLok® anchor is then fully seated into the pilot hole using a mallet and the sutures are individually micro-tensioned.

After disengaging the red safety lever, the anchor is deployed by squeezing the trigger. The PopLok anchor's audible "POP" confirms fixation of the suture and disengagement of the driver.

ROTATOR CUFF SURGICAL TECHNIQUES

Arthroscopic **Single Row Rotator Cuff Repair** using the 3.5 or 4.5mm PopLok® Knotless Suture Anchor – *Three Individual Technique Options*







Simple Stitch

Pass individual sutures through the torn tendon. Place paired suture limbs through a 4.5mm PopLok anchor and insert near the lateral edge of the footprint of the supraspinatus. For smaller tears only a single anchor and 2 sutures are used (a total of 4 suture limbs). For larger tears, 2 or even 3 anchors can be used. The anchors should be spaced at least 9mm apart. Suture organization is enhanced by utilizing individual (suture) stab incisions.

Horizontal Mattress (Inverted) Stitch

A horizontal mattress (inverted) stitch provides broader compression of the tendon onto the greater tuberosity. A single horizontal mattress stitch should be used with each 3.5mm anchor to prevent bunching of the tendon, which can occur when attempting to bring too many different parts of the tendon down to a single anchor point.

Mattress Loop Stitch

For this stitch, the sutures are passed similarly to a horizontal mattress except both ends of the suture are brought back through the trailing loop. The loop closes around the tendon as the suture is tensioned in the anchor. Only a single pair of suture limbs (or single loop) should be used per 3.5mm anchor to prevent bunching of the tendon during tensioning.

NOTE: For these single row techniques, see page 5 for additional information on "Seat, Tension, POP and GO!"



ROTATOR CUFF SURGICAL TECHNIQUES

Arthroscopic Double Row Rotator Cuff Repair using the 3.5 or 4.5mm PopLok® Knotless Suture Anchor



Suture Loading and Punch

With two CrossFT[™] suture anchors placed medially and mattress sutures tied, a single suture limb from each anchor is loaded into the first PopLok anchor. A hemostat is clipped on the suture limbs and the loaded PopLok anchor hangs from the cannula. To facilitate insertion of the anterior anchor through the lateral working portal, the shoulder may be externally rotated, so the perfect location can be achieved. The PopLok punch is then inserted down to the laser line for the anterior anchor.



Anterior Anchor

Using same insertion angle as the punch, the PopLok anchor is inserted down to the laser line. With counter pressure applied, each suture is tensioned individually, then locked into the cleat on the handle. After the red safety lever is disengaged, the anchor is deployed by squeezing the trigger, making the audible "POP".



Posterior Anchor

The shoulder is internally rotated to allow the posterior anchor to be inserted through the lateral working portal. The second suture from each medial row anchor is loaded into the second PopLok anchor. The same technique - hemostat, punch, anchor insertion, tension, "POP"- is repeated for this posterior anchor.



Final Construct

The Katana® Suture Cutter cuts the remaining suture tails and the completed repair can be viewed from both the posterior and lateral portals, confirming compression of the tendon to the footprint and appropriate placement of the anchors. A final view of the repair from the articular side is also recommended, to ensure medial row compression.

QUICK REFERENCE: Remember these four simple steps in every Technique



SEAT ANCHOR

Sutures are loaded into the anchor and the anchor seated into bone tunnel.

NOTE: A maximum of two #2 Hi-Fi® suture limbs should be used with the 3.5mm anchor, and four suture limbs with the 4.5mm anchor. The 2.8 and 3.3mm anchors will accept two #2 or four #0 limbs of HiFi® Suture.





TENSION

Sutures are pulled to proper tension and locked into cleat on handle.



POP

Trigger is squeezed until audible "POP", which locks the inner sleeve into the outer sleeve and deploys the wings.



GO

The sutures are locked between the sleeves. The inserter automatically detaches from the anchor. Suture tails are cut with the Katana® Suture Cutter.



ORDERING INFORMATION

SPECTRUM® II SET		POPLOK® KNOTLESS SUTURE ANCHOR	
Spectrum II Handle	C6350	INSTABILITY REPAIR	
Spectrum II Sterilization Tray		2.8mm PopLok w/one #2 Hi-Fi Suture	GKP-2801
Spectrum II Roller Wheel Replacement Kit		2.8mm PopLok w/two #0 Hi-Fi Sutures	
		3.3mm PopLok w/one #2 Hi-Fi Suture	
LIMITED REUSE SUTURE HOOKS		3.3mm PopLok w/two #0 Hi-Fi Sutures	
Suture Hook 45° Right	C6360	PopLok 2.8/3.3mm Drill Bit	
Suture Hook 45° Left		Drill Guide	
Suture Hook 60° Right		PopLok Obturator	
Suture Hook 60° Left			
Suture Hook 90° Right		ROTATOR CUFF REPAIR	
Suture Hook 90° Left		3.5mm PopLok Anchor (no suture)	CKP-3500
Suture Hook CorkScrew, Right		3.5mm PopLok Anchor w/one #2 Hi-Fi Suture	
Suture Hook CorkScrew, Left		4.5mm PopLok Anchor (no suture)	
Suture Hook Straight	C6368	4.5mm PopLok Anchor w/two #2 Hi-Fi Sutures	
Suture Hook Crescent,		3.5mm PopLok Punch	
Small, 3.0 x 15.0mm		4.5mm PopLok Punch	
Suture Hook, Crescent,		SRS Instrument Tray	RCR-TRAY
Medium, 4.0 x 20.0mm	C6370		
Suture Hook, Crescent,		HI-FI® HIGH STRENGTH SUTURE	
Large, 6.0 x 25.0mm	C6371	(STERILE, 12 PER BOX)	
		#2, 36 in. single strand, (blue and white co-braid) no needle	H5120
DISPOSABLE SUTURE HOOKS		#2, 36 in. single strand, (white) no needle	H5130
Suture Hook, 45° Right (Red)	C6380	#2, 36 in. single strand, (white and green co-braid) no needle	H5140
Suture Hook, 45° Left (Blue)	C6381	#2, 36 in. single strand, (black and white co-braid) no needle	
Suture Hook, 60° Right (Orange)		#0, 36 in. single strand, (black and white co-braid) no needle	H5301
Suture Hook, 60° Left (Yellow)		#0, 36 in. single strand, (white) no needle	H5302
Suture Hook, Straight (Pink)	C6384		
Suture Hook, Crescent, Small,		DRY-DOC® CANNULA	
3.0 x 15.0mm (White)	C6385	Dry-Doc 5x85mm	C7350
Suture Hook, Crescent, Medium,		Dry-Doc 6x85mm	
4.0 x 20.0mm (Teal)	C6386	Dry-Doc 7x85mm	
Suture Hook, Crescent, Large,		/	
6.0 x 25.0mm (Purple)	C638/	Dry-Doc 8x75mm	
ACCESSORIES		Dry-Doc 8x85mm	
Super Shuttle® Suture Passer (8/box)		REUSABLE CANNULATED METAL OBTURATOR	
Loop Handle Knot Pusher		5.0mm x 85mm	C7380
Crochet HookGrasping Forceps, 3.4mm	С6105	6.0mm x 85mm	
Grasping Forceps, 3.4mm Diameter, Straight with Ratchet	11 1001	7.0mm x 85mm	
Suture Retrieval Forceps, 3.4mm	. 11.1001	7.0mm x 95mm	C7382
Diameter	16 1018		
Diameter		8.0mm x 85mm	C7395
KATANA® HIGH-STRENGTH SUTURE	CUTTER		

orthroscopic Repairs Ising the PopLok[®] Inotless Suture Anchor



11311 Concept Boulevard Largo, Florida 33773-4908 Customer Service (800) 237-0169 USA Fax (727) 399-5256 International Fax (727) 397-4540

linvatec.com srs.linvatec.com