


UNIQUE BONE VOID FILLER

Tribio™ Backfill Plugs



Delivered sterile and ready to use

Osteoconductive scaffold rapidly absorbs the surrounding blood and fluids to aid in new bone formation



A GLOBAL EXTREMITY COMPANY

UNIQUE BONE VOID FILLER

Tribio™ Backfill Plugs

Tribio™ Backfill Plugs are a proprietary material delivered sterile and ready for use in any bony voids or gaps of the skeletal system. Consisting of hydroxyapatite, tri-calcium phosphate and bioactive glass in a collagen matrix, the Backfill Plugs are sized to easily fit in

surgically created osseous defects; 5.5mm dia. (for 6mm void) and 7.5mm dia. (for 8mm void) by 40mm in length. The osteoconductive scaffold rapidly absorbs the surrounding blood and fluids to aid in new bone formation during the healing process.

CAT NO.	DESCRIPTION	HOLE DIA.	IMPLANT SIZE
M80 SB006	Tribio™ Backfill Plugs	6mm	5.5mm x 40mm
M80 SB008	Tribio™ Backfill Plugs	8mm	7.5mm x 40mm





Tribio™ Backfill Plug Technique




1 After use of the CoLink® Bone Graft Harvester, flush and suction the subcutaneous tissue at the lateral heel incision to make sure no residual small pieces of graft are remaining.

2 Use the corresponding sized backfill plug to fill the void left by the drill or trephine (5.5mm Tribio™ Implant for a 6mm void, and 7.5mm dia. for an 8mm void).



3 Insert the Tribio™ Implant into the harvested hole with a slight twisting motion into the tunnel.

4 With the use of a hemostat or bone tamp, drive the Tribio™ Implant fully into the bone until it is flush, or slightly sub-flush



with the cortex.

Confirm the plug is not sitting proud and trim any excess material above the bone surface.

The material will rapidly absorb the surrounding blood and other fluids and will be replaced by bone over time. Repeat the flush and suction of the soft tissue to make sure no granules of the Tribio™ Implant are left behind.

Also available from In2Bones...

CoLink® Bone Graft Harvester


The In2Bones CoLink® Bone Graft Harvester is a single-use, pre-assembled bone graft harvesting device provided sterile that may be used to harvest bone from various sites in the body including the calcaneus, iliac crest, proximal tibia, distal tibia,

distal radius, and distal femur.

The Bone Graft Harvester is available in two sizes, 6mm and 8mm outer diameter.

CoLink® Bone Graft Harvester / Sterile, Single-use

CAT NO.....	DESCRIPTION	SIZE
G05 S1006....	CoLink® Bone Graft Harvester.....	6mm
G05 S1008....	CoLink® Bone Graft Harvester.....	8mm



All content contained herein is furnished for informational purposes only. In2Bones does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgment and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and precautions are listed in the implant package insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure. Availability of these products might vary from a given country or region to another as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.



A GLOBAL EXTREMITY COMPANY

Corporate Headquarters

In2Bones Global, Inc. • Memphis, TN • USA
844. 602. 6637 • Info@i2b-USA.com

International Office

In2Bones SAS • Lyon • France
+33 (0)4 72 29 26 26

FOLLOW US...



In2Bones.com

Tribio, CoLink, the In2Bones name and logo are trademarks of In2Bones or its affiliates. In2Bones USA Memphis TN, 38119 USA / In2Bones SAS, 69130 Ecully, France © 2018 In2Bones USA, Memphis, TN • All rights reserved • SP-16005-02 Rev D, 0920