

TriVay®
TibioTaloCalcaneal (TTC)
Arthrodesis System





TABLE OF CONTENTS

Indications / Contraindications2
System Overview / Specifications3
Instrumentation Design Features3
Surgical Technique Overview
Distal Screws First5
Proximal Screws First6
Surgical Technique
Patient Positioning / Bone Prep
Anterior Approach7
Trans-fibular Approach8
Plantar entry point9
Guide Wire insertion9
Intramedullary bone preparation 10
Targeting Jig assembly11
Nail Fixation and Compression
Distal Screws placement first13
Subtalar Compression Screw
Anterior approach15
Posterior approach16
Talar Screw placement17
Tibio-Talar Joint Compression18
Proximal Screw placement18
Proximal Screws placement first 21
Closing and postoperative protocol22
Revision technique22
Ordering information23
Implants
Instruments
Regulatory information24

INDICATIONS

The TriWay® Nail Ankle Arthrodesis System is intended for use in tibiotalocalcaneal arthrodesis and treatment of trauma to the hindfoot and distal tibia. Examples include:

- Post-traumatic and degenerative arthritis involving both ankle and subtalar joints
- Rheumatoid arthritis with severe deformity
- Revision of failed ankle arthrodesis with subtalar involvement or with insufficient talar body
- Revision of failed total ankle arthroplasty with subtalar intrusion
- Talar deficiency conditions, including avascular necrosis (requiring a tibiocalcaneal arthrodesis)
- Neuropathy or neuromuscular deformity or other neuromuscular disease with severe deformity or instability of the ankle, including Charcot foot
- Severe pilon fractures with trauma to the subtalar joint
- Malunited tibial pilon fractures

The addition of an IBS[®] 6.5mm Compression Screw through the subtalar joint and through the Nail is required.

CONTRAINDICATIONS

The implant should not be used in a patient who has currently, or who has history of:

- Acute or chronic, systemic inflammations,
- Active infections,
- Sensitivity/allergies to the implant materials.
- · Bone pathologies that may compromise the rigidity of the implant fixation (examples include: osteoporosis, acute cystic developments, acute osteopenia, bone tumor, etc...)

In2Bones[®] as the manufacturer of this device, does not practice medicine. The surgeon who performs any implant procedure is responsible for determining and using the appropriate surgical techniques for implanting the device in each patient. This Surgical Technique Manual is furnished for information purposes, as an aid to properly use the device and its dedicated instruments.

TriWay° TibioTaloCalcaneal Arthrodesis System SYSTEM OVERVIEW / SPECIFICATIONS / DESIGN RATIONALE















IMPLANT SPECIFICATIONS

Solid Nail Construct

Multi-planar fixation

- Medio-lateral fixation
- Postero-anterior fixation
- Oblique screw in subtalar joint line with a slight valgus to match subtalar anatomical alignment

7° posterior offset Nail to calcaneocuboid joint line

- Designed for increased purchase in the calcaneus
- Straight Nail philosophy for bone preparation surgical steps

Anatomical Design

- Right & left versions
- 3 diameters (10-11-12mm)
- 3 lengths (160-200-250mm)

SCREWS DESIGN RATIONALE

5.0 mm Cotter Screws

- Available in 25mm to 90mm lengths, 5mm increments
 - For tibia, talus and calcaneus fixation
 - Solid construct and works as a rod and allows for adequate fixation in bone
 - No screw head impingement
 - No screw thread / Nail impingement


5.0 mm Headed Screws

- Available in 20mm to 35mm lengths, in 2.5mm increments
- Headed option for proximal tibial fixation


6.5mm IBS C & C+ Screws

- Can be placed either with an anterior or posterior approach
- Compression screw provides subtalar compression
- When distal screws are placed first, there is no need to anticipate compression effect upon final position of the Nail
- Placed across the subtalar joint and through the Nail
- 6.5mm IBS C Screws: distal thread is 15mm long and is recommended when subtalar screw is preferred from posterior approach
- 6.5mm IBS C+ Screws: distal thread is 30mm long and is recommended when subtalar screw is preferred from anterior approach


INSTRUMENTATION DESIGN FEATURES



SURGICAL TECHNIQUE OVERVIEW...with Distal Screws First



TriWay® TTC Arthrodesis System SURGICAL TECHNIQUE OVERVIEW... with Distal Screws First



6. Distal screws insertion

5. Compression

SURGICAL TECHNIQUE: Anterior Approach







1. Patient Positioning and Bone Preparation with Anterior Approach

If an anterior approach is preferred, patient can be placed in supine position on the operating table.

A SUPINE PATIENT POSITIONING

The patient is placed in the supine position on a radiolucent operating table with the affected extremity elevated in neutral alignment. If the surgeon wishes to insert the subtalar compression screw with an anterior approach, this is the preferred position. General or spinal anesthesia is required and a thigh tourniquet is used.

3 SURGICAL EXPOSURE

An anterior incision is performed at the tibiotalar joint line level. Care should be taken with neurovascular structures and tendons.

A lateral approach is done in the sinus tarsi area to access the subtalar joint.

BONE PREPARATION

Osteophytes and cartilage are removed from the tibiotalar joint line with the use of curette, chisel, bone rongeur or sawblade, according to surgeon preference.

The Tibio Talar Cutting Guide (NO2 00023) can be used to perform parallel cuts over the tibiotalar joint line. Place the Cutting Guide block over the joint line and the 2.5x200mm K-wire (K10 NS200) into the proximal part of the Guide to align the Guide over the tibial axis. Once the Guide is well positioned, fix the Guide on the tibia and on the talus with 4 2.5x100mm K-wires (K10 NS251). Perform the tibial and talar parallel cuts with the use of a sawblade directly within the Guide slots. Remove the Guide once the cuts are completed and remove bone fragments.

The subtalar joint line is also prepared to ensure proper subtalar fusion.

SURGICAL TECHNIQUE... with Trans-fibular Approach

2. Patient Positioning and Bone Preparation with Trans-fibular Approach

(A) If a trans-fibular approach is preferred, patient can be positioned in the lateral position on the operating table. If a supine position is preferred, please refer to the supine position description page.

LATERAL PATIENT POSITIONING

The patient is positioned in the lateral decubitus position on a radiolucent operating table with the affected extremity accessible. If the surgeon wishes to address the subtalar compression screw with a posterior approach, this is the ideal position. General or spinal anesthesia is required and a thigh tourniquet is used.



A curvilinear incision is made over the fibula, curving distally along the peroneus tendon. Care should be taken with neurovascular structures and tendons. The incision is then extended distally to the sinus tarsi area to visualize the subtalar joint and facilitate access for bone preparation.



The distal part of the fibula is removed in an oblique fashion in order to expose the tibiotalar and the subtalar joints lines. The distal part of the fibula shall be used either for bone graft, or shall be repositioned at the end of the procedure to reproduce the mortise of the ankle joint.

BONE PREPARATION

Osteophytes and cartilage are removed from the tibiotalar joint line with a curette, chisel, bone rongeur or sawblade, according to surgeon preference. The subtalar joint line is also prepared to ensure good subtalar fusion.








SURGICAL TECHNIQUE



3. Plantar Entry Point

A The plantar entry point is visually localized with a K-wire placed in the tibial axis and in the alignment of the second ray. Due to the posterior offset design of the implant, the entry point is moved 3 to 5mm posteriorly. A longitudinal incision, approximately 2-3cm should be made at this intersection. Careful dissection is then utilized to gain access to the plantar surface of the calcaneus.



4. Guide Wire Insertion

A Insert the 3.2x350mm Guide Wire (N02 00014) through the center of the lateral column of the calcaneus up to the center of the talar dome.

Once the 3.2x350mm Guide Wire is located in the center of the talar dome, care must be taken to position the ankle in 5° dorsiflexion before Guide Wire insertion in the tibial canal shaft. Insert the Guide Wire in the tibia.

Confirm the Guide Wire position with a C-Arm. Guide Wire must be aligned with the talar lateral process.




Take care to position the ankle in 5° dorsiflexion to avoid equinus position after Nail insertion.

SURGICAL TECHNIQUE

5. Calcaneal, Talar and Tibial Bone Preparation

A Place the Soft Tissue Protector (N02 00006) over the 3.2x350mm Guide Wire (N02 00014) to avoid any soft tissue impingement. Insert the 7.0mm cannulated Drill (NO2 00027) over the Guide Wire and reamuntil the laser marking is aligned with the base of the Soft Tissue Protector.

Remove the 7.0mm cannulated Drill. Take care to keep the 3.2mm Guidewire in place during the removal process. B Insert the 14.0mm Reamer (N02 00016) over the 3.2mm Guidewire. Use the Soft Tissue Protector to avoid any soft tissue impingement. Ream untill the laser marking is aligned with the base of the Soft Tissue Protector.




TIBIAL CANAL REAMING

Tibial canal must be prepared with cannulated flexible reamers. Cannulated flexible Reamers and a Bead-Tipped Guide Wire (C01 S105) can be provided upon request by In2Bones as an option to the general instrument set.



Cannulated flexible reamers must not be used over the 3.2x350mm Guide Wire (N02 00014). Choice of the reamer guide must be done following the flexible reamers' manufacturer recommendation. Be sure to follow manufacturer's recommendation while using the flexible reamers.

• Ream the tibial canal 1mm above the desired Nail diameter size. Failure to do so would create excessive resistance to nail insertion. Nail diameter and length will be selected according to surgeon preference.



6. Targeting Jig Assembly

A Screw in the Index Lock (N02 00020) into the Index Tab (N02 00005) until its threads appear in the index tab window.

③ Slide the Index Tab, cylindrical part first, into the Targeting Jig (N02 00001) with the Index Lock on the opposite side of the Locked / Free laser-marks, until the laser line on the Index Tab is aligned in the "Free" position.

⊙ Engage the Compression Sleeve (N02 00010) into the Compression Wheel slot (N02 00009) with the domed surface up.



• Align the flat surfaces of the Compression Sleeve and the Connector (N02 00003) facing up.



Screw the Sleeve and the Connector together. Make sure the flat surfaces face up during this step.

● Attach the Nail assembly to the Connector and lock the Fixation Screw (N02 00004) with the non cannulated T25 Screwdriver (G01 01641).

The Nail must be tightened firmly.


- Insert the Nail assembly within the Targeting Jig. Make sure the Index Tab is in the "Free" position.
- 1 Push and maintain the Index Lock.
- 2 Place the Nail assembly within the Index Tab/Targeting Jig as shown.
- 3 Pull up slightly on the Nail assembly until a visible and audible click of the Index Tab is noticed to confirm proper alignment.



After the assembly, the Compression Wheel should be near the bottom of the Jig.




Confirm that the Nail is angled toward the Jig.



SURGICAL TECHNIQUE

G To lock, push the Index Tab in the Targeting Jig until the "Locked" position is reached with the Index Tab.

1 To modify the Nail position, pull the Index Tab until the line faces "Free" on the Targeting Jig.




• Rotate the Nail by 1/4 of a turn to obtain the desired position



POSTERIOR POSITION

LATERAL AND MEDIAL POSITION

- 1 Push the Index Tab back into "Locked".
- 2 Screw the Hitting Plate (N02 00011) onto the Targeting Jig.



7. Nail Fixation and Compression with Distal **Screws First**

NAIL INSERTION

Make sure the 3.2x350mm Guide Wire is removed before inserting the Nail.

A Place the Nail in the back position on the Targeting Jig with the Index Tab on the "Free" position, by aligning the Targeting Jig arm with the back position as illustrated.




Lock the position by pushing the Index Tab to the "Locked" position.

The Nail is introduced manually in the prepared bone surfaces through the plantar incision. Confirm the insertion depth with fluoroscopy.

3 A radiopaque demarcation is visible to visualize the distal part of the Nail. Be sure the Nail's distal part is flush to the calcaneus plantar aspect. If needed, the Nail can be further impacted with a mallet on the Hitting Plate or with the Slap Hammer Guide (N02 00024) and the Slap Hammer Mass (N02 00025) assembled and threaded into the Hitting Plate.

Confirm the final positioning of the Nail with fluoroscopy.





DISTAL SCREWS INSERTION

(NO2 00013) and the Drill Guide (NO2 00007) into the calcaneus screw hole of the Targeting Jig. Insert the K-wire Guide (N02 00008) into the subtalar screw hole of the Targeting Jig.

SUBTALAR SCREW PREPARATION


Insert the 2.0x300mm K-wire (N02 00017) into the K-wire Guide until the K-wire protrudes from the talar neck.



Make sure to select the appropriate hole of the Targeting Jig, according to the operative side.




Make sure that the Sleeve and Drill Guide are in contact with the bone for accurate measurement.



SURGICAL TECHNIQUE

CALCANEAL SCREW INSERTION

• Insert the Trocar Tip (N02 00026) into the Protecting Sleeve and Drill Guide and puncture the calcaneus cortical wall by impacting the Trocar Tip with a mallet. Remove the Trocar Tip.




(3) Insert the 4.3mm non cannulated Drill (N02 00015) in the Protecting Sleeve and the Drill Guide. Drill the calcaneus until the calcaneo-cuboid joint line. Control the Drill position with fluoroscopy.




Care must be taken to stop drilling before the drill enters the cuboid bone.

Read the screw length directly on the Drill over the Drill Guide with the Drill in place in the bone. Make sure the Protecting Sleeve is pressed firmly to the bone before reading screw length.

Remove the Drill and the Drill Guide from the Protecting Sleeve.



G Use the non cannulated T25 Screwdriver (G01 01641) to insert the appropriate length Cotter Screw in the calcaneus.



(3)

SUBTALAR COMPRESSION SCREW INSERTION WITH ANTERIOR APPROACH

The subtalar compression screw can either be inserted with a posterior approach or an anterior approach.

- ① In case an **anterior approach** is preferred, be sure the 2.0mm K-wire already inserted in the subtalar joint protrudes from the talar neck in order to insert the cannulated 4.3x275mm Drill (G01 01651). Drill with the cannulated Drill from the talus to the calcaneus over the K-wire.
- Direct visualization of the screw length can be read on the cannulated Drill at the bone contact with the cannulated Drill. Remove the Drill.



Remove 5mm from the measurement to select the appropriate screw length as the compression screw will compress the subtalar joint line during its insertion.


A 6.5mm IBS C+ Screw is preferred in order to have a longer threaded portion in the calcaneus.

Insert the appropriate 6.5mm IBS Screw with the cannulated T25 Screwdriver (G01 00411) over the K-wire. Once the screw placement is confirmed with fluoroscopy, remove the K-wire.





Make sure to select the appropriate hole of the Targeting Jig, according to the operative side.







SURGICAL TECHNIQUE

SUBTALAR COMPRESSION SCREW INSERTION WITH POSTERIOR APPROACH

The subtalar compression screw can either be inserted with a posterior approach or an anterior approach.

In case a **posterior approach** is preferred, remove the K-wire Guide while keeping the K-wire in place. Insert the Protecting Sleeve and the Drill Guide into the subtalar screw hole of the Targeting Jig. Take care to have the K-wire located within the Drill Guide.

Insert the cannulated 4.3x275mm Drill (G01 01651) into the Drill Guide and over the K-wire. Drill with the cannulated Drill from the calcaneus to the talus over the K-wire.

Direct visualization of the screw length can be read on the cannulated Drill at the Drill Guide contact with the cannulated Drill. Remove the Drill and the Drill Guide.



Remove 5mm from the measurement to select the appropriate screw length as the compression screw will compress the subtalar joint line during its insertion.


A 6.5mm IBS C Screw is preferred in order to have a shorter threaded portion in the talus.


M Insert the appropriate 6.5mm IBS Screw with the cannulated T25 Screwdriver (G01 00411) over the K-wire. Once the screw placement is confirmed with fluoroscopy, remove the K-wire.





Make sure to select the appropriate hole of the Targeting Jig, according to the operative side.





TALAR SCREW INSERTION

The talar screw can either be inserted from a lateral or medial approach based on surgeon preference.

Rotate the Targeting Jig 90° to the desired side by pulling the Index Tab to the « Free » position. Once set, lock the position by pushing the Index Tab to the « Locked » position.

① Insert the Protecting Sleeve and the Drill Guide into the talar screw hole of the Targeting Jig.



Make sure to select the appropriate hole of the Targeting Jig, according to the operative side.




Make sure that the Sleeve and Drill Guide are in contact with the bone for accurate measurement.

Place the Trocar Tip in the Drill Guide (N02 00026) and puncture the cortex by impacting the Trocar Tip with a mallet. Remove the Trocar Tip.

The screw length can either be read over the 4.3mm non cannulated Drill or with the Depth Gauge (N02 00018). If the Depth Gauge is preferred, remove the 4.3mm non cannulated Drill and the Drill Guide before inserting the Depth Gauge.

- Drill with the 4.3mm non cannulated Drill in the Drill Guide. The talar screw must be placed bicortically. Remove the 4.3mm non cannulated Drill.
- Remove the Drill Guide from the Protecting Sleeve and insert the Screwhead Countersink. Drill the Countersink to prepare the cortex.
- Select the appropriate length Cotter Screw and insert it with the non-cannulated T25 Screwdriver within the Protecting Sleeve.



SURGICAL TECHNIQUE

TIBIO-TALAR JOINT LINE COMPRESSION

Manual compression is achieved to compress the tibio-talar joint when the distal screws are inserted first.

Manual compression of the tibiotalar joint can be achieved directly by pushing or hammering the Hitting Plate or by using the Slap Hammer. Hammering motion needs to be applied in alignment with the nail, failure to do so might lead to breakage of the hitting plate. Hammering motion needs to be applied in alignment with the nail, failure to do so might lead to breakage of the hitting plate.

TIBIAL SCREWS INSERTION

The Targeting Jig arm is positioned on the lateral or medial side of the tibia according to surgeon preference.

Rotate the Targeting Jig 90° to the desired side by pulling the Index Tab to the « Free » position. Once set, lock the position by pushing the Index Tab to the « Locked » position.

The posterior holes are always used to target the Nail.



both of the Targeting Jig arm tibial holes. Be sure to select the approriate holes according to the size of the Nail inserted in the tibia. For 160mm length Nails use the 2 most distal holes. For the 200mm length Nails use the 2 most proximal holes. For the 250mm Nails free hand targeting will be performed.

M Place the Trocar Tip in the most distal Drill Guide and puncture the cortex by impacting the Trocar Tip with a mallet. Remove the Trocar Tip.




Make sure that the Sleeve and Drill Guide are in contact with the bone for accurate measurement.



non cannulated Drill in the Drill Guide. Screw length can be read over the 4.3mm non cannulated Drill.

With the Drill still in place in the distal Drill Guide, place the Trocar Tip in the proximal Drill Guide and puncture the cortex by impacting the Trocar Tip with a mallet. Remove the Trocar Tip.



SURGICAL TECHNIQUE

2 Remove Drill Bit from the most distal hole. Drill the proximal tibial screw hole with the 4.3mm non cannulated Drill in the Drill Guide. The tibial screws must be placed bicortically. The screw length can be read over the 4.3mm non cannulated Drill.

AA Remove the Drill Guide from the Protecting Sleeve of the proximal tibial hole and insert the Screwhead Countersink. Drill the Countersink to prepare the cortex.


BB Select the appropriate length Cotter or Headed Screw and insert into in the proximal tibial screw hole within the Protecting Sleeve using the non cannulated T25 Screwdriver.

Remove the Drill Guide from the Protecting Sleeve of the proximal tibial hole and insert the Screwhead Countersink. Drill the Countersink to prepare the cortex.

D Select the appropriate length Cotter or Headed Screw and insert it into the distal tibial screw hole within the Protecting Sleeve with the non cannulated T25 Screwdriver.



The tibial screws must be placed bicortically.




8. Nail Fixation and Compression with Proximal Screws Insertion First


⚠ Insert the K-wire Guide (N02 0008) into the more distal proximal tibial screw hole (oblong) within the Targeting Jig. Insert the K-wire into the Guide.

Depth of insertion should account for the amount of compression (up to 35mm) to be applied after placing the most proximal Cotter or Headed Screw through the nail.

Note: Leave the K-wire in place to keep a closed frame between the nail and the arm of the Targeting Jig.

- (3) Insert the Protecting Sleeve and the Drill Guide within the most proximal tibia screw hole within the arm of the Targeting Jig. Place the Trocar Tip in the proximal Drill Guide and puncture the cortex by impacting the Trocar Tip with a mallet. Remove the Trocar Tip.
- © Drill the most proximal tibia screw hole with the 4.3mm non-cannulated Drill in the Drill Guide. The screw length can be read over the 4.3mm non-cannulated drill.
- Remove the Drill Guide from the Protecting Sleeve of the most proximal tibial screw hole and insert the Screwhead Countersink. Drill the Countersink to prepare the cortex.
- © Select the appropriate length Cotter or Headed Screw and insert it into the most proximal tibial screw hole through the Protecting Sleeve with the non-cannulated T25 Screwdriver.
- Remove the K-wire, K-Wire Guide, and Protecting Sleeve from the Targeting Jig.
- **G** Insert the Protecting Sleeve and Drill Guide into the more distal of the proximal tibial screw holes. Following the same process, drill and countersink the hole, then fill with the appropriate Cotter or Headed Screw.
- ① Turn the wheel clockwise to apply compression in the subtalar joint and the tibio-talar joint. Be sure the Nail is not protruding into the calcaneal plantar aspect after compression is achieved.





SURGICAL TECHNIQUE

9. Closing and Postoperative Protocol

A Confirm hardware placement with AP and lateral X-rays, before detaching the Nail from the Jig. Remove the Jig with the Screwdriver by unscrewing the internal holding Fixation Screw. The wound is closed in layers and the patient is treated with additional fixation if appropriate.

The patient should remain non-weight-bearing for six weeks or until there is early radiographic evidence of consolidation at the arthrodesis sites.

Further protection with a walking boot or brace may help ease the transition to weight bearing.



X-rays provided by Stephen Kearns, MD, FRCS (Tr & Orth), Galway University Hospitals


10. Removal Technique


Connect the Slap Hammer Guide (N02 00024) to the distal part of the Nail.

A Remove all the Cotter and Headed Screws with the non cannulated T25 Screwdriver, starting from proximal to distal. In case a Screw is locked in the bone, use the Screw Extractor (N02 00019). Screw the Screw Extractor over the head of the Screw and pull on the Screw Extractor, or use the non cannulated T25 Screwdriver to remove the Screw.

Confirm with fluoroscopy if the subtalar screw has been inserted from anterior or posterior. Insert the 2.0mm K-wire in the subtalar screw. Remove the subtalar 6.5mm IBS Screw with the cannulated T25 Screwdriver.

13 Once all screws have been removed, place the Slap Hammer Mass (N02 00025) over the Slap Hammer Guide and gently extract the Nail by hitting on the Slap Hammer.





ORDERING INFORMATION

TriWay® Artl	nrodesis Nails / Sterile
PART#	DESCRIPTION RIGHT / LEFT
N20 ST101	10mm Lg., 160mm., R
N20 ST201	10mm Lg 160mm L
N20 ST102	10mm Lg., 200mmR
N20 ST202	10mm Lg., 200mm., L
N20 ST111	. 11mm Lg 160mm R
	. 11mm Lg 160mm L
N20 ST112	. 11mm Lg 200mm R
	. 11mm Lg 200mm L
	. 11mm Lg 250mm R
	. 11mm Lg 250mm L
	. 12mm Lg 160mm R
N20 ST221	. 12mm Lg 160mm L
N20 ST122	. 12mm Lg 200mm R
	. 12mm Lg 200mm L
	12mm Lg., 250mm., R
	. 12mm Lg 250mm L
N20 ST223	12mm Lg., 250mm., L

TriW	Vay [®] Cotter Scre	ws / Sterile
PART		
N20	ST5255.0mm	Lg25mm
N20	ST5305.0mm	Lg30mm
N20	ST5355.0mm	Lg35mm
N20	ST5405.0mm	Lg40mm
N20	ST5455.0mm	Lg45mm
N20	ST5505.0mm	Lg50mm
N20	ST555 5.0mm	Lg55mm
N20	ST5605.0mm	Lg60mm
N20	ST5655.0mm	Lg65mm
N20	ST5705.0mm	Lg70mm
N20	ST5755.0mm	Lg75mm
N20	ST5805.0mm	Lg80mm
N20	ST5855.0mm	Lg85mm
N20	ST5905.0mm	Lg90mm

TriWay [®] Headed Screws / Sterile			
PART#	DESCRIPTION		
N20 ST720	.5.0mm Lg20mm		
N20 ST722	. 5.0mm Lg22.5mm		
N20 ST725	. 5.0mm Lg25mm		
N20 ST727	. 5.0mm Lg27.5mm		
N20 ST730	. 5.0mm Lg30mm		
N20 ST732	. 5.0mm Lg32.5mm		
N20 ST735	. 5.0mm Lg35mm		
IBS _® C+ 6.5mm Sterile			
Compression Screw			

Compression Screw			
PART#	DESCRIPTION		
	6.5mm Lg60mm		
	6.5mm Lg65mm		
	6.5mm Lg70mm		
	6.5mm Lg75mm		
	6.5mm Lg80mm		
	6.5mm Lg85mm		
S65 ST590	6.5mm Lg90mm		


Compression Screw

PART#	DESCRIPTION
S65 ST160	6.5mm Lg 60mm
S65 ST165	6.5mm Lg 65mm
S65 ST170	6.5mm Lg 70mm
S65 ST175	6.5mm Lg 75mm
	6.5mm Lg 80mm
S65 ST185	6.5mm Lg 85mm
S65 ST190	6.5mm Lg 90mm
S65 ST195	6.5mm Lg 95mm
S65 ST300	6.5mm Lg100mm

IBS® C 6.5mm Sterile

INSTRUMENT SET

PART	NUMBER	DESCRIPTION QUANTITY
1	N02 00013	Protecting Sleeve2
2	N02 00003	Connector 1
3	N02 00004	Fixation Screw1
4	N02 00008	Guide for 2.0mm Dia K-wire2
5	N02 00024	Slap Hammer Guide1
6	N02 00025	Slap Hammer Mass1
7	N02 00011	Hitting Plate1
8	N02 00023	Tibio Talar Cutting Guide1
9	G01 00411	T25 Screwdriver Cannulated1
		(use over 2.0mm K-wire)
9	G01 01641	T25 Screwdriver Non Cannulated1
10	G01 01391	"T" Handle For Medium Screwdriver1
11	N02 00009	Compression Wheel1
12	N02 00001	Targeting Jig1
13	N02 00020	Index Lock1
14	N02 00006	Soft Tissue Protector1
15	N02 00019	Screw Extractor1
16	N02 00010	Compression Sleeve1
17	N02 00005	Index Tab1
18	G01 40041	K-wire Tube1
18	N02 00017	2.0mm Dia K-wire Lg 300mm3
18	K10 NS200	2.5mm Dia K-wire Lg 200mm 1 Sharp Tip
		Non Sterile1
18	K10 NS251	2.5mm Dia K-wire Lg 100mm 1 Sharp Tip
		Non Sterile4



PART	NUMBER	DESCRIPTION QUANTITY
19	N02 00026	Trocar Tip 4.3mm Dia1
20	N02 00027	7.0mm Dia Cannulated Drill
21	N02 00016	14.0mm Dia Cannulated Reamer1
22	N02 00014	3.2mm X 350mm Guide Wire
23	N02 00015	4.3mm Non Cannulated Drill
24	G01 01651	4.3mm x 275mm Cannulated Drill (use over 2.0mm K-wire) 1
25	N02 00018	Depth Gauge1
26	G01 01591	T25 Screwdriver Tip Cannulated K-wire 2.0mm 1 (use over 2.0mm K-wire)
27	N02 00021	Screwhead Countersink1
28	N02 00007	4.3mm Drill Guide

ADDITIONAL INSTRUMENTS

PART NUMBER C01 S1025	DESCRIPTION Bead-tipped Guide Wire, 5mm x 500mm,	QUANTITY Sterile1	
FLEXIBLE REAMERS			
BJ1103B5	Flexible Reamer, 9.0mm	1	
BJ1103B6	Flexible Reamer, 9.5mm	1	
BJ1103B7	Flexible Reamer, 10.0mm	1	

Flexible Reamer, 10.5mm1

Flexible Reamer, 11.0mm1

PART NUMBER	DESCRIPTION	QUANTITY
BJ1103B10	Flexible Reamer, 11.5mm	1
BJ1103B11	Flexible Reamer, 12.0mm	1
BJ1103B12	Flexible Reamer, 12.5mm	1
BJ1103B13	Flexible Reamer, 13.0mm	1

BJ1103B8

BJ1103B9



Regulatory Information

RECOMMENDATION

It is recommended to carefully read the instructions for use available in the package insert.

REIMBURSEMENT

Reimbursement may vary from country to country. Check with local authorities.

MANUFACTURER

In2Bones SAS 28, chemin du Petit Bois 69130 Ecully — FRANCE Tél: +33 (0)4 72 29 26 26 Fax: +33 (0)4 72 29 26 29

DOCUMENT

Reference: ST-DIG-TRIWAY-US-EN-082023

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

Always refer to the appropriate instructions for use for complete clinical instructions. Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.

Except otherwise specified, all commercial trademarks used in this document, In2Bones and the In2Bones logo are trademarks or registered trademarks of In2Bones or its subsidiaries.

CONMED Foot & Ankle

In2Bones USA 6000 Poplar Ave., Suite 115 Memphis, TN 38119 844 602 6637 www.conmed.com/in2bones In2Bones SAS 28, chemin du Petit Bois 69130 Ecully France +33 (0)4 72 29 26 26 CONMED Corporation 11311 Concept Blvd Largo, FL 33773 USA

[™] and ® denote Trademarks and Registered Trademarks of CONMED or its affiliates. Patent pending. ©2023 CONMED Corporation. All Rights Reserved.

All content contained herein is furnished for informational purposes only. In2Bones does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgment and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and precautions are listed in the implant package insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure. Availability of these products might vary from a given country or region to another as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.