



Enhance[®]

Pre-Shaped Allograft Wedges

The ENHANCE[®] line provides pre-shaped calcaneal lengthening and osteotomy wedge allograft solutions for foot deformity correction procedures

**MTF
ADVANTAGE**

No Terminal Radiation
No Harsh Chemicals

Over 9 million grafts distributed



In Partnership with:



HALL[®] MICROFREE[®]
CORDLESS POWER SYSTEM



ENHANCE[®]
INSTRUMENT TRAY



ENHANCE[®]
BONE TAPS



To learn more about these and other innovative products, call **1-866-4CONMED** or visit **www.CONMED.com**.

PRE-SHAPED ALLOGRAFT WEDGE SOLUTIONS FOR DEFORMITY CORRECTION



Enhance Instrument Tray

Allograft Wedges

Pre-Shaped Wedges from CONMED in partnership with MTF Biologics provide surgeons with convenient and reliable allografts for flatfoot deformity correction.

- Dense cancellous bone provides an excellent reincorporation environment
- Two-side cortical rim provides structural strength
- Pre-sized design eliminates the need for intraoperative shaping
- Cotton and Evans technique specific shapes

Note: Allograft wedges may be trimmed from the cancellous sides to customize length and width.

EVANS WEDGE



H	H1	L	W
6mm	4mm	22mm	22mm
8mm	5mm	22mm	22mm
10mm	5mm	22mm	22mm
12mm	5mm	22mm	22mm
14mm	6mm	22mm	22mm

COTTON WEDGE



H	L	W
5mm	22mm	16mm
6mm	22mm	16mm
7mm	22mm	16mm
8mm	22mm	16mm

ARCH WEDGE



H	L	W
5mm	22mm	16mm
6mm	22mm	16mm
7mm	22mm	16mm
8mm	22mm	16mm

Instrumentation

EVANS TRIAL



BONE TAMPS



INSTRUMENT SET



DISTRACTOR



COTTON TRIAL



Calcaneal Lengthening Procedure



1 Opening osteotomy prepared with sagittal saw and osteotome.



2 Converging K-wires (max. of 2.8mm) are utilized with pin distractor.



3 Trial Sizers represent allograft wedge dimensions.



4 Reconstitute the graft utilizing distractor, implant pre-shaped wedge.
Note: Graft length and width may be modified for improved fit



5 Final graft position can be manipulated by gently tamping the graft.



6 Completed procedure with proper graft placement
Surgeon discretion to be used for graft fixation.



Processed by:
 Musculoskeletal Transplant Foundation
 125 May Street
 Edison, New Jersey 08837
 www.mtfbiologics.org

Ordering Information

Description	Catalog Number	Description	Catalog Number	Description	Catalog Number
Evans Wedge		Cotton Wedge		Arch Wedge	
6mm	403106	5mm	4013205	5mm	403305
8mm	403108	6mm	4013206	6mm	403306
10mm	403110	7mm	4013207	7mm	403307
12mm	403112	8mm	4013208	8mm	403308
14mm	403114				
Allograft Wedge Instrumentation					
Instrument Set, Full					903000
	Osteotome, 16mm	Osteotome, 10mm	Mallet		
	Distractor	Bone Tamp, Large	Bone Tamp, Small		
	Evans Trial, 6mm	Evans Trial, 8mm	Evans Trial, 10mm		
	Evans Trial, 12mm	Cotton Trial, 5mm	Cotton Trial, 6mm		
	Cotton Trial, 8mm				
Instrument Tray, Empty					903001

Tissue Safety and Integrity

MTF is a non-profit organization founded in 1987. Since its inception, MTF has recovered over 139,000 donors and distributed over 9 million grafts for transplantation while maintaining an exemplary safety record. Because of all of the safety steps instituted by the Board of Directors, Medical Board of Trustees, and Donation Board of Trustees, MTF has become the number one tissue bank in the nation and is one of the largest providers of grafts in the world.

- MTF employs the latest technologies such as Nucleic Acid Testing and conducts serological tests on all donors.
- MTF leads the industry in donor selection criteria with standards that exceed those set by the AATB.

To learn more about this and other innovative products, call **1-866-4CONMED** or visit www.CONMED.com.

Toll Free: 1-866-4CONMED
 International: 1-727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com