



## Latex Information

To keep our customers aware of any possible or potential reactions and/or side effects from the use of our products, we provide the following list of devices and device components that are known to contain natural rubber latex. These substances may cause allergic reactions in some individuals.

Product	CAT #	Latex Component
Steigmann-Goff™ Endoscopic Ligator Kit	000221, 000227, 200220	“O” ring
Steigmann-Goff™ Endoscopic Ligator	100225	“O” ring
DuraGlide3™ Stone Removal Balloon	007123-007125, 007143-007145	Balloon
Arthroscopy Tube Set (83000 Pump)	83100	Tube

While there are no alternative latex free tubing sets able to function with the 83000 Arthroscopy Pump System, our other pumps (Apex® Universal Irrigation System, 87000 Pump System, 10K® Pump system, 24K® Pump System, and CrystalView™ Pro Fluid Management System) utilize tubing sets that are not made with natural rubber latex.

The following kits contain latex ONLY in the plunger of the syringe used for infusing lidocaine, sterile water, medications, or feeding formula:

Product	CAT #
Monopty®	000402, 000439, 000441

Please contact us if you have any questions.