




DRSHD™ 1080p

High Definition Image Management Solution

Product Specifications

Tablet Compatible



Still Capture:	1280 x 720 pixels, BMP, JPEG and PNG file format	Environmental: (Operating)	Temperature: 16.7 to 26.7° Temperature: 62° to 80° F Relative Humidity: 30 to 60% Atmospheric: 700 to 1060hPa
Motion Video Capture:	1080p (25/30fps); 720p (50/60fps); 640 x 480 25/30 fps AVCHD H.264 .mp4 files	Environmental: (Storage)	Temperature: -6.6 to 50° Temperature: 20° to 122° F Relative Humidity: 10 to 90% Atmospheric: 500 to 1060hPa
Hard Drive Capacity:	2TB Total capacity, 1.77TB for storage	Overall Dimensions:	CPU: 6in H x 14.0in W x 14.0in D 15.9cm H x 35.5cm W x 35.5cm D
Inputs:	DVI-I x 1: 720p, 1080i, 1080p, 50/59.94/60hz S-Video x 1, PAL/NTSC Composite x 2, PAL/NTSC Audio: Mic and Line level Foot Pedal: RS-232 Hand Switch: RS-232	Touch Screen:	14in H x 13.8in W x 8in D (Panel and base) 35.5 cm H x 35 cm W x 20.3 cm D
Outputs:	Audio Output: Line level	Weights:	CPU: 9.3kg / 20.0 lb Touch Screen: 7.3kg / 16.1lb
Data I/O:	USB 2.0: Front: Type 1 x 1, Rear: Type 1 x 4 Ethernet: Gigabit RJ45 x 1 Keyboard/Mouse I/O: USB Operating System: Windows 7 Embedded GPIO: RJ-11		
Power Input:	CPU: 100 - 240 VAC, 50/60 Hz		

To learn more about this and other innovative products,
call **1-800-237-0169** or visit **ConMed.com**.



DRS02

UP-DR80MD



VP8500
Video Cart with
VP8520
Secondary Arm

Ordering Information

Description

Catalog Number

CPU, High Definition Digital Recording System	DRSHD 1080p
Touch Screen, 15in, High Definition Digital Recording System	DRS02
Footswitch, 2 Pedal, Record and Still Capture	DRS03
GPIO Cable (for video monitor recording indicator)	DRS28
Power Supply for DRS02 Touch Screen	DRS11
Extended Cable Kit for DRSHD Touch Screen, 10ft (3m)	DRS19
Keyboard, USB, for DRSHD	DRS29
Cable, DVI-D, 6ft	DRS25
WiFi Kit with Installation Card	DRS31
Secondary Arm (for VP8500 Video Cart)	VP8520
Sony Digital Color Printer	UP-DR80MD
Print Pack, for UP-DR80MD, Letter Size for U.S. and Canada	UPC-R81MD
Print Pack, for UP-DR80MD, A4 Size for Europe	UPC-R80MD

DRSHD 1080p
DRS02
DRS03
DRS28
DRS11
DRS19
DRS29
DRS25
DRS31
VP8520
UP-DR80MD
UPC-R81MD
UPC-R80MD

Local: 727-392-6464
Toll Free: 800-237-0169

ConMed.com
customer_service@livatec.com