



Complete Arthroscopic Capital Solutions

FLUID MANAGEMENT, RESECTION, ABLATION, VISUALIZATION

From performance to simplicity, this portfolio delivers a wide range of benefits for healthcare providers and the patients they treat.


CRYSTALVIEW™ PRO



EDGE®



LOOKING GLASS®



TOUCH GLASS



ARTHROSCOPIC
RESECTION SYSTEM



Imagine, what we
could do together!

Contact us at customerexperience@CONMED.com,
internationalorders@CONMED.com or visit www.CONMED.com

Ablation

CONMED has the Edge®



Edge® Bipolar Arthroscopic RF System (AES-1)

Offers a wide range of performance benefits ranging from optimized power settings to joint temperature monitoring. The Edge® System delivers efficient and effective resection, ablation, and coagulation performance across all disciplines. The strong and stable plasma field efficiently and effectively removes tissue debris that can impair a surgeon's field of view. Bubbles are effectively removed by the optimized suction slot in the electrode. The edge of the large non embedded electrode can be used for sculpting and fine tuning with maximum control.

Benefits

- Plug & play performance
- Faster bulk tissue ablation with less risk
- Blade and sculpt with one versatile probe
- Real-time intraarticular joint temperature monitoring

Fluid Management

Focus on the Patient, not the Pump



CrystalView® Pro Fluid Management System (CV4000)

Combines high performance inflow and outflow for optimal distension and clear visualization with intuitive features and customizable surgeon presets to help simplify the OR. CrystalView™ Pro is designed to let you focus on your patient, not on the pump.

Benefits

- Plug & play joint preset
- Simple & intuitive touch screen
- Reliable fluid management
- Programmable surgeon settings

Resection



Arthroscopic Resection System (D4000A)

CONMED's next generation of powered instrument consoles for arthroscopic resection, the D4000A has intuitive features that deliver new levels of simplicity, versatility, and performance. It supports 31 handpieces including CONMED Shavers and Hall® Power Tools.

Benefits

- Plug & play performance
- All joints in one hand piece
- All your blade & bur requirements

Looking Glass®

The Performance You Want, the Durability You Need.

Integrated Surgical Visualization System (AV4132)

Has a new look and feel compared to traditional monitors. Remember this is not just a monitor, it's an integrated system housing the camera controller and power source. MAGIC is in the monitor.

Touch Glass Image Capture (AV41SIC)

Touch Glass is designed to be scalable, giving you options for functionality and cost.

Benefits

- 4K Ultra-High Definition (UHD) image quality
- Plug & play performance
- Reduced OR footprint
- Ultra durable eLuma™ electronic light guide (ELG)

CAMERA HEAD



ELUMA™





Ordering Information

RESECTION

Description	Qty	Catalog No.
Consoles		
D4000 Drive System	1	D4000
D4000A Drive System with Irrigation	1	D4000A
Blades and Burs		
Stealth™ Diamond-Like Carbon (DLC) Coated Blades	1	SDPS-C002
Lanza™ Tapered Speciality Shaver Blades	1	9391A
Mimix™ Small Joint Blades and Burs	1	LE-29FRP
Signature Series™ Blades and Burs	1	HPS-C001
Ultra Series™ Dual Purpose Shaver Blades	1	DPS-C001
Standard Blades and Burs	1	HPS-CB01
Handpieces		
Ergo™ 2-Buttons Shaver Handpiece	1	D4240
Ergo™ No-Button Shaver Handpiece	1	D4200
Foot Pedal	1	C9863

ABLATION

Description	Qty	Catalog No.
Consoles		
Arthroscopic Energy System	1	AES-1
Probes		
Arthroscopic Energy 90° Probe with Suction	1	AES-90SN
Arthroscopic Energy 50° Probe with Suction	1	AES-50S
Arthroscopic Energy 30° Probe	1	AES-30
Arthroscopic Energy 90° Probe with Suction, Extended Length	1	AES-90SL
Arthroscopic Energy 50° Probe with Suction, Extended Length	1	AES-50SL
Arthroscopic Energy Straight Hook Probe	1	AES-HK00
Arthroscopic Energy 20° Hook Probe	1	AES-HK20

FLUID MANAGEMENT

Description	Qty	Catalog No.
Pumps		
CrystalView™ Pro Fluid Management System	1	CV4000
CrystalView™ Pro Remote	1	CV4040
10K Pump	1	10K
24K Pump	1	24K
Tubing		
EcoFlow™ Day-Use Pump Inflow Tube Set	1	10K600
EcoFlow™ Single-Use Patient Inflow Tube Set	1	10K605
Inflow-Only Tube Set	1	10K100
Inflow-Only Tube Set with Suction Tube	1	10K150
Outflow, Shaver Suction Tube Set	1	24K100
PuddleVac™ Floor Suction Device	1	9321



Ordering Information

LOOKING GLASS®

Description	Qty	Catalog Number
Video System		
Looking Glass® 4K Integrated Visualization System	1	AV4132
Looking Glass® Camera Head – 4K 3MOS EYECUP	1	AV41CHE
Looking Glass® Camera Head – 4K 3MOS EYECUP – Locking	1	AV41CHE-L*
Looking Glass® Cable Connector Interface	1	AV41CON
Touch Glass Standard Image Capture	1	AV41SIC†
Touch Glass Content Management System	1	AV41CMS†
USB Keyboard for Touch Glass	1	DRS29
Light System		
ELG – Electronic Light Guide	1	AV41LG
ELG – Dyonics/Wolf Adapter	1	AVELG51
ELG – Adapter ACMI	1	AVELG52
ELG – CONMED Threaded Adapter	1	AVELG53
Scope to Threaded Adapter	1	AVELG53SA
Direct LED Light Source	1	LS8000
Translator LS8000	1	AV41TLS8
Video Cart		
Video Cart, 120VAC	1	AV8500
Primary Articulating Arm	1	AV8570
Secondary Arm with Keyboard Tray	1	AV8571
Shelf Looking Glass® VP8500 Adapter	1	AV8541
Mounting Bracket Looking Glass® CCI	1	AV8542
Storage Drawer	1	VP8576
Footswitch Bin	1	VP8578
Camera Head Holder	1	VP8577



* Cable Connector Interface

† These items are classified under low voltage directive

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com