



UHD 4mm Arthroscopes with QuickLatch™

CONMED's UHD 4mm QuickLatch™ Arthroscopes and accessories provide increased clarity and detail that allows for more vivid images than our previous HD scopes. The UHD scopes and our fully redesigned QuickLatch™ Sheaths and Bridges are a perfect combination for improved imaging and performance.



LOOKING GLASS™
4K 3MDS EYECUP CAMERA HEAD



LOOKING GLASS™
4K INTEGRATED VISUALIZATION SYSTEM



ARTHROSCOPE
STERILIZATION TRAY



To learn more about these and other innovative products, call **1.866.4CONMED** or visit **CONMED.com**.

ULTRA-HIGH
DEFINITION
QUICKLATCH™
ARTHROSCOPES



ULTRA-HIGH DEFINITION
(UHD) ARTHROSCOPES

Redefining Definition

The better surgeons can see, the better they can treat their patients. CONMED's UHD 4mm QuickLATCH™ Arthroscopes and accessories are designed to set a new standard in vivid clinical visualization.

Ultra-High Definition (UHD) Scopes

CONMED offers higher quality optics throughout the whole imaging chain for a quality image in Orthopedics and ENT.

- **More Consistent Clarity and Brightness**
 - Featuring an increase in the optical quality, improvements in light transmission, and an overall increase in clarity and brightness from the center of the image all the way to the perimeter.
- **Realistic Representation**
 - Improved resolution and depth of field bring attention to critical details and provide a realistic perception of the working cavity.
- **Shorter Working Length**
 - The redesign of the 4mm arthroscope shortens the working length from 170mm to 157mm.
- **Autoclavable**
 - For faster turnaround time, increased efficiency, and lower sterilization costs.
- **Sapphire Lens**
 - Sapphire is used on the distal tip of the lens because it is more scratch resistant than conventional optical glass.

To improve the way you see your patients, schedule an evaluation today by contacting your local CONMED Sales Representative or visit CONMED.com for more information.

QuickLatch™ 4.0mm Sheaths

The QuickLatch™ design increases durability and provides distal tip protection. Tangent position of the scope and proprietary stopcock design provide excellent flow rates in both the 5.5mm and 5.9mm sizes.



QUICKLATCH™
4.0MM SHEATHS

QuickLatch™ 4.0mm Bridge System

The QuickLatch™ Bridge for the 4mm scope is available in both single and double stopcock versions. The QuickLatch™ bridge redesign allows for lavaging, which is key for many arthroscopic cases post-operatively to flush loose debris.



SINGLE STOPCOCK
QUICKLATCH™ BRIDGE

Bridge Cannulas

CONMED offers multiple sizes for multiple requirements. Cannulated obturators are compatible with 3.175mm switching sticks and portal rods.

Inner Diameter	4.6mm	5.8mm	7.6mm
Standard	OC4601	OC5801	OC7601
Standard 30/70°	OC4606	OC5806	OC7606
Trocar	OC4603	OC5803	OC7603
Conical Obturator	OC4604	OC5804	OC7604
Cannulated Obturator	OC4605	OC5805	OC7605

Technology to help you see today, durability to help you see tomorrow

Technology to help you see today, durability to help you see tomorrow. Whether you're looking for a high resolution system with truer colors, features that drive OR efficiency, or reliable performance including a 7-year warranty against defects arising from autoclaving, Looking Glass™ is the solution for your OR.



LOOKING GLASS™

Ordering Information

Description	Catalog Number
UHD Scopes	
Ultra HD 4mm, 0° Scope, 157mm WL	UHD4000
Ultra HD 4mm, 30° Scope, 157mm WL	UHD4300
Ultra HD 4mm, 30° Scope, Reverse Light Post, 157mm WL	UHD4302
Ultra HD 4mm, 70° Scope, 157mm WL	UHD4700
Sheath Systems	
5.5mm, Single Stopcock Arthroscope Sheath System	UQL5531
5.5mm, Single Stopcock Arthroscope Sheath System, Reverse Light Post	UQL5531-R
5.5mm, Double Stopcock Arthroscope Sheath System	UQL5532
5.5mm, Double Stopcock Arthroscope Sheath System, Reverse Light Post	UQL5532-R
5.9mm, Single Stopcock Arthroscope Sheath System	UQL6131
5.9mm, Single Stopcock Arthroscope Sheath System, Reverse Light Post	UQL6131-R
5.9mm, Double Stopcock Arthroscope Sheath System	UQL6132
5.9mm, Double Stopcock Arthroscope Sheath System, Reverse Light Post	UQL6132-R
QuickLatch™ Bridge Systems	
Single Stopcock QuickLatch™ Bridge	UQL6111
Single Stopcock QuickLatch™ Bridge, Reverse Light Post	UQL6111-R
Double Stopcock QuickLatch™ Bridge	UQL6112
Double Stopcock QuickLatch™ Bridge, Reverse Light Post	UQL6112R
4.0mm, Blunt Obturator	UQL6101
4.0mm, Conical Obturator	UQL6102
4.0mm, Sharp Trocar	UQL6103
4.0mm, Cannulated Obturator	UQL6104
28 cm Switching Stick	T4906-11
36 cm Switching Stick	T4906
Portal Rod	T4907
Stopcock Locks	T4914