




# TruPass® SW (Short-Wire) Triple Lumen Papillotome

Compatible with short-wire ERCP techniques, CONMED's TruPass® SW Papillotomes are designed to deliver guaranteed orientation.



**GORE® VIABIL®**  
BILIARY ENDOPROSTHESIS



**ACE™**  
BILIARY GUIDEWIRE



**FLEXXUS®**  
ENDOSCOPIC BILIARY STENT



To learn more about these and other innovative products, call **1-866-4CONMED** or visit **www.CONMED.com**.

**Guaranteed Orientation: TruPass® Papillotomes are guaranteed to orient between 11 and 1 o'clock while approaching the papilla.**

**Smooth, Atraumatic Radiopaque Tip:** The smooth tapered tip reduces tissue trauma to the papilla. Tips are also radiopaque to enable fluoroscopic visualization.

**Color-Coded Cutting Wire Markers:** Easy visibility for determining papillotome depth and cutting wire location. TruPass now available with both the FXWire and ACE straight tip guidewires.

▶ **For a complete ERCP Solution, use TruPass® SW Triple Lumen Papillotomes with FIT™ SW Stone Balloons and FX Wire® Advanced Measurement Guidewire.**

## Ordering Information

Description	Catheter Length	Tip O.D.	Catheter O.D.	Cutting Wire Length	Cutting Wire Type	Guidewire Compatibility	Qty	Catalog Number
<b>NEW</b> TruPass® Triple Lumen Papillotome								
TruPass®	190cm	4.5 F	7-4.5 F	20mm	Monofilament	035" /.89mm	1/box	T2043F
TruPass®	190cm	4.5 F	7-4.5 F	30mm	Monofilament	035" /.89mm	1/box	T3043F
<b>TruPass® SW Triple Lumen Papillotome with Guidewires</b>								
TruPass® with FX Wire®	190cm	4.5 F	7-4.5 F	20mm	Monofilament	035" /.89mm	1/box	P2043F
TruPass® with FX Wire®	190cm	4.5 F	7-4.5 F	30mm	Monofilament	035" /.89mm	1/box	P3043F

Color Coded Guidewire Port, Exit Marker Band, Cutting Wire Markers.

To learn more about this and other innovative products, call **1-866-4CONMED** or visit [www.CONMED.com](http://www.CONMED.com).

Toll Free: 1-866-4CONMED  
International: 1-727-214-3000  
[www.CONMED.com](http://www.CONMED.com)

[customerexperience@conmed.com](mailto:customerexperience@conmed.com)  
[internationalorders@conmed.com](mailto:internationalorders@conmed.com)