



Fixation Products IV Securement

A complete line of quality fixation products to meet your IV securement needs. One company with one goal – to satisfy the customer in the most cost-effective method possible.



STAT 2[®]
CONTROLLER



STAT 2[®] PUMPETTE
CONTROLLER



REUSABLE SpO₂
SENSOR



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com**.



FIXATION
PRODUCTS

Veni-Gard® IV Stabilization Dressing

The ultimate in fixation products, providing stabilization, anchoring, and securement in one product.

Key Features Include:

- Aggressive adhesive foam border system holds objects securely in place and reduces the need for additional securement devices.
- Clear membrane window allows visualization of the site while fixation is maintained.
- Window has a high Moisture Vapor Transmission Rate (MVTR) for breathability.
- Water resistant allowing patients to maintain ADLs.
- Design decreases chance of runaways or catheter pistoning.
- One hand application provides for ease of use.
- Supplied with stretchable sterile tape strips for added stabilization of infusion lines.
- Available in a variety of shapes and sizes for a multitude of fixation applications.
- Sterile.

Veni-Gard® IV Stabilization Dressing

	IV Catheter Securement	Peripheral Line Securement	Jugular Line Securement	Central Line Securement	PICC Line Securement	Femoral Line Securement	Epidural Line Securement	Pediatric Line Securement	Aggressive Securement
705 - 4431		■	■						
730 - 4432								■	
710 - 2730	■			■					
750 - 4181			■	■					■
740 - 1440					■				
745 - 1441			■			■	■		
740 - 1000									■
745 - 1000					■		■		

Product Configurations:



705 - 4431
Veni-Gard TM



730 - 4432
Veni-Gard TM Jr.



710 - 2730
Veni-Gard SP



750 - 4181
Veni-Gard ML

FIXATION
PRODUCTS



740-1440
Veni-Gard AP



745 - 1441
Veni-Gard TPN



740 - 1000
Neoderm



745 - 1000
Neoderm



525 French Road, Utica, NY 13502

FIXATION
PRODUCTS

Ordering Information

Description	Size	Pouch	Box	Case	Catalog Number
Veni-Gard® IV					
Veni-Gard® TM, Traditional omega shape	–	1	100	500	705-4431
Veni-Gard® TM Jr, Pediatric omega shape	–	1	100	500	730-4432
Veni-Gard® SP, Omega shape with IV notch	–	1	100	500	710-2730
Veni-Gard® ML, Elongated omega shape	–	1	25	150	750-4181
Veni-Gard® AP, 2" x 3" Narrow Bordered	50.8mm x 76.2mm	1	100	500	740-1440
Veni-Gard® TPN, 4" x 5" Narrow Bordered	102mm x 127mm	1	25	150	745-1441
Neoderm®					
Neoderm®, 2" x 3" Wide Bordered	50.8mm x 76.2mm	1	100	500	740-1000
Neoderm® Plus, 4" x 5" Wide Bordered	102mm x 127mm	1	100	500	745-1000
Tape					
Strip-N-Tape™, Non-sterile foam tape strip	–	5	1000	5000	701-0510

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.CONMED.com