

# PRESSURE INFUSER IRRIGATION PUMPS

- A full line of single 1-liter and 3-liter systems to meet all of your laparoscopic irrigation pressure source needs.
- The Power Pack Compressor for the Infuser Irrigation Pumps — Super Quiet, Super Convenient!
- All of the tubing and connectors needed to adapt to your OR.



# Pressure Infuser Irrigation Pumps

## The 1-Liter Pressure Infuser Irrigation Pump

- Takes standard irrigation bag, 500 – 1000ml
- IV pole mountable (¾" to 1½")
- Locking pressure adjustment knob
- Comes with 12' of gas source tubing to fit connectors below



## The Dual 1-Liter Pressure Infuser Irrigation Pump

- Takes two bags, 500 – 1000ml — **GREAT for irrigators with dual spikes!**
- IV pole mountable (¾" to 1½")
- Individual locking pressure adjustment knob for each bag chamber
- Comes with 12' of gas source tubing to fit connectors below



## The 3-Liter Pressure Infuser Irrigation Pump

- Most versatile — Takes single bag, 500 – 3000ml
- IV pole mountable (¾" to 1½")
- Locking pressure adjustment knob
- Comes with 12' of gas source tubing to fit connectors below
- Also can take 2 bags, 500 – 1000ml



## The Power Pack — Pressure Infuser Irrigation Pump

- AC powered air compressor for 270178, 270123, and 270165 Infuser Irrigation Pumps
- 110V version available
- IV pole mountable (¾" to 1½")
- UL recognized

## Connectors for Pump Systems



Ohmeda



Chemtron



Puritan-Bennett



Oxequip



Schrader — Non-swivel



Schrader — Swivel

## ORDERING INFORMATION

### Pressure Infuser Irrigation Pump

Ref. No.	Description
270178	The 1-Liter Pressure Infuser Irrigation Pump
270123	The Dual 1-Liter Pressure Infuser Irrigation Pump
270165	The 3-Liter Pressure Infuser Irrigation Pump
270179	The Power Pack — Infuser Pump Compressor

### Gas Source Quick Connect Couplers

Company	Oxygen	Air	Nitrous Oxide
Ohmeda	CD984	CD984A	CD984N
Chemtron	CD983	CD983A	CD983N
Puritan-Bennett	CD989	CD989A	CD989N
Schrader (non-swivel)	CD987	CD987A	CD987N
Schrader (swivel)	CD988	CD988A	CD988N
Oxequip	CD985	CD985A	CD985N

#### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
Utica, NY 13502  
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com