

Proforma™
HF 4.5 CANNULA



Precision Access

FEATURES

- The Only 4.5F Tip, .035" Guidewire Compatible Multi-Lumen Cannula
- Smooth Tapered Tip and Small OD
- Distal Radiopaque Marker Band
- Color Coded Guidewire and Contrast Port
- Ergonomically-designed, Multi-lumen Handle
- Provides Optimal Flow Rate of Contrast Through Dedicated Lumen
- Endoscopic Measurement System
- Available in Straight or Curved Configurations



Performance Guaranteed



43% GREATER DISTAL CATHETER TAPER LENGTH THAN LEADING COMPETITION

Designed for improved deep cannulation performance

4.5F SMOOTH TAPERED TIP AND SMALL OD

Reduces risk of trauma to the papilla and may provide enhanced cannulation performance.

DISTAL RADIOPAQUE MARKER BAND

Enables fluoroscopic visualization of tip location.

COLOR CODED GUIDEWIRE AND CONTRAST PORT

Visual indicator eliminates confusion.

ENHANCED COLUMN STRENGTH

Unique shaft design eliminates the need for a stylet.

Ordering Information

Sterile and disposable

NEW

ProForma™ HF 4.5 Cannulas

Catheter OD 7-5F/1.7mm, Catheter Length 190 cm

Catalog Number	Tip Type	Tip O.D.	Guidewire Compatibility
050907	Curved	4.5F/1.5mm	.035"/.89mm
050906	Straight	4.5F/1.5mm	.035"/.89mm

Pre-curved ProForma™ Cannulas

Catheter OD 5F/1.7mm, Catheter Length 200 cm

Catalog Number	Tip Type	Tip O.D.	Guidewire Compatibility
050310	Standard	5.0F/1.7 mm	.035"/.89mm
050320	Long Taper	4.5F/1.5 mm	.035"/.89mm
050330	Short Taper	4.5F/1.5 mm	.035"/.89mm
050370	Metal Ball	5.5F/1.8 mm	.035"/.89mm

Straight ProForma™ Cannulas

Catheter OD 5F/1.7 mm, Catheter Length 200 cm


Catalog Number	Tip Type	Tip O.D.	Guidewire Compatibility
050300	Standard	5.0F/1.7 mm	.035"/.89mm
050315	Long Taper	4.5F/1.5 mm	.035"/.89mm
050325	Short Taper	4.5F/1.5 mm	.035"/.89mm

Ultra Tapered Tip Cannulas


Catheter OD 5.5F/1.8 mm, Catheter Length 200 cm, *Non-Sterile

Catalog Number	Tip Type	Tip O.D.	Guidewire Compatibility
000156	Tapered, Ultra*	4.0F/1.3 mm	.025"/.89mm
000183	Tapered, 5-4-3*	3.5F/1.2 mm	.018"/.89mm

Superior Pushability



Enhanced Cannulation Design



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



129 Concord Road, Building #3
Billerica, MA 01821-7031
(800) 225-1332

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com