

FXWire™

Advanced Measurement Guidewire

Applied Ingenuity™

CONMED™
ENDOSCOPIC TECHNOLOGIES

An Everyday Guidewire with YOU in Mind

- Hydrophilic, Flexible (Floppy), Atraumatic Tip
 - Enhanced stricture access
 - Ease of cannulation
- Proprietary REGILIAN™ Nitinol Core
 - Optimal column strength for enhanced stricture access
 - Kink resistant
 - Ease of device exchange
 - Superior torque response
- Lubricious Silicone Coating
 - Smooth device exchange
- Light Blue PTFE Jacket with Black Markings
 - Reduced endoscopic glare
 - Enhanced visualization
 - Meets ANSI/AAMI HF-18 standard for cutting over the guidewire



The Only Fully Integrated Measurement Guidewire Designed to Meet ALL of Your Needs




An Advanced Measurement Guidewire Combining Accuracy & Simplicity

- Radiopaque 5 cm Distal Tip with 2 Radiopaque Marker Bands at 10 cm and 15 cm
 - Facilitates accurate fluoroscopic stricture measurement
- Endoscopic Measurement System
 - Individual 1 cm spaced markings
 - Multiple markings every 5 cm for ease of measurement
 - Helps detect inadvertent movement
- Proprietary Proximal Marker Measurement System
 - 1 cm spaced marker bands provide improved stricture length measurement without changing current habits (naked eye)



FXWire™ FEATURES



ORDERING INFORMATION

Sterile, Disposable

Item Number	Wire Type	Guidewire O.D.		Length cm	Tip Type	Torque Vise	Unit
		inches	millimeters				
FX35450S	Standard wire	.035"	.89 mm	450cm	Straight tip	No	3/box
FX35450A	Standard wire	.035"	.89 mm	450cm	Torquable angled tip	Yes	3/box
FX35260S	Standard wire	.035"	.89 mm	260cm	Straight tip	No	3/box
FX35260A	Standard wire	.035"	.89 mm	260cm	Torquable angled tip	Yes	3/box

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
 Utica, NY 13502
 (315) 797-8375
 Customer Service: 1-800-448-6506
 FAX: 1-800-438-3051
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
 email: info@mail.conmed.com
 www.conmed.com