



NAKAO



SPIDER-NET™

RETRIEVAL DEVICE



Inspired by Nature
to Capture & Retrieve

Applied Ingenuity™

 **CONMED™**
ENDOSCOPIC TECHNOLOGIES

Nakao Spider-Net™ vs. Leading Competitor

LOSS PREVENTION DESIGN


- Pouch closes around specimen during retraction reducing risk of loss

TRANSLUCENT NET MATERIAL

- When wet, net material becomes translucent for enhanced endoscopic visibility

SPECIMEN FRIENDLY DESIGN

- Atraumatic net configuration minimizes potential for minced or crushed artifact




GREATER VISIBILITY WHEN WET VS. LEADING COMPETITIVE NET

NAKAO SPIDER-NET™ CREATES A POUCH FOR OPTIMAL SPECIMEN CAPTURE AND RETENTION



SPECIMEN PROTECTION COLLAR REDUCES THE RISK OF SPECIMEN LOSS OR DAMAGE.



NAKAO SPIDER-NET™ RETRIEVAL DEVICES

Non-Sterile, Disposable

Catalog Number	Description	Net Size	Approximate Length	Unit
0230AM	Standard Nakao Spider-Net™ Multi-Pack	2.5cm x 6cm	230cm	20/Box
00230A	Standard Nakao Spider-Net™	2.5cm x 6cm	230cm	10/Box
0160MG	Standard Nakao Spider-Net™ - EGD	2.5cm x 6cm	160cm	5/Box

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375
Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com