

flexxus[®]
biliary stent



Superior Design, Visible Performance


Applied Ingenuity™

CONMED[™]
ENDOSCOPIC TECHNOLOGIES

flexxus[®]


biliary stent

*Superior Design,
Visible Performance*



Unique single handed pistol grip delivery

- One-click control provides unsurpassed placement control and accuracy
- Reduces the potential for premature stent deployment



*For the palliative treatment
of patients with malignant
biliary strictures*

Features






- Highest expansion force available to help maintain lumen patency
- Offering the best fluoroscopic visibility, even under difficult conditions
- Minimal shortening helps ensure accurate placement
- Highly flexible Nitinol stent design conforms even to tortuous anatomy
- Flared stent ends provide excellent stent anchoring within the duct, significantly reducing the potential for migration.
- "One-click single handed pistol grip delivery easy-to-use"

*Data on File

Stricture prior to Flexxus®
stent placement


Flexxus® stent immediately
after deployment

Post deployment



flexxus[®]

biliary stent



Ordering Information

Sterile, Single Patient Use, **Recommended** Guidewire: .035" (FXWire – FX35450S)

Catalogue Number	Stent Diameter x Length	Catheter Diameter x Length	Unit
LXB08040	8mm x 40mm	7.5 F x 190cm	1/Box
LXB08060	8mm x 60mm	7.5 F x 190cm	1/Box
LXB08080	8mm x 80mm	7.5 F x 190cm	1/Box
LXB08100	8mm x 100mm	7.5 F x 190cm	1/Box
LXB10040	10mm x 40mm	7.5 F x 190cm	1/Box
LXB10060	10mm x 60mm	7.5 F x 190cm	1/Box
LXB10080	10mm x 80mm	7.5 F x 190cm	1/Box
LXB10100	10mm x 100mm	7.5 F x 190cm	1/Box

*Warning: The safety and effectiveness of this device for use in the vascular system has not been established.
MRI compatible*

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY

Distributed by:



525 French Road
 Utica, NY 13502
 (315) 797-8375
 Customer Service: 1-800-448-6506
 FAX: 1-800-438-3051
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
 email: info@mail.conmed.com
 www.conmed.com

Manufacturer:
 Angiomend GmbH & Co.
 Medizintechnik KG
 Subsidiary of C.R. Bard, Inc.
 Wachhausstraße 6
 76227 Karlsruhe
 Germany
 Tel: ++49 721 9445 - 9
 Fax: ++49 721 9445 - 111