

# Vizeon™ EUS Fine Needle Aspiration System

**NEW**

25 gauge  
needle design



Designed to provide optimal visualization  
and enhance procedural yield

Applied Ingenuity™

**CONMED™**  
ENDOSCOPIC TECHNOLOGIES

# Vizeon™ EUS Fine Needle Aspiration System

## Enhancing Visibility

### PROPRIETARY VIZEOLYTE TREATED TIP

Provides optimal needle visualization

### ADJUSTABLE SHEATH LENGTH

Ensures scope compatibility

### NITINOL STYLET

Resists kinking and enhances ease of reinsertion

### STEEL LUER-LOCKING HUB

Maintains firm scope connection throughout the procedure

### LOCKING ASPIRATION SYRINGE AND VALVE

Enhances procedural efficiency and yield

### SPRING COIL OR TEFLON SHEATH MATERIAL

Enhances distal tip flexibility

### LOCKABLE STOP ON NEEDLE PISTON

Ensures controlled needle penetration up to 8.5cm

NUMBERED REFERENCE MARKERS – PROVIDE EASY POINT OF REFERENCE

ADJUSTABLE SHEATH PISTON – ENSURES SCOPE COMPATIBILITY

NEEDLE TIP TREATMENT – PROVIDES OPTIMAL NEEDLE VISUALIZATION

RESISTS KINKING AND ENHANCES EASE OF REINSERTION

NITINOL STYLET

MAINTAINS FIRM SCOPE CONNECTION THROUGHOUT PROCEDURE

STEEL LUER-LOCKING HUB

## VIZEON™ EUS FINE NEEDLE ASPIRATION SYSTEM

5 Per Box

Catalog Number	Gauge	Stylet	Needle Diameter	Outer Diameter	Sheath Material
US1822R	22	Round	0.7mm / .028"	1.8mm	Stainless Steel
US1822EB	22	Extended Bevel	0.7mm / .028"	1.8mm	Stainless Steel
US2722R	22	Round	0.7mm / .028"	2.7mm	Stainless Steel
US2119R	19	Round	1.0mm / .039"	2.1mm	Stainless Steel
US1825R (NEW)	25	Round	0.5mm/0.20"	1.8mm	Teflon
US2725R (NEW)	25	Round	0.5mm/0.20"	2.7mm	Teflon

The needle is packaged sterile and includes an aspiration syringe and valve.

### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com