

Optimizer[®] Polypectomy Snares



Measured Performance

Applied Ingenuity™

 **CONMED™**
ENDOSCOPIC TECHNOLOGIES

Optimizer[®] Polypectomy Snares

Optimizer[®] Key Feature

GRADATIONS

- Assists in estimating the amount of tissue ensnared in the loop

REDUCED THROW HANDLE

- Shortened distance to operate the snare handle

FIRM OR SOFT WIRE TYPES

- Both types available in four loop sizes

**GRADATIONS
HELP INDICATE
RESECTION SIZE
& PROGRESS**



Sheath O.D.: 2.3 mm; Sheath Length: 230 cm

Latex^{FREE}

<i>Sterile, Disposable</i>		Loop		Active		Sheath		
Catalog Number	Size		Shape	Cord	O.D.	Length		Unit
Firm Loop Wire								
000472	Large		Oval	Olympus	2.3mm	230cm		10/box
000473	Large		Oval	Microvasive	2.3mm	230cm		10/box
000474	Medium		Oval	Olympus	2.3mm	230cm		10/box
000475	Medium		Oval	Microvasive	2.3mm	230cm		10/box
000476	Small (Mini)		Oval	Olympus	2.3mm	230cm		10/box
000477	Small (Mini)		Oval	Microvasive	2.3mm	230cm		10/box
000478	X-Small		Oval	Olympus	2.3mm	230cm		10/box
000479	X-Small		Oval	Microvasive	2.3mm	230cm		10/box
Soft Loop Wire								
000460	Large		Oval	Olympus	2.3mm	230cm		10/box
000461	Large		Oval	Microvasive	2.3mm	230cm		10/box
000462	Medium		Oval	Olympus	2.3mm	230cm		10/box
000463	Medium		Oval	Microvasive	2.3mm	230cm		10/box
000464	Small (Mini)		Oval	Olympus	2.3mm	230cm		10/box
000465	Small (Mini)		Oval	Microvasive	2.3mm	230cm		10/box
000466	X-Small		Oval	Olympus	2.3mm	230cm		10/box
000467	X-Small		Oval	Microvasive	2.3mm	230cm		10/box

Boston Scientific Endoscopy is a trademark of Boston Scientific Corporation. Olympus is a registered trademark of Olympus America, Inc.

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



Subsidiary of ConMed Corporation

525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

© 2007 CONMED Corporation, 4/2007, Control #12631 Rev A

MBR9126