

Singular[®] Polypectomy Snares



Short Throw Performance
For Nurse Comfort

Applied Ingenuity[™] |

 **CONMED[™]**
ENDOSCOPIC TECHNOLOGIES

Singular[®] Polypectomy Snares

Singular[™] Key Feature

**PATENTED TRUE 2 TO 1
PULLY SYSTEM DELIVERS
UNIQUE SHORT THROW
PERFORMANCE**

SHORT THROW HANDLE

- Facilitates one-handed operation and user comfort
- Provides enhanced tactile feel

FULL RANGE OF LOOP SHAPES AND SIZES

- Four Oval Sizes
- Medium Crescent
- Medium Hexagonal

FIRM OR SOFT WIRE TYPES



Sheath O.D.: 2.3 mm; Sheath Length: 230 cm

Sterile, Disposable

Latex^{FREE}

Catalog Number	Size	Loop	Shape	Active Cord	O.D.	Sheath Length	Unit
Firm Loop Wire							
000957	Medium		Crescent	Olympus	2.3mm	230cm	10/box
000958	Medium		Crescent	Microvasive	2.3mm	230cm	10/box
000959	Medium		Hexagonal	Olympus	2.3mm	230cm	10/box
000960	Medium		Hexagonal	Microvasive	2.3mm	230cm	10/box
000963	Large		Oval	Olympus	2.3mm	230cm	10/box
000964	Large		Oval	Microvasive	2.3mm	230cm	10/box
000965	Medium		Oval	Olympus	2.3mm	230cm	10/box
000966	Medium		Oval	Microvasive	2.3mm	230cm	10/box
000978	Small (Mini)		Oval	Olympus	2.3mm	230cm	10/box
000979	Small (Mini)		Oval	Microvasive	2.3mm	230cm	10/box
000991	X-Small		Oval	Olympus	2.3mm	230cm	10/box
000992	X-Small		Oval	Microvasive	2.3mm	230cm	10/box
Soft Loop Wire							
000977	Large		Oval	Olympus	2.3mm	230cm	10/box
000980	Large		Oval	Microvasive	2.3mm	230cm	10/box
000983	Medium		Oval	Olympus	2.3mm	230cm	10/box
000984	Medium		Oval	Microvasive	2.3mm	230cm	10/box
000985	Small (Mini)		Oval	Olympus	2.3mm	230cm	10/box
000986	Small (Mini)		Oval	Microvasive	2.3mm	230cm	10/box
000987	X-Small		Oval	Olympus	2.3mm	230cm	10/box
000988	X-Small		Oval	Microvasive	2.3mm	230cm	10/box

Boston Scientific Endoscopy is a trademark of Boston Scientific Corporation. Olympus is a registered trademark of Olympus America, Inc.

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



Subsidiary of ConMed Corporation

525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

© 2007 CONMED Corporation, 4/2007, Control #12613 Rev A

MBR9127