

Precisor[®] EXL[™] Disposable Biopsy Forceps



Superior CAM Design
Technology Designed to Cut Tissue Not Tear

Applied Ingenuity[™] |

 **CONMED[™]**
ENDOSCOPIC TECHNOLOGIES

Precisor® EXL™ Disposable Biopsy Forceps

■ Superior CAM Design Technology

- Designed to cut tissue; not tear larger jaw sample capacity.
- Greatest cutting force as jaws close
- Helps maintain channel integrity

■ Jaw Sample Capacity

- Precisor® jaws have a greater capacity, and more defined edge

■ Jaw Length

- Precisor® jaws have a shorter cup length for enhanced pushability

■ Jaw Opening Span

- Precisor® jaws have a wide jaw span

■ Three Ring Handle Design

- Precisor® cold and hot forceps handles provide greater control and tactile feel

Boston Scientific Endoscopy is a trademark of Boston Scientific Corporation.

Olympus is a registered trademark of Olympus America, Inc.



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506
FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

Precisor® EXL™ Coated Disposable Biopsy Forceps

Featuring Channel-Guard™ Low Friction Sheath

Latex **FREE**

Sterile, Coated

Catalog Number	Configuration	O.D.	Length	Color	Unit
Standard GI					
000381	Alligator Cup	2.3mm	160cm	Yellow	10/box
000382	Alligator Cup with Needle	2.3mm	160cm	Yellow	10/box
000383	Oval Cup	2.3mm	160cm	Yellow	10/box
000384	Oval Cup with Needle	2.3mm	160cm	Yellow	10/box
000385	Alligator Cup	2.3mm	230cm	Blue	10/box
000386	Alligator Cup with Needle	2.3mm	230cm	Blue	10/box
000389	Oval Cup	2.3mm	230cm	Blue	10/box
000388	Oval Cup with Needle	2.3mm	230cm	Blue	10/box

Jumbo GI

00820C	Oval Cup	3.1mm	230cm	Blue	10/box
00821C	Oval Cup with Needle	3.1mm	230cm	Blue	10/box
00822C	Alligator Cup	3.1mm	230cm	Blue	10/box
00823C	Alligator Cup with Needle	3.1mm	230cm	Blue	10/box

Pediatric GI

Sterile, Coated

Catalog Number	Configuration	O.D.	Length	Color	Unit
Standard GI					
100512	Oval Cup	1.8mm	160cm	Yellow	10/box
100513	Alligator Cup	1.8mm	160cm	Yellow	10/box

Precisor Broncho® Disposable Biopsy Forceps

Sterile, Disposable

Catalog Number	Description	Length	O.D.	Unit
Standard GI				
100502	Oval Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100503	Alligator Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100506	Oval Cup Forceps	115cm	1.8mm	10/box
100507	Alligator Cup Forceps	115cm	1.8mm	10/box

Precisor® Disposable Biopsy Forceps

Sterile, Uncoated

Catalog Number	Configuration	O.D.	Length	Unit
Standard GI				
000837	Oval Cup	2.3mm	160cm	10/box
000838	Oval Cup with Needle	2.3mm	160cm	10/box
000840	Alligator Cup	2.3mm	160cm	10/box
000841	Alligator Cup with Needle	2.3mm	160cm	10/box
000816	Oval Cup	2.3mm	230cm	10/box
000817	Oval Cup with Needle	2.3mm	230cm	10/box
000818	Alligator Cup	2.3mm	230cm	10/box
000819	Alligator Cup with Needle	2.3mm	230cm	10/box
Jumbo GI				
000820	Oval Cup	3.1mm	230cm	10/box
000821	Oval Cup with Needle	3.1mm	230cm	10/box
000822	Alligator Cup	3.1mm	230cm	10/box
000823	Alligator Cup with Needle	3.1mm	230cm	10/box

Pediatric GI

100516	Oval Cup	1.8mm	160cm	10/box
100517	Alligator Cup	1.8mm	160cm	10/box

Precisor® Hot Disposable Biopsy Forceps

Sterile, Coated

Catalog Number	Configuration	O.D.	Length	Coating	Cord	Unit
Standard GI						
000852	Oval Cup	2.3mm	230cm	Red	Olympus®	10/box
000853	Oval Cup	2.3mm	230cm	Red	Microvasive®	10/box
000854	Alligator Cup	2.3mm	230cm	Red	Olympus®	10/box
000856	Alligator Cup	2.3mm	230cm	Red	Microvasive®	10/box