



FLEXITIP™

Disposable Sclerotherapy Needle



Flexible Tip, Visible Performance

Applied Ingenuity™

CONMED™
ENDOSCOPIC TECHNOLOGIES

FLEXITIP™

Disposable Sclerotherapy Needle

FlexiTip™ Key Feature

OPTIC YELLOW
TIP FOR
VISIBILITY



OPTIC YELLOW TIP

- Highly visible in bloody field

ULTRA FLEXIBLE TIP

- Easily passed over deflected elevator of duodenoscope

RETRACTABLE, LOCKABLE NEEDLE

- Minimizes needle bounce
- Optimizes control during procedure

ERGONOMIC HANDLE DESIGN

- Provides control
- Convenient introduction of injection media

FLEXITIP™ DISPOSABLE SCLEROTHERAPY NEEDLES

Latex **FREE**

Sterile, Disposable

Catalog Number	Description	Needle		Sheath			Unit
		Length	Gauge	Length	O.D.		
000215	Optic Yellow Tip Needle	4mm	25 g	160cm	2.3mm	5/box	
000216	Optic Yellow Tip Needle	5mm	25 g	160cm	2.3mm	5/box	
000217	Optic Yellow Tip Needle	6mm	25 g	160cm	2.3mm	5/box	
000218	Optic Yellow Tip Colonoscope Needle	5mm	25 g	230cm	2.3mm	5/box	
000198	Standard Needle	4mm	25 g	160cm	2.3mm	5/box	
000199	Standard Needle	5mm	25 g	160cm	2.3mm	5/box	
000200	Standard Needle	6mm	25 g	160cm	2.3mm	5/box	
000138	Standard Colonoscope Needle	5mm	25 g	230cm	2.3mm	5/box	

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



Subsidiary of ConMed Corporation

525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

© 2007 CONMED Corporation, 5/2008, Control #17255 Rev A

MBR 9133 Rev 1