

# Sure Shot<sup>®</sup>



**Metal Sheathed Injection Needle**

Applied Ingenuity™

**CONMED™**  
ENDOSCOPIC TECHNOLOGIES



## Sure Shot® Key Feature

### STAINLESS STEEL SHEATH

- Eliminates Kinking
- Puncture Resistant
- Enhances Pushability
- Maximizes Needle Actuation While Retroflexed


### RETRACTABLE, LOCKABLE NEEDLE

- Minimizes Needle Bounce
- Optimizes Control During Procedure

### ERGONOMIC HANDLE DESIGN

- Provides Control
- Convenient Introduction of Injection Media

### METAL SHEATH MAXIMIZES NEEDLE ACTUATION



## SURE SHOT® METAL SHEATHED INJECTION NEEDLE

Latex **FREE**

*Sterile, Disposable*

Catalog Number	Description	Needle		Sheath			Unit
		Length	Gauge	Length	O.D.		
100218	Metal Sheathed Injection Needle	5mm	25 g	230cm	2.5mm	5/box	

#### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



Subsidiary of ConMed Corporation

525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com  
 © 2007 CONMED Corporation, 9/2007, Control #12741 Rev A