


# Apollo<sup>®</sup>

Papillotomes

Featuring

Advanced Cannulation Triple Lumen  
**Apollo<sup>3</sup> AC<sup>™</sup>**  
Papillotome

Advanced Cannulation  
Triple Lumen Papillotome



Get Oriented

Applied Ingenuity<sup>™</sup>

 **CONMED<sup>™</sup>**  
ENDOSCOPIC TECHNOLOGIES



Advanced Cannulation Triple Lumen

**Apollo**<sup>3 AC</sup>™  
Papillotome


# Advanced Cannulation Triple Lumen Papillotome

**4.5F Tip** provides .035" compatible triple lumen papillotome.

- Smooth taper reduces trauma to the papilla
- 4.5F tip designed to improve cannulation performance



## Papillotome Tip



- CONMED®
- Existing Technologies



# Advanced Cannulation Double Lumen Papillotome




## Guaranteed Orientation

The Apollo® papillotomes are guaranteed to orient between 11 and 1 o'clock.

## Other Advantages

- **Radiopaque distal tip** enables fluoroscopic visualization of tip location.
- **Color coded cutting wire markers** provide endoscopic visibility for determining papillotome depth and cutting wire location.
- **Exit marker on the shaft of the papillotome** provides a visual cue to the clinician as to when the papillotome is exiting the scope.
- **Color coded guidewire port.** Visual indicator eliminates confusion in a dark procedure room.
- **Full range of tip shapes** accommodate a variety of procedure needs.



## Ordering Information

### Tapered Tip Apollo® Triple Lumen Papillotomes

Catalog Number	Tip Type	Tip Length	Tip OD	Cutting Wire Length	Cutting Wire Type	Catheter Length	Catheter OD	Guidewire Compatibility
007104	Short Nose	8 mm	5F	20 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
007105	Short Nose	8 mm	5F	30 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
7104AC	Short Nose	5 mm	4.5F	20 mm	Monofilament	190 cm	7-4.5F	.035"/.89 mm
7105AC	Short Nose	5 mm	4.5F	30 mm	Monofilament	190 cm	7-4.5F	.035"/.89 mm
007106	Long Nose	20 mm	5F	20 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
007107	Long Nose	20 mm	5F	30 mm	Monofilament	190 cm	7-5F	.035"/.89 mm

### Tapered Tip Apollo® Double Lumen Papillotomes

Catalog Number	Tip Type	Tip Length	Tip OD	Cutting Wire Length	Cutting Wire Type	Catheter Length	Catheter OD	Guidewire Compatibility
007100	Short Nose	8 mm	5F	20 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
007101	Short Nose	8 mm	5F	30 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
7100AC	Short Nose	5 mm	4.5F	20 mm	Monofilament	190 cm	7-4.5F	.035"/.89 mm
7101AC	Short Nose	5 mm	4.5F	30 mm	Monofilament	190 cm	7-4.5F	.035"/.89 mm
007102	Long Nose	20 mm	5F	20 mm	Monofilament	190 cm	7-5F	.035"/.89 mm
007103	Long Nose	20 mm	5F	30 mm	Monofilament	190 cm	7-5F	.035"/.89 mm

### Beveled Tip Apollo® Triple Lumen Papillotomes

Catalog Number	Tip Type	Tip Length	Tip OD	Cutting Wire Length	Cutting Wire Type	Catheter Length	Catheter OD	Guidewire Compatibility
050071	Short Nose	8 mm	5.5F	20 mm	Monofilament	190 cm	7-5.5F	.035"/.89 mm
050073	Short Nose	8 mm	5.5F	30 mm	Monofilament	190 cm	7-5.5F	.035"/.89 mm

### Beveled Tip Apollo® Double Lumen Papillotomes

Catalog Number	Tip Type	Tip Length	Tip OD	Cutting Wire Length	Cutting Wire Type	Catheter Length	Catheter OD	Guidewire Compatibility
050063	Short Nose	8 mm	5.5F	20mm	Monofilament	190cm	7-5.5	.035"/.89 mm
050064	Short Nose	8 mm	5.5F	30mm	Monofilament	190cm	7-5.5	.035"/.89 mm
050053	Short Nose	8 mm	5.5F	20mm	Braided	190cm	7-5.5	.035"/.89 mm
050054	Short Nose	8 mm	5.5F	30mm	Braided	190cm	7-5.5	.035"/.89 mm

All Apollo papillotomes accommodate an Olympus® active cord and are packaged 1/box.

#### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
Utica, NY 13502  
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com