


ELIMINATOR[®]

Biliary Balloon Dilators




Stronger than the stricture.

Applied Ingenuity™

CONMED™
ENDOSCOPIC TECHNOLOGIES

ELIMINATOR[®]

Biliary Balloon Dilators



■ Patented High Pressure PET Design


- Does not waist or dogbone.

■ Radiopaque Marker Bands

- Provide fluoroscopic identification of balloon working length.

■ .035 Guidewire Compatibility

Performance Guaranteed!



Ordering Information

Eliminator Biliary Balloon Dilators

Catalog Number	Balloon Size	Catheter OD	Length	Guidewire Compatibility
000347	12F/4 mm x 2 cm	5F	180 cm	.035" / .89 mm
000348	18F/6 mm x 2 cm	5F	180 cm	.035" / .89 mm
000349	24F/8 mm x 3 cm	5F	180 cm	.035" / .89 mm

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com