

NOT YOUR ORDINARY BLADE...

- EFFICIENT
- EFFECTIVE
- VALUE-ADDED




EFFICIENT, EFFECTIVE, VALUE-ADDED

ConMed Linvatec offers a unique line of dual-function blades that can effectively resect both soft tissue and bone. These high-performance resection devices will save you time and money by minimizing the need to switch between a blade and a bur.

Ultracut®


High performance cutting surface on the inner tube constructed using ultra-hard Trimrite steel, allows this dual-function blade to be used for the effective resection of soft tissue and bone. The large cutting window helps reduce clogging.



3.5mm Ultracut.....	35ULT-ST-ZZ
4.2mm Ultracut.....	C9405A
4.2mm 15° Pre-Bent Ultracut	42ULT-ZZ-MT
4.2mm Extended Length Ultracut.....	EL9405
5.5mm Ultracut.....	C9415

Mako™



The double window tooth pattern of the split inner tube combined with its resilient cutting edges enables the Mako blade to maximize bone and tissue resection safely and effectively. The cutting edges spring inward under heavy loads, minimizing debris generation.



4.2mm Mako	9298A
5.5mm Mako	C9298A

Lightning®

With new enhanced design features, the Lightning Bur is more powerful and effective than ever! Its unique fluted inner tube further breaks down the bone and tissue as it passes through the bur directly to the suction port, reducing clogging. The Lightning bur is ideal for ultra-aggressive resection of bone and tissue.



4.5mm Lightning	H9133
-----------------------	-------

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com