



## High Definition Scopes and QuickLatch®

- 4.0mm HD Scope Systems
- Quicklatch Sheaths (4.0mm scopes)
- True autoclavability
- Eyepiece and Cartridge available



*CLEARLY TRUE HIGH DEFINITION™*

True **HD**

**CONMED**  
LINVATEC


# High Definition Scopes

ConMed Linvatec brings another first to market with the HD arthroscope. Continuing our drive of optical optimization, the enhanced lens design ensures maximized transmission of high contrast light from the arthroscope into the True HD camera head. ConMed Linvatec can now offer the highest quality optics throughout the whole Imaging chain for a True HD image in ENT and Orthopedics. Our HD scope range is compatible with your existing QuickLatch® sheath range ensuring you don't have to settle for low flow to operate in high definition.

- Improved light transmission
- High contrast to the picture edge
- Improved depth of field
- High contrast to the back of the joint

	EYECUP	CARTRIDGE
0 Degree	HD4000	HD4001
30 Degree	HD4300	HD4301
70 Degree	HD4700	HD4701
30 Degree Rev	HD4302*	

\*For use with QLR6131 or QLR6132 sheath sets only.



# Sheath Systems

## QuickLatch® 4.0mm Sheaths

QuickLatch design increases durability and provides excellent distal tip protection. Tangent position of scope and proprietary stopcock design provide excellent flow rates in both the 5.5mm and 5.9mm sizes

	<b>SINGLE</b>	<b>DUAL</b>
5.5mm Arthroscopy	QL5531	QL5532
5.9mm High Flow	QL6131	QL6132
5.9mm Two-Connect	QL6231	QL6232

See Back Panel for Trocar and Accessory Items



## Bridge Cannulas



Multiple sizes for multiple requirements. Cannulated obturators are compatible with 3.175mm switching sticks and portal rods.

	<b>4.6mm</b>	<b>5.8mm</b>	<b>7.6mm</b>
	Inner Diameter		
Standard	OC4601	OC5801	OC7601
Standard 30/70°	OC4606	OC5806	OC7606
Trocar	OC4603	OC5803	OC7603
Conical Obturator	OC4604	OC5804	OC7604
Cannulated Obturator	OC4605	OC5805	OC7605

## QuickLatch 4.0mm Bridge System

The QuickLatch Bridge for 4mm scopes is available in both single and double stopcock versions. Two-Connect adapter allows for "bridging" with ConMed Linvatec's pressure sensing pumps (4.6mm Cannulas are not compatible).



Single Stopcock	QL6111
Dual Stopcock	QL6112
Two-Connect	QL6210*

\*With QL6230, T5101 or T5102

# High Definition 4.0mm Scopes



CLEARLY TRUE HIGH DEFINITION™  
True HD

## Scopes

Outside Diameter  
Field of View  
Angle of View  
Working Length

## High Definition 4.0mm Scope

4mm  
105°  
0°, 30°, 70°  
0°, 172mm  
30°, 170mm  
N/A  
70°, 173mm

## Sheaths

Outside Diameter  
  
Working Length

## QuickLatch Sheaths (4.0mm Scopes)

5.5mm Arthroscopy Sheath, 5.5mm  
5.9mm High Flow Sheath, 5.9mm  
5.9mm Two-Connect Sheath, 5.9mm  
5.5mm Arthroscopy Sheath, 142mm  
5.9mm High Flow Sheath, 142mm  
5.9mm Two-Connect Sheath, 125mm



## ACCESSORIES

### Sheaths Description

Blunt Obturator  
Conical Obturator  
Sharp Trocar  
Cannulated Obturator  
28cm Switching Stick  
36cm Switching Stick  
Portal Rod  
Stopcock Locks

### QuickLatch Sheaths (4.0mm Scopes)

QL6101  
QL6102  
QL6103  
QL6104  
T4906-11  
T4906  
T4907  
T4914



**HIGH DEFINITION Monitors**



**IM4000 High Definition Camera**



### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard  
Largo, FL 33773-4908  
(727) 392-6464  
Customer Service: 1-800-237-0169  
FAX: (727) 399-5256  
International FAX: +1 (727) 397-4540  
email: customer\_service@linvatec.com  
www.linvatec.com