

HIGH DEFINITION Monitor



ConMed Linvatec is proud to offer the VP4726/F High Definition Monitor. This 26 inch 1080p medical grade monitor incorporates a 16:9 aspect ratio and TrueHD picture quality. Available with or without fiber optic input. You can now maximize the benefits of High Definition technology at the cutting edge of surgical HDTV.


- Fully HD and Digital Compatible
- Unrivalled 1920 x 1200 Resolution
- Sealed Front Panel in Light-weight Design with VESA Mount
- Wide-screen 16:9 aspect ratio
- PIP functionality



HD 1080p



HIGH DEFINITION Monitor



HD 1080p

VP4726/F Specifications:

LCD Panel:	Type TFT Active Matrix LCD with an AR-coated protection panel
Resolution (W x H)	1920 x 1200 (WUXGA)
Picture Size (W x H)	Approximately 21.7 x 13.6 inches (551 x 345mm)
Picture Size (Diagonal)	26 inches (660mm)
Aspect Ratio	16:9
Contrast Ratio	1000:1
Luminance	400 cd/m ²
Viewing Angle	89° /89° /89° /89°
Dimensions (W x H x D)	
Inches	24.5 x 16.6 x 4.0
mm	622 x 422 x 102
Weight	18.2 lbs. (8.5 Kg)



The VP4726/F enable the wide-screen viewing of High Definition images using HD-SDI, DVI, or Fiber Optic (VP4726F only) video inputs.



Catalog Number	Description
VP4726	HD LCD Monitor
VP4726F	HD LCD Monitor, with Fiber Optic Input

Input/Output

Inputs	
SDI, HD-SDI	BNC (x1)
S-video, Y/C	BNC (x2)
R/G/B/S, YPbPr	BNC (x5)
Composite	BNC (x1)
DVI	DVI-d
VGA	HD-15
Fiber Optic	4 pin (VP4726F only)
Control Inputs	
Serial Remote	D-sub 9pin, RS-232C
GPIO	RJH-4
Outputs	
SDI, HD-SDI	BNC (x1)
DVI	DVI-d

General

Power Consumption	
Maximum	120W
Stand-by Mode	20W
Power Requirement	
Display	24 VDC, 6.25A Max
DC Power Supply	100 to 240 VAC, 50/60Hz
Temperature	
Operating	0 to +40° C (+32 to +104° F)
Storage	-20 to +60° C (-4 to 140° F)
Humidity	
Operating	30 to 75% non-condensing
Storage	10 to 85% non-condensing
Mounting	100mm VESA

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com