


NEW

LINVATEC
SRS

SHOULDER RESTORATION SYSTEM

Simple | Secure | Versatile



CONMED
LINVATEC

Simple solutions for complex procedures

- **Simple**
- **Secure:** Maximize footprint compression for enhanced healing
- **Versatile:** Intra-operative choice

Primary or Medial Row Fixation



SuperRevo®-FT and ThRevo®-FT

Titanium fully threaded anchor
5.0mm Diameter

- **Proven Revo® technology**
- **Double or triple loaded with suture**
- **Self-drilling**

CrossFT™

All PEEK-OPTIMA® polymer fully threaded anchor

4.5, 5.5, 6.5mm Diameter

- **Cancellous and cortical threads provide exceptional fixation strength**
- **Patent pending broaching punch minimizes cortical stress fractures**
- **Cannulated to channel growth factors**
- **Double or triple loaded with suture**



EM

*for single and double row,
SLAP, and instability repairs*

NEW

Primary or Lateral Row Fixation



PopLok™

All PEEK-OPTIMA® polymer knotless anchor
3.5, 4.5mm Diameter

- Provides ability to tension sutures individually and separately from anchor seating
- Suture locks within the anchor
- Sub-cortical wings provide fixation
- Tactile and audible "POP" confirms fixation

Seat, tension, pop and go!


SLAP and Instability Repairs



Presto™ Suture Anchor

100% Poly Lactic Acid (PLA) anchor
2.9mm Diameter

- Press fit style
- Single or double loaded
- Purple in color for easy visualization
- Pre-loaded on an arthroscopic disposable inserter



Paladin™

High Strength 96L / 4D Self-Reinforced™ Poly Lactic Acid
5.0, 6.5mm Diameter

- Low profile head of only 2.75mm
- Proximal and Distal laser etch indicating eyelet orientation
- Distal etch insertion depth guides

ORDERING INFORMATION

PopLok™ Anchor

CKP-4500	4.5mm PopLok Anchor
CKP-4502	4.5mm PopLok Anchor w/two #2 HiFi® Sutures
CKP-3500	3.5mm PopLok Anchor
CKP-3501	3.5mm PopLok Anchor w/one #2 HiFi® Suture

CrossFT™ Anchor

CFP-4502	4.5mm CrossFT , 2 sutures
CFP-4503	4.5mm CrossFT , 3 sutures
CFP-5502	5.5mm CrossFT , 2 sutures
CFP-5503	5.5mm CrossFT , 3 sutures
CFP-6502	6.5mm CrossFT , 2 sutures
CFP-6503	6.5mm CrossFT , 3 sutures

Paladin™ Rotator Cuff Anchor

C5070H	Paladin 5.0mm Driver with Anchor
C6570H	Paladin 6.5mm Driver with Anchor

Presto™ Anchor

C6271H	Single Loaded Presto™ Anchor
C6272H	Double Loaded Presto™ Anchor
C6179	Instability Twist Drill Bit

Super Revo®- FT & ThRevo®- FT Anchors

CF6140H	5.0mm Super Revo-FT
CF6160H	5.0mm ThRevo-FT

PopLok™ Instrumentation

PKL-35M	3.5mm PopLok Punch
PKL-45M	4.5mm PopLok Punch

CrossFT™ Instruments

PFT-00M	Universal CrossFT Broaching Punch
PFT-45M	4.5mm CrossFT Broaching Punch
PFT-55M	5.5mm CrossFT Broaching Punch
PFT-65M	6.5mm CrossFT Broaching Punch
TFT-45M	4.5mm CrossFT Tap
TFT-55M	5.5mm CrossFT Tap
TFT-65M	6.5mm CrossFT Tap

Bio-Instability Instrumentation for Presto™


C6171	Bio-Instability Drill Guide, Fishmouth
C6172	Bio-Instability Drill Guide, Fishmouth
C6173	Blunt Obturator
C6174	Sharp Trocar
C6176	Bio-instability Bone Punch, 2.1mm
C6178	Bio-instability Sterilization tray
C6179	Bio-instability Twist Drill Bit, 2.1mm

Paladin™ Rotator Cuff Instruments

C5075	5.0mm Paladin Tap
C6575	6.5mm Paladin Tap

Linvatec-SRS Sterilization Tray

RCR-TRAY	Linvatec Shoulder Restoration System Tray
----------	---



Visit SRS.Linvatec.com today to learn more.

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com



CrossFT™ and PopLok™ are comprised of PEEK-OPTIMA® polymer from Invibio® Biomaterial Solutions. PEEK-OPTIMA and INVIBIO are registered trademarks of Invibio Ltd. All Rights are Reserved.



This product contains Dyneema® Purity. Dyneema® Purity is a registered trademark of Royal DSM N.V.