

*One company, one goal — to satisfy the customer
in the most cost-effective method possible.*






Fixation Products

**A Complete Line of Quality Fixation
Products to Meet All of Your
IV Securement Needs**

Veni-Gard® IV Stabilization Dressing

The ultimate in fixation products, providing stabilization, anchoring, and securement in one product.

- Aggressive adhesive foam border system holds objects securely in place and reduces the need for additional securement devices.
- Clear membrane window allows visualization of the site while fixation is maintained.
- Window has a high MVTR (Moisture Vapor Transmission Rate) for breathability.
- Water resistant allowing patients to maintain ADL's.
- Design decreases chance of runaways or catheter pistoning.
- One hand application provides for ease of use.
- Supplied with stretchable sterile tape strips for added stabilization of infusion lines.
- Available in a variety of shapes and sizes for a multitude of fixation applications.



Clear-Gard™ Transparent IV Dressing

A completely transparent fixation device for IV cannulation.

- Completely transparent for full visualization of IV site and surrounding skin.
- Designed for patients with fragile skin.
- High moisture vapor transmission rate.
- Support wings hold IV cannula in place while still allowing total access to hub.



ClearSite® TM Transparent Membrane Dressing

A completely transparent fixation device available in a variety of sizes for varying applications.



- Stabilization carrier prevents dressing from curling and self-adhesion upon patient application.
- Unique delivery system allows for true one-handed application.
- High moisture vapor transmission rate (MVTR) ensures breathability and patient comfort.

| ClearSite® TM | 3M Tegaderm™ | S&N OpSite™ |
|---------------|--------------|-------------|
| 741-1022 | 1622W | 4654 |
| 741-1023 | 1624W | 4007 |
| 741-1034 | 1635 | — |
| 741-1044 | 1630 | — |
| 741-1040 | 1627 | 4542 |
| 741-1045 | 1626W | 4008 |
| 741-1068 | 1628 | 4631 |
| 741-1080 | 1629 | 4987 |

1. CONMED does not manufacture Tegaderm™ or OpSite™ products.
2. Tegaderm™ is a registered trademark of 3M.
3. OpSite™ is a registered trademark of Smith & Nephew.

SureTape™ Sterile Tape Strips

A sterile tape application system designed to prevent cross contamination.

- Patented “third hand” delivery system allows for sterile application by single clinician.
- Can be easily applied with gloves.
- Unique material tears straight each time.
- Prevents cross contamination.
- Available in a variety of sizes.



Fixation Product Line

ORDERING INFORMATION

| Product | Description | Packaging | | |
|----------------------------|-------------------------------|-----------|------|------|
| | | Pouch | Box | Case |
| Veni-Gard® | | | | |
| 701-0510 Strip-N-Tape® | Non-sterile foam tape strip | 5 | 1000 | 5000 |
| 705-4431 Veni-Gard® TM | Traditional omega shape | 1 | 100 | 500 |
| 705-4432 Veni-Gard® TM Jr. | Pediatric omega shape | 1 | 100 | 500 |
| 710-2730 Veni-Gard® SP | Omega shape with IV notch | 1 | 100 | 500 |
| 750-4181 Veni-Gard® ML | Elongated omega shape | 1 | 25 | 150 |
| 740-1440 Veni-Gard® AP | 2" X 3" Narrow bordered | 1 | 100 | 500 |
| 745-1441 Veni-Gard® TPN | 4" X 5" Narrow bordered | 1 | 25 | 150 |
| 740-1000 Neoderm® | 2" X 3" Wide bordered | 1 | 100 | 500 |
| 745-1000 Neoderm® Plus | 4" X 5" Wide bordered | 1 | 100 | 500 |
| Clear-Gard™ | | | | |
| 755-2002 Clear-Gard™ | Winged shape | 1 | 100 | 400 |
| ClearSite® TM | | | | |
| 741-1022 ClearSite® TM | 2" X 2" Transparent membrane | 1 | 100 | 400 |
| 741-1023 ClearSite® TM | 2" X 3" Transparent membrane | 1 | 100 | 400 |
| 741-1034 ClearSite® TM | 3" X 4" Transparent membrane | 1 | 100 | 400 |
| 741-1044 ClearSite® TM | 4" X 4" Transparent membrane | 1 | 100 | 400 |
| 741-1045 ClearSite® TM | 4" X 5" Transparent membrane | 1 | 50 | 200 |
| 741-1040 ClearSite® TM | 4" X 10" Transparent membrane | 1 | 50 | 200 |
| 741-1068 ClearSite® TM | 6" X 8" Transparent membrane | 1 | 25 | 100 |
| 741-1080 ClearSite® TM | 8" X 10" Transparent membrane | 1 | 50 | 200 |
| SureTape™ | | | | |
| 800-0403 SureTape™ | 4" X 3" Sterile tape strips | 4 | 100 | 400 |
| 800-0406 SureTape™ | 4" X 6" Sterile tape strips | 4 | 50 | 200 |
| 800-0410 SureTape™ | 4" X 10" Sterile tape strips | 4 | 50 | 200 |

FIXATION PRODUCTS APPLICATION CHART

| | IV Catheter Securement | Peripheral Line Securement | Jugular Line Securement | Central Line Securement | PICC Line Securement | Femoral Line Securement | Epidural Line Securement | Pediatric Line Securement | Fragile Skin Securement | Post-Op Wound Closure | Aggressive Securement |
|----------------------|---------------------------|-------------------------------|----------------------------|----------------------------|-------------------------|----------------------------|-----------------------------|------------------------------|----------------------------|--------------------------|--------------------------|
| Veni-Gard® | | | | | | | | | | | |
| 705-4431 | | X | X | | | | | | | | |
| 730-4432 | | | | | | | | X | | | |
| 710-2730 | X | | | X | | | | | | | |
| 750-4181 | | | X | X | | | | | | | X |
| 740-1440 | | | | | X | | | | | | |
| 745-1441 | | | X | | | X | X | | | | |
| 740-1000 | | | | | | | | | | | X |
| 745-1000 | | | | | X | | X | | | | |
| ClearGard™ | | | | | | | | | | | |
| 755-2002 | X | X | | | | | | | | | |
| ClearSite® TM | | | | | | | | | | | |
| 741-1022 | X | | | | | | | X | | | |
| 741-1023 | X | | | | | | | | X | | |
| 741-1034 | X | X | | | | X | | | X | | |
| 741-1044 | | X | X | X | X | X | | | X | | |
| 741-1045 | | | | X | X | | | | | X | |
| 741-1040 | | | | | | | | | | X | |
| 741-1068 | | | | | | | X | | | X | |
| 741-1080 | | | | | | | | | | X | |
| SureTape™ | | | | | | | | | | | |
| 800-0403 | X | | | | X | | | | | X | X |
| 800-0406 | | | | | X | | | | | X | X |
| 800-0410 | | | | | | | | | | X | X |

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • IMAGING • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS



525 French Road
 Utica, NY 13502
 (315) 797-8375
 Customer Service: 1-800-448-6506
 FAX: 1-800-438-3051
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
 email: info@mail.conmed.com
 www.conmed.com

