

**Designing, developing, and marketing the most complete line of multifunction electrodes.**

- ConMed's Direct Connect electrodes are comparable to Zoll®, Medtronic® QUIK-COMBO™, Philips Codemaster, and Philips Heartstart electrodes.
- ConMed electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.
- With ConMed's radiotransparent, radiotransparent sterile, and electrodes with carbon leadwires, The Difference is Clear.®
- Radiotranslucent pediatric, and mini-pediatric electrodes are available.
- Electrodes with Pre-Connect packaging are available.
- ConMed offers adapters or ConMed cables, allowing multiple brands of defibrillators to use the same electrodes, with either the Anderson Connection or R2® connection.



**PADPRO®**

\* See reverse side for trademark information.

**CONMED™**  
PATIENT CARE



# Adult and Pediatric Multifunction Electrodes for Defibrillation/Pacing/Cardioversion/Monitoring

**Size A**



**Size A/B**



**Size C**



**Size D**



**Pre-Connect Packaging**



ConMed electrodes have a 24-month shelf life.

## ORDERING INFORMATION

Product Description	Pouch Color	Background Label Color	Size	Anderson Connector <sup>1</sup>	Zoll <sup>®</sup> <sup>2</sup>	Medtronic <sup>®</sup> <sup>3</sup>	Philips Barrel-Style <sup>4</sup>	Philips Plug-Style <sup>5</sup>	R2 <sup>®</sup> Connector <sup>6</sup>
Radiotransparent, Adult	Silver	White	A	2001	2001Z	2001M	2001P	2001H	2001R
Radiotransparent, Adult, w/carbon leadwires	Silver	White	A	2001-C	2001Z-C	2001M-C	2001P-C	2001H-C	N/A
Radiotransparent, Adult, w/Pre-Connect Packaging	Silver	White	A	2001-PC	2001Z-PC	2001M-PC	2001P-PC	N/A	N/A
Radiotransparent, Adult, w/large-posterior pad	Silver	White	A/B	2001-EPS	N/A	N/A	N/A	N/A	N/A
Radiotransparent, Adult w/Sterile Packaging	Silver	White	A	2502	2502Z	2502M	2502P	N/A	N/A
Radiotranslucent, Adult	White	White	A	2516	2516Z	2516M	2516P	2516H	2516R
Radiotranslucent, Adult, w/Pre-Connect Packaging	White	White	A	2516-PC	2516Z-PC	2516M-PC	2516P-PC	N/A	N/A
Radiotranslucent, Pediatric	White	Pink	D	2603	2603Z	2603M	2603P	2603H	2603R
Radiotranslucent, Mini Pediatric	White	Blue	C	2602	2602Z	2602M	2602P	N/A	N/A

**Sizes:**

A = 6.0" (15.2cm) X 4.25" (10.8cm)

B = 9.0" (22.9cm) X 6.0" (15.2cm)

C = 2.5" (6.4cm) X 2.125" (5.4cm)

D = 4.25" (10.8cm) X 2.875" (7.3cm)

\* Product and company names listed are trademarks or trade names of their respective companies.

<sup>1</sup> Electrodes with the Anderson connector can be utilized in environments with multiple brands of defibrillators (via PadPro<sup>®</sup> adapters or PadPro<sup>®</sup> cables).

<sup>2</sup> Zoll models are comparable to Zoll<sup>®</sup> multifunction electrodes.

<sup>3</sup> Medtronic models are comparable to Medtronic Physio Control<sup>®</sup> QUIK-COMBO<sup>™</sup> electrodes.

<sup>4</sup> Philips Barrel-style models are comparable to Philips Codemaster (Barrel-Style) electrodes.

<sup>5</sup> Philips Plug-Style models are comparable to Philips Heartstart (Plug-Style) electrodes.

<sup>6</sup> Electrodes with the R2<sup>®</sup> connector can be utilized in environments with multiple brands of defibrillators (via R2<sup>®</sup> adapters or R2<sup>®</sup> cables).

**CONMED CORPORATION PRODUCT AREAS:**

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com

