



Large Hub Arthroscopy Blade and Bur Chart

For use with C9820, C9824, C9828, C9964, MC9820, MC9828, D9820 and D9824 handpieces


	Specifications	3.5mm	3.7mm	4.0mm	4.2mm	4.5mm	4.8mm	5.5mm	6.0mm
 Great White™ The ultimate blade for aggressive soft tissue resection, including removal of the ACL stump.	Cat. No.: Length: Max Speed:	3.5mm 9399A 13cm 6,000rpm			4.2mm 9299A 13cm 6,000rpm		4.8mm 9297A 13cm 6,000rpm		
 Mako™ Ideal for procedures such as subacromial decompression, acromioplasty and any other procedure where soft tissue and/or bone cutting is desired.	Cat. No.: Length: Max Speed:				4.2mm 9298A 13cm 6,000rpm			5.5mm C9298A 13cm 6,000rpm	
 Cuda® All purpose, aggressive blade. Ideal for extensive tissue removal.	Cat. No.: Length: Cat. No.: Length: Max Speed:	3.5mm C9253 13cm			4.2mm C9254 13cm EL9356 19cm 6,000rpm		4.8mm C9258 13cm	5.5mm C9255 13cm	
 Gator® Ideal for aggressive meniscal trimming, joint debridement, plica and synovium removal, and bone resection.	Cat. No.: Length: Max Speed:	3.5mm C9264 13cm 6,000rpm			4.2mm 9263A 13cm 6,000rpm		4.8mm C9262 13cm 6,000rpm	5.5mm 9260A 13cm 6,000rpm	
 Ultracut™ Excellent for general debridement applications including soft tissue and bone resection.	Cat. No.: Length: Max Speed:				4.2mm C9405 13cm 6,000rpm				
 Full-Radius Resector Excellent for all general debridement applications including soft tissue and bone resection.	Cat. No.: Length: Max Speed:	3.5mm C9248 13cm 6,000rpm			4.2mm 9247A 13cm 6,000rpm		4.8mm C9245 13cm 6,000rpm	5.5mm 9246A 13cm 6,000rpm	
 Dragon® End Cutter Multi-purpose blade ideal for tough tissue resections.	Cat. No.: Length: Max Speed:	3.5mm C9290 13cm 6,000rpm			4.2mm C9291 13cm 6,000rpm			5.5mm C9292 13cm 6,000rpm	
 End Cutter Useful for all general debridement applications.	Cat. No.: Length: Max Speed:	3.5mm C9283 13cm 6,000rpm			4.2mm C9284 13cm 6,000rpm				
 Meniscus Cutter Useful for all meniscal trimming and joint debridement of the patella and trochlear chondromalacia.	Cat. No.: Length: Max Speed:	3.5mm C9760 13cm 6,000rpm			4.2mm 9761A 13cm 6,000rpm			5.5mm C9762 13cm 6,000rpm	
 Slotted Whisker Non-aggressive, for debridement of the fine synovium, fatty tissue and articulating surfaces.	Cat. No.: Length: Max Speed:				4.2mm C9274 13cm 6,000rpm			5.5mm C9275 13cm 6,000rpm	
 Oval Ideal for resecting bone and cartilage during acromioplasty and notchplasty.	Cat. No.: Length: Max Speed:			4.0mm H9101 H9101RH* 13cm 10,000rpm					6.0mm H9102 13cm 10,000rpm
 Spherical Useful for bony site preparation, cartilage and osteochondral debridement, and osteophyte resection.	Cat. No.: Length: Max Speed:	3.5mm H9110 13cm 10,000rpm				4.5mm H9111 13cm 10,000rpm		5.5mm H9112 H9112RH* 13cm 10,000rpm	
 Cyclone® Useful for aggressive bony site preparation, cartilage and osteochondral debridement, and osteophyte resection.	Cat. No.: Length: Max Speed:					4.5mm H9118 H9118RH* 13cm 10,000rpm		5.5mm H9119 13cm 10,000rpm	
 Tapered Excellent for shaping and resecting bone and cartilage during acromioplasty and notchplasty.	Cat. No.: Length: Max Speed:			4.0mm H9121 H9121RH* 13cm 10,000rpm					6.0mm H9122 13cm 10,000rpm
 Vortex® Router Ultra-aggressive. Ideal for extensive bone and cartilage removal in the subacromial space.	Cat. No.: Length: Max Speed:					4.5mm H9131 13cm 10,000rpm			6.0mm H9132 13cm 10,000rpm
 Lightning™ Ideal for procedures such as subacromial decompression, acromioplasty and any other procedure where both soft tissue and bone cutting is desired.	Cat. No.: Length: Max Speed:					4.5mm H9133 13cm 10,000rpm			
 Great White™ The ultimate blade for aggressive soft tissue resection, including removal of the ACL stump.	Cat. No.: Angle: Length: Max Speed:				4.2mm C9299A Pre-Bent 15° 13cm 6,000rpm				
 Cuda® All purpose aggressive blade. Ideal for extensive tissue removal.	Cat. No.: Angle (max): Length: Max Speed: Cat. No.: Angle:		3.7mm P9355 Bendable 30° 13cm 6,000rpm		4.2mm P9356 Bendable 30° 13cm 6,000rpm C9256 Pre-Bent 15°		4.8mm P9357 Bendable 20° 13cm 6,000rpm		
 Gator® Ideal for aggressive meniscal trimming, joint debridement, plica and synovium removal, and bone resection.	Cat. No.: Angle (max): Length: Max Speed: Cat. No.: Angle:		3.7mm P9365 Bendable 30° 13cm 6,000rpm		4.2mm P9366 Bendable 30° 13cm 6,000rpm C9266 Pre-Bent 15°		4.8mm P9367 Bendable 20° 13cm 6,000rpm		
 Full-Radius Resector Excellent for all general debridement applications including soft tissue and bone resection.	Cat. No.: Angle (max): Length: Max Speed: Cat. No.: Angle:		3.7mm P9375 Bendable 30° 13cm 6,000rpm		4.2mm P9376 Bendable 30° 13cm 6,000rpm C9241 Pre-Bent 15°		4.8mm P9377 Bendable 20° 13cm 6,000rpm		

Very Aggressive

Aggressive

Least Aggressive

Burs and Routers



For use with Merlin® bendable blades P9305



Linvatec
A CONMED® Company

11311 Concept Boulevard,
Largo, Florida 33773-4908
Phone: (727) 392-6464
Customer Service: (800) 237-0169
USA Fax: (727) 399-5256
International Fax: (727) 397-4540

Sterling