

**PRECISION
AND FLEXIBILITY
IN HIP ARTHROSCOPY**



HIP ARTHROSCOPY DISPOSABLE KIT

The ConMed Linvatec Hip Arthroscopy Disposable Kit includes everything needed to access the hip joint during hip arthroscopy procedures. The kit includes the following components:

- (2) Hip Arthroscopy Cannulae, 8.0 x 120mm
- Distention Needle, 17 Gauge 1.47 x 152.4mm
- Disposable Syringe, 30mL
- (2) Nitinol Guidewires, 1.14 x 355.6mm
- Cannulated Switching Stick, 3.48 x 317.5mm
- Disposable Arthroknife®, Banana Blade, 5.59 x 254mm

Note: Hip Arthroscopy Cannulae are latex free.

ORDERING INFORMATION

Description	Cat. No.
Hip Arthroscopy Disposable Kit, 1ea (Hip Arthroscopy Reusable Obturator sold separately)	HA-222DK
Accessories	
Hip Arthroscopy Reusable Obturator, 1ea	9837
Straight Blades	
4.2mm Great White® Extended Length Non-bendable Blade, 19cm	EL9299
4.2mm Gator® Extended Length Non-bendable Blade, 19cm	EL9263A
5.5mm Gator Extended Length Non-bendable Blade, 19cm	EL9260A
Pre-Bent Blades	
4.2mm Great White Extended Length Pre-Bent 15 degree Blade, 19cm (window facing down)	EL9299D
4.2mm Great White Extended Length Pre-Bent 15 degree Blade, 19cm (window facing up)	EL9299U
4.2mm Full-Radius Resector Extended Length Pre-Bent 15 degree Blade, 19cm (window facing down)	EL9247D
4.2mm Full-Radius Resector Extended Length Pre-Bent 15 degree Blade, 19cm (window facing up)	EL9247U
Bendable Blades	
4.2mm Poly Cuda® Extended Length Bendable Blade, 19cm	EL9356
Burs	
4.0mm Oval Extended Length Bur, 19cm	EL9101
4.5mm Spherical Extended Length Reduced Hood Bur, 19cm	EL9111RH



Hip Arthroscopy Reusable Obturator



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
 www.linvatec.com