

With a variety of shapes, sizes, types, gels, materials, adhesion levels, and applications to offer, CONMED can meet all of your ECG electrode monitoring needs. One company with one goal — to satisfy our customer in the most cost-effective manner possible.



Adult ECG Electrodes

CONMED has got you covered...

CONMED
PATIENT CARE

ADULT ECG ELECTRODES — TAPE

CONDUCTIVE ADHESIVE GEL ELECTRODES



Cleartrace®

- Clear Neoderm® tape allows visual monitoring of the patient's skin condition at all times.
- Perforated tape permits skin ventilation and moisture evaporation for increased patient comfort.
- Conforms to patient's body without pulling on the skin.
- Will not absorb fluids.



Cleartrace® 2



Cleartrace® 2 LT

- Square/Oval shape adds "staying" power to the design — up to seven days of use.
- Clear Neoderm® tape allows visual monitoring of the patient's skin condition at all times.
- Neoderm® tape is waterproof yet allows moisture to escape.
- Available with radiotranslucent (Cleartrace® 2) or metal (Cleartrace® 2 LT) snap.

WET GEL ELECTRODES



Invisatrace®



Ultratrace®

- Wet gel formulation for strong signals and stable baselines.
- Available with radiotranslucent (Invisatrace®) or metal (Ultratrace®) snap.
- Tab for easy application and removal.
- Clear Neoderm® tape gives patient comfort and fluid-proof performance.
- Ability to monitor skin condition at all times.

ADULT ECG ELECTRODES — CLOTH



Plia-Cell®



Plia-Cell® Diaphoretic

- Oval shape allows for convenient V-lead placement.
- ClearSite® solid adhesive gel is gentle to the patient.
- Radiotranslucent snap.
- Low profile cloth conforms to body for increased patient comfort.
- Available with a normal (Plia-Cell®) or a diaphoretic (Plia-Cell® Diaphoretic) adhesive.

SPECIALTY ELECTRODES



Special Studies Electrode



ECG Backpads



ECG Strippers

- The SSE Special Studies Electrode is a unique radiotransparent electrode providing uninterrupted views of the chest, when utilized with unique radiotransparent leadwires.
- The ECG Backpad is an integrated array of electrodes designed for placement on the back for maximum chest exposure.
- ECG Strippers are designed for ease and speed of set up while eliminating tangled leadwires.

ADULT ECG ELECTRODES — FOAM

CONDUCTIVE ADHESIVE GEL ELECTRODES



Suretrace®



Positrate®

- Closed cell foam keeps fluids from penetrating.
- Pull tab makes it easy to apply and remove the electrode.
- Available in normal (Suretrace®) and diaphoretic (Positrate®) adhesive.



Silvon®



Silvon® Diaphoretic

- Low profile foam conforms to body for increased patient comfort.
- Oval shape allows for convenient V-lead placement.
- ClearSite® solid adhesive gel is gentle to the patient.
- Radiotranslucent snap.
- Available in normal (Silvon®) and diaphoretic (Silvon® Diaphoretic) adhesive.



Accutac®




Accutac® Diaphoretic




High Demand® II

- Transparent film repels fluid that might otherwise cause loss of contact.
- Non-irritating ClearSite® gel adheres gently to skin.
- Radiotranslucent snap.
- Available in normal (Accutac®), diaphoretic (Accutac® Diaphoretic), and diagnostic (High Demand® II) adhesive.

WET GEL ELECTRODES



Exactrace



Instatrace®

- Closed cell foam keeps fluids from penetrating.
- Pull tab makes it easy to apply and remove the electrode.
- Available in normal (Exactrace) and diaphoretic (Instatrace®) adhesive.
- Wet gel formula provides quick signal retention.



Dyna/Trace®



Dyna/Trace® Stress



Dyna/Trace® Mini

- Wet gel formula provides quick signal retention.
- Available in normal (Dyna/Trace®) and aggressive (Dyna/Trace® Stress) adhesion for all levels of application.
- Dyna/Trace® Mini is designed to be used along the V-leads thereby reducing the need to overlap.
- Radiotranslucent snap.

RESTING EKG ELECTRODES



Fastrace® 4

- Silver/Silver Chloride provides stable baselines and excellent electrical characteristics.
- Conductive adhesive gel maintains contact with the skin, even on hairy patients.
- Enhanced substrate extends shelf life (36 months).



Snaptrace® II

- The small design is ideal for adults and pediatrics and offers easy placement for V-leads.
- Solid adhesive hydrogel provides staying power, but gentle release.
- Radiotranslucent snap.

Adult ECG Electrodes

Adult ECG Electrodes — Tape

Product	Product Number	Packaging		
		Pouch	Box	Case
Cleartrace®	1700-001	1	30	600
	1700-003	3	30	600
	1700-005	5	50	600
	1700-030	30	—	600
Cleartrace® 2	2700-001	1	30	600
	2700-003	3	30	600
	2700-005	5	50	600
	2700-030	30	—	600
Cleartrace® 2 LT	2710-001	1	30	600
	2710-003	3	30	600
	2710-005	5	50	600
	2710-030	30	—	600
Invisatrace®	1680-001	1	30	600
	1680-003	3	30	600
	1680-005	5	50	600
	1680-030	30	—	600
Ultratrace®	1690-001	1	30	600
	1690-003	3	30	600
	1690-005	5	50	600
	1690-030	30	—	600

Adult ECG Electrodes — Cloth

Plia-Cell®	01-3810	1	30	900
	01-3820	25	—	250
	01-3830	3	30	900
Plia-Cell® Diaphoretic	01-3310	1	30	900
	01-3320	25	—	250
	01-3330	3	30	900

Specialty Electrodes

SSE	01-5010	1	30	300
	01-5050	50	—	500
Backpad	01-3130	1	—	50
Strippers (5-lead)	38-402-1	1	25	100

Resting EKG Electrodes

Fastrace® 4	1915-100	100	—	1000
Snaptrace® II	1921-100	100	—	1000

Adult ECG Electrodes — Foam

Product	Product Number	Packaging		
		Pouch	Box	Case
Suretrace®	1800-001	1	30	600
	1800-003	3	30	600
	1800-005	5	50	600
	1800-030	30	—	600
Positrac® (Diaphoretic)	1870-001	1	30	600
	1870-003	3	30	600
	1870-005	5	50	600
	1870-030	30	—	600
Silvon®	01-3510	1	30	900
	01-3530	3	30	900
	01-3533	30	—	600
	01-3540	4	40	1200
	01-3550	5	50	1500
Silvon® Diaphoretic	01-3610	1	30	900
	01-3630	3	30	900
	01-3633	30	—	600
	01-3640	4	40	1200
	01-3650	5	50	1500
Accutac®	01-8050	50	—	500
Accutac® Diaphoretic High Demand® II	01-8150	50	—	500
	01-6150	50	—	500
Exactrace	1510-001	1	30	600
	1510-003	33	30	600
	1510-005	5	50	600
	1510-030	30	—	600
Instatrace® (Diaphoretic)	1590-001	1	30	600
	1590-003	3	30	600
	1590-005	5	50	600
	1590-030	30	—	600
Dyna/Trace®	1500-001	1	30	600
	1500-003	3	30	600
	1500-004	4	40	480
	1500-005	5	50	600
	1500-030	30	—	600
Dyna/Trace® Stress	1550-001	1	30	600
	1550-005	5	50	600
Dyna/Trace® Mini	1555-005	5	50	600

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

© 2006 CONMED Corporation, 03/2006, Control #8443 Rev C

Ag/AgCl