


With a variety of shapes, sizes, types, gels, materials, adhesion levels, and applications to offer, CONMED can meet all of your ECG electrode monitoring needs. One company with one goal — to satisfy our customer in the most cost-effective manner possible.



Adult ECG Electrodes

CONMED has got you covered...



ADULT ECG ELECTRODES — TAPE

CONDUCTIVE ADHESIVE GEL ELECTRODES



Cleartrace™



Cleartrace™ Diagnostic

- Clear Neoderm® tape allows visual monitoring of the patient's skin condition at all times.
- Perforated tape permits skin ventilation and moisture evaporation for increased patient comfort.
- Conforms to patient's body without pulling on the skin.
- Will not absorb fluids.
- Available in general purpose (Cleartrace™) and diagnostic (Cleartrace™ Diagnostic) adhesive.

LONG TERM MONITORING ELECTRODES



Cleartrace™ 2



Cleartrace™ 2 LT

- Square/Oval shape adds "staying" power to the design — up to seven days of use.
- Clear Neoderm® tape allows visual monitoring of the patient's skin condition at all times.
- Neoderm® tape is waterproof yet allows moisture to escape.
- Available with radiotranslucent (Cleartrace™ 2) or metal (Cleartrace™ 2 LT) snap.

WET GEL ELECTRODES



Invisatrace®



Ultratrace®

- Wet gel formulation for strong signals and stable baselines.
- Available with radiotranslucent (Invisatrace®) or metal (Ultratrace®) snap.
- Tab for easy application and removal.
- Clear Neoderm® tape gives patient comfort and fluid-proof performance.
- Ability to monitor skin condition at all times.

REPOSITIONABLE ECG ELECTRODES

VINYL TAPE



Cleartrode™



Cleartrode™ RTL

- Transparent tape allows visual monitoring of the skin under the electrode.
- Conductive adhesive gel does not dry out and allows for extended wear without irritating the baby's skin.
- Cleartrode® electrodes are re-wettable for cost effective long term use and ease of application and removal.
- Available with durable nickel-plated brass snap or carbon-based radiotranslucent (RTL) snap.
- Latex free to protect against possible latex sensitivity issues.
- PVC free.

SOFTFACE® CLOTH



Softface® Small



Softface® Small RTL



Softface® Medium



Softface® Medium RTL



Softface® Large



Softface® Large RTL

- Soft fabric material and mild conductive adhesive provide secure, comfortable electrode contact.
- Humidity has little effect on Softface® electrodes' ability to stick or to be removed without damaging fragile skin.
- Conductive adhesive gel does not dry out and allows for extended wear without irritating the baby's skin.
- Softface® electrodes are re-wettable for cost effective long term use and ease of application and removal.
- Available with durable nickel-plated brass snap or carbon-based radiotranslucent (RTL) snap.
- Available in three sizes (small, medium and large) for optimal placement and application for all patient sizes.
- Latex free to protect against possible latex sensitivity issues.
- PVC free.

ADULT ECG ELECTRODES — FOAM

CONDUCTIVE ADHESIVE GEL ELECTRODES



Suretrace®



Positrace®

- Closed cell foam keeps fluids from penetrating.
- Pull tab makes it easy to apply and remove the electrode.
- Available in general purpose (Suretrace®) and diaphoretic (Positrace®) adhesive.



Silvon®



Silvon® Diaphoretic

- Low profile foam conforms to body for increased patient comfort.
- Oval shape allows for convenient V-lead placement.
- ClearSite® solid adhesive gel is gentle to the patient.
- Radiotranslucent snap.
- Available in general purpose (Silvon®) and diaphoretic (Silvon® Diaphoretic) adhesive.

WET GEL ELECTRODES



Exactrace®



Instatrace®

- Closed cell foam keeps fluids from penetrating.
- Pull tab makes it easy to apply and remove the electrode.
- Available in general purpose (Exactrace®) and diaphoretic (Instatrace®) adhesive.
- Wet gel formula provides quick signal retention.



Dyna/Trace®



Dyna/Trace® Stress



Dyna/Trace® Mini

- Wet gel formula provides quick signal retention.
- Available in general purpose (Dyna/Trace®) and diagnostic (Dyna/Trace® Stress) adhesion for all levels of application.
- Dyna/Trace® Mini offers easy placement for V-leads.
- Radiotranslucent snap.

SPECIALTY ELECTRODES



Special Studies Electrode



ECG Backpads



ECG Strippers

- The SSE Special Studies Electrode is a unique radiotransparent electrode providing uninterrupted views of the chest, when utilized with unique radiotransparent leadwires.
- The ECG Backpad is an integrated array of electrodes designed for placement on the back for maximum chest exposure.
- ECG Strippers are designed for ease and speed of set up, elimination of tangled leadwires, and minimization of cross contamination.

RESTING EKG ELECTRODES



Fastrace® 4

- Silver/Silver Chloride provides stable baselines and excellent electrical characteristics.
- Conductive adhesive gel maintains contact with the skin, even on hairy patients.
- PVC free.



Snaptrace® II

- The small design is ideal for adults and pediatrics and offers easy placement for V-leads.
- Solid adhesive hydrogel provides staying power, but gentle release.
- Radiotranslucent snap.
- PVC free.

Adult ECG Electrodes

Adult ECG Electrodes — Tape

Product	Product Number	Packaging		
		Pouch	Box	Case
Cleartrace™	1700-001	1	30	600
	1700-003	3	30	600
	1700-004	—	—	600
	1700-005	5	50	600
	1700-030	30	—	600
	1700-050	50	—	1000
Cleartrace™ Diagnostic	1750-001	1	30	600
	1750-003	3	30	600
	1750-005	5	50	600
	1750-030	30	—	600
Cleartrace™ 2	2700-001	1	30	600
	2700-003	3	30	600
	2700-005	5	50	600
	2700-030	30	—	600
Cleartrace® 2 LT	2710-001	1	30	600
	2710-003	3	30	600
	2710-005	5	50	600
	2710-030	30	—	600
Invisatrace®	1680-001	1	30	600
	1680-003	3	30	600
	1680-005	5	50	600
	1680-030	30	—	600
Ultratrace®	1690-001	1	30	600
	1690-003	3	30	600
	1690-005	5	50	600
	1690-030	30	—	600

Repositionable ECG Electrodes

Cleartrode™	1720-003	3	—	300
Cleartrode™ RTL	1720C-003	3	—	300
Softrace® Small	2320-003	3	—	300
Softrace® Small RTL	2320C-003	3	—	300
Softrace® Medium	2330-003	3	—	300
Softrace® Medium RTL	2330C-003	3	—	300
Softrace® Large	2350-003	3	—	300
	2350-005	5	—	600
Softrace® Large RTL	2350C-003	3	—	300
	2350C-005	5	—	600

Adult ECG Electrodes — Foam

Product	Product Number	Packaging			
		Pouch	Box	Case	
Suretrace®	1800-001	1	30	600	
	1800-003	3	30	600	
	1800-005	5	50	600	
	1800-030	30	—	600	
	Positrace® (Diaphoretic)	1870-001	1	30	600
		1870-003	3	30	600
1870-005		5	50	600	
1870-030		30	—	600	
Silvon®	01-3510	1	30	900	
	01-3530	3	30	900	
	01-3533	30	—	600	
	01-3540	4	40	1200	
	01-3550	5	50	1500	
Silvon® Diaphoretic	01-3610	1	30	900	
	01-3630	3	30	900	
	01-3633	30	—	600	
	01-3640	4	40	1200	
	01-3650	5	50	1500	
Exactrace	1510-001	1	30	600	
	1510-003	33	30	600	
	1510-005	5	50	600	
	1510-030	30	—	600	
Instatrace® (Diaphoretic)	1590-001	1	30	600	
	1590-003	3	30	600	
	1590-005	5	50	600	
	1590-030	30	—	600	
Dyna/Trace®	1500-001	1	30	600	
	1500-003	3	30	600	
	1500-004	4	40	480	
	1500-005	5	50	600	
	1500-030	30	—	600	
Dyna/Trace® Stress	1550-001	1	30	600	
	1550-005	5	50	600	
Dyna/Trace® Mini	1555-005	5	50	600	

Specialty Electrodes

SSE	01-5010	1	30	300
	01-5050	50	—	500
Adult Backpad	01-3130	1	—	50
Strippers (5-lead)	38-402-1	1	25	100
Strippers (4-lead)	38-404-1	1	25	100

Resting EKG Electrodes

Fastrace® 4	1915-100	100	—	1000
Snaptrace® II	1921-100	100	—	1000

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com

© 2008 CONMED Corporation, 01/2008, Control #8443 Rev D

Ag/AgCl  

MBR 9074 Rev. 1