



DURAglide[™]

Stone Removal Balloons



Applied Ingenuity[™] |

CONMED[™]
ENDOSCOPIC TECHNOLOGIES

DURAglide™

Stone Removal Balloons

DURAglide 3™



.035"
COMPATIBLE


Triple Lumen Stone Removal Balloons

Ordering Information

Catalog No.	Balloon		Catheter		Guidewire Compatibility
	O.D.(mm)	VOL.(cc)	O.D.	Length	
Injection Above Balloon					
007123	8.5	.75	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007124	11.5	1.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007125	15	3.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
Injection Below Balloon					
007143	8.5	.75	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007144	11.5	1.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
007145	15	3.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm

Performance Guaranteed!

Triple Lumen



DURAglide



.035"
COMPATIBLE



Double Lumen Stone Removal Balloons

Ordering Information

Catalog No.	Balloon		Catheter		Guidewire Compatibility
	O.D.(mm)	VOL.(cc)	O.D.	Length	
Standard					
000637	11.5	1.5	7F/2.3 mm	200 cm	.035"/.89 mm
000638	15	3.5	7F/2.3 mm	200 cm	.035"/.89 mm
Tapered					
000635	8.5	.75	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
000646	11.5	1.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm
000647	15	3.5	7-5F/2.3-1.7 mm	200 cm	.035"/.89 mm

Performance Guaranteed!

Double Lumen



- Increased Injection Lumen Capacity
Increased Flow Reduces Hand Fatigue
- Color Coded Guidewire Port
Visual Indicator Eliminates Confusion
- 3 Radiopaque 1cm Spaced Distal Marker Bands
Enables Fluoroscopic Visualization of Balloon Position
- Pebax Catheter Shaft
Enhances Pushability



- Standard and Tapered Catheter Options
Accommodates a Variety of Duct or Stone Sizes
- 3 Radiopaque 1cm Spaced Distal Marker Bands
Enables Fluoroscopic Visualization of Balloon Position
- Pebax Catheter Shaft
Enhances Pushability

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375
Customer Service: 1-800-448-6506
FAX: 1-800-438-3051
International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235
email: info@mail.conmed.com
www.conmed.com