


# Reusable Electrosurgical Pencils

# Reusable Electrosurgical Pencils



## Ordering Information

Catalog Number	Product Description	Packaging
130316	Rocker Switch, Disposable Blade Electrode, 10' (3.05m) Cable, Sterile validated to 20 uses when utilized in accordance with the sterilization instructions provided in the package insert	10/case 1/pkg
130316-15	Rocker Switch, Disposable Blade Electrode, 15' (4.57m) Cable, Sterile validated to 20 uses when utilized in accordance with the sterilization instructions provided in the package insert	10/case 1/pkg
130317	Push Button, Disposable Blade Electrode, 10' (3.05m) Cable, Sterile validated to 20 uses when utilized in accordance with the sterilization instructions provided in the package insert	10/case 1/pkg
130317-15	Push Button, Disposable Blade Electrode, 15' (4.57m) Cable, Sterile validated to 20 uses when utilized in accordance with the sterilization instructions provided in the package insert	10/case 1/pkg
130384	PowerPoint® Plus, Push Button, Disposable Blade Electrode, 15' (4.57m) Cable, Non-Sterile validated to 50 uses when utilized in accordance with the sterilization instructions provided in the package insert	10/case 1/pkg
60-0815-001	Foot Control Pencil, 10' (3.05m) Cord, Non-Sterile, Sold without an electrode	1 each
60-0235-001	Reusable Safety Holster	10/case 1/pkg
137668	Disposable Safety Holster	40/case 1/pkg
138029	Tip Cleaner, Single Use, Sterile	40/case 1/pkg

Please contact your local ConMed Territory Manager for information on our complete product line and accredited training programs.

### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: [info@mail.conmed.com](mailto:info@mail.conmed.com)  
[www.conmed.com](http://www.conmed.com)