



ABCFlex™ Argon Beam GI Probe

## ABCFlex™ Argon Beam GI Probes

ConMed Electrosurgery provides advanced technology for the gastroenterologist with ABCFlex™ Probes and the System 7550™ ABC® ESU or System 7500® ABC® ESU. Argon Beam Coagulation® technology enables the physician to achieve non-contact hemostasis with minimal thermal injury of adjacent or underlying tissue.

ABC® technology is ideally suited for GI applications and is a valuable tool for the gastroenterologist. Like standard electrosurgery, energy is delivered to tissue to create hemostasis, but the argon beam uses a focused column of ionized argon gas at room temperature to deliver the energy. The resulting cool argon beam moves blood away from the bleeding tissue and coagulates directly to the mucosal surface which is intended to limit thermal spread and may result in less necrosis. It can be used to create superficial ablation of mucosal tissue over a large area, or remedy Arteriovenous Malformations (AVMs).

Two versions of the ABCFlex™ Probes are available. The standard version (#133023 and #133270) is designed to work with ConMed's System 7550™ and System 7500®. Additional versions of the ABCFlex™ Probes (#133323 and #133327) are designed to work exclusively with the System 7550™. The additional ABCFlex™ Probes are constructed to provide better initiation as well as eliminate leaks when attached to the System 7550™.

## Intended Use

Products #133023 and #133270 are intended for use with the System 7500® and System 7550™ as flexible catheters for non-contact superficial ablation and hemostasis. These devices are used exclusively through an endoscope.

Products #133323 and #133327 are intended for use with the System 7550™ as flexible catheters for non-contact superficial ablation and hemostasis. These devices are used exclusively through an endoscope.

Please consult product labels and inserts for any indications, contradictions, hazards, warnings, precautions and instructions for use.



## Ordering Information

Product	Product Description	Quantity
133023	ABCFlex™ Probe 2.3mm (7Fr.) x 220cm, Original	10/case
133270	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, Original	10/case
133323	ABCFlex™ Probe 2.3mm (7Fr.) x 220cm, System 7550™	10/case
133327	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, System 7550™	10/case
133000	Extension Adapter for use with ABCFlex™ Probes	1 each

Please contact your ConMed Territory Manager for information on our complete product line and accredited training program.

### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road  
 Utica, NY 13502  
 (315) 797-8375  
 Customer Service: 1-800-448-6506  
 FAX: 1-800-438-3051  
 International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235  
 email: info@mail.conmed.com  
 www.conmed.com

