

THE CORE ENTREE® II REPOSABLE BODY ACCESS SYSTEM

- **Modular Reposable System** provides choice with flexibility and value.
- **Disposable One-Hole Reducer and Valve** accepts 5mm – 12mm instruments.
- **Disposable Audible Dilating Trocar** minimizes the trocar wound defect while dilating the remaining tissue.
- **Reusable Low-Profile Metal Cannula** enables ease of instrument exchange and includes a lifetime warranty.
- **Disposable Procedure-Specific Trocar Kits** provide convenience and additional cost savings




THE CORE ENTREE® II


Disposable Feel and Performance...
Reusable Price

The Core Entree® II Reposable Body Access System

ENTREE® II Cannula and Cannula with Fascia Threads, Reusable




Cannula




Cannula with Fascia Threads

ORDER CODE		DESCRIPTION
Smooth	Threaded	
CD1605	CD1605FH	5mm, Without Stopcock, SHORT Length
CD1606	CD1606FH	5mm, SHORT Length
CD1608	CD1608FH	5mm, Without Stopcock, LONG Length
CD1609	CD1609FH	5mm, LONG Length
CD1607	CD1607FH	7/8mm, 90mm Length
CD1610	CD1610FH	10mm, 110mm Length
CD1612	CD1612FH	10/12mm, 100mm Length

Trocar Procedure Kit, 5mm, Sterile



ORDER CODE	DESCRIPTION
CD7100G	1 ea, 5mm CORE Audible Dilating Trocar, SHORT Length 2 ea, 5mm ENTREE II Cannula Valve



CD7150G	1 ea, 5mm CORE Audible Dilating Trocar, LONG Length 2 ea, 5mm ENTREE II Cannula Valve
---------	--

Trocar Procedure Kit, 7/8mm, Sterile




ORDER CODE	DESCRIPTION
CD7308G	1 ea, 7/8mm CORE Audible Dilating Trocar, 90mm Length 1 ea, 5-12mm ENTREE II Valve/Reducer

Trocar Procedure Kit, 10mm, Sterile




ORDER CODE	DESCRIPTION
CD7465G	1 ea, 10mm CORE Audible Dilating Trocar, 110mm Length 1 ea, 5-12mm ENTREE II Valve/Reducer

Trocar Procedure Kit, 10/12mm, Sterile




ORDER CODE	DESCRIPTION
CD7565G	1 ea, 10/12mm CORE Audible Dilating Trocar, 100mm Length 1 ea, 5-12mm ENTREE II Valve/Reducer

ENTREE II Valve/Reducer, Sterile



ORDER CODE	DESCRIPTION
CD775	5-12mm

ENTREE II Cannula Valve, Sterile



ORDER CODE	DESCRIPTION
CD650	5mm

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



525 French Road
Utica, NY 13502
(315) 797-8375

Customer Service: 1-800-448-6506

FAX: 1-800-438-3051

International Sales: +1(315) 797-8375 • FAX: +1(315) 735-6235

email: info@mail.conmed.com

www.conmed.com