

SmartScrew® II

Proprietary Self-Reinforced
96L/4D copolymer

Clover shape screw head fits
easily with new center driver

Self retaining screw — screw fits
to the driver and will not drop

Low profile head design

Fully and partially
threaded implants

Cannulated and Non-cannulated

Depth markings in cannulated
drills (drilling deepness can be
estimated from drill bit when
drill/tap guide used)

 Material
Advantage™

 CONMED
LINVATEC

ORDERING INFORMATION:



Instruments for SmartScrew® II Non-Cannulated and Cannulated

Screw Size	Drill Bit	Drill/Tap Guide	Countersink	Depth Gauge	Bone Tap	Drill Bit for Lag Technique	Screw Driver
Non-Cannulated SmartScrew II							
2.0mm	DB5015 (1.5 mm)	GE5020	CS5020	DG5001	BT5020 (2.0 mm)	DB5020 (2.0 mm)	SD5020
2.7mm	DB5020 (2.0 mm)	GE5027	CS5027	DG5001	BT5027 (2.7 mm)	DB5027 (2.7 mm)	SD5027
3.5mm	DB5025 (2.5mm)	GE5035	CS5035	DG5001	BT5035 (3.5 mm)	DB5035 (3.5mm)	SD5035
4.5mm	DB5035 (3.5 mm)	GE5045	CS5045	DG5000	BT5045 (4.5 mm)	DB5045 (4.5mm)	SD5045
Cannulated SmartScrew II*							
3.5mm	DRC5125 (2.5 mm)	GE5035	CSC5135	DG5002	BTC5135 (3.5 mm)	DRC5135 (3.5 mm)	SDC5135
4.5mm	DRC5135 (3.5mm)	GE5045	CSC5145	DG5002	BTC5145 (4.5 mm)	DRC5145 (4.5mm)	SDC5145

3.5 mm K-Wire KWC1200 (1.2 mm)

4.5 mm K- Wire KWC1600 (1.6 mm)

*Cannulated Screws available in 3.5 and 4.5mm diameter only



CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com