




Infinity™ Quad Tendon

Harvesting System




Engineered for consistent
quad graft harvesting

Scan for more information.

Imagine what we could do together.

Designed for Performance

The Infinity™ Quad Tendon Harvesting System supports consistent harvesting of the quad tendon for ACL reconstruction using a simple, repeatable technique. The System includes the following next-generation components:



Infinity™ Disposable Double Blade



Infinity™ Tendon Stripper (left- and right-handed configurations)



Infinity™ Proximal Tendon Cutter

Infinity™ Quad
Tendon Retractor



Imagine what we could do together.

A New Standard in Harvest Visualization

The Infinity™ Quad Tendon Harvesting System features a retractor designed to accommodate arthroscopic visualization during tendon harvest.



Infinity™ Quad
Tendon Retractor



- The adjustable retractor provides **controlled, variable skin retraction**
- The retractor maintains elevation of the surgical incision **without requiring additional instrumentation or assistance**, while a dedicated slot accommodates a scope for **continuous visualization** during tendon harvest



Imagine what we could do together.


Precision, Perfected

The Infinity™ Quad Tendon Harvesting System ensures consistent performance and ease-of-use for every procedure.


Infinity™ Disposable Double Blades

The system utilizes a proximal-to-distal tendon cutting technique, enabling controlled harvesting of the tendon.

The reverse double blades feature a depth limiter that enhances patient safety and delivers consistent 6 mm or 8 mm harvest depth.



The Infinity™ Tendon Stripper features a blunt cutting edge, designed to separate layers and adhesions to help preserve graft integrity.



The Infinity™ Proximal Tendon cutter enables consistent and reliable resection of the proximal end of the quad tendon.

Imagine what we could do together.

Ordering Information

Product Name	Use	SKU
Infinity™ Disposable Double Blade, 6 mm x 9 mm	Disposable	DB-0609
Infinity™ Disposable Double Blade, 6 mm x 10 mm	Disposable	DB-0610
Infinity™ Disposable Double Blade, 6 mm x 11 mm	Disposable	DB-0611
Infinity™ Disposable Double Blade, 8 mm x 8 mm	Disposable	DB-0808
Infinity™ Disposable Double Blade, 8 mm x 9 mm	Disposable	DB-0809
Infinity™ Disposable Double Blade, 8 mm x 10 mm	Disposable	DB-0810
Infinity™ Proximal Tendon Cutter and Retractor Kit	Disposable	QTCR
Infinity™ Attachable Double Blade 8 mm x 8 mm	Disposable	ADB-0808
Infinity™ Attachable Double Blade 8 mm x 9 mm	Disposable	ADB-0809
Infinity™ Attachable Double Blade 8 mm x 10 mm	Disposable	ADB-0810
Infinity™ Attachable Double Blade 6 mm x 9 mm	Disposable	ADB-0609
Infinity™ Attachable Double Blade 6 mm x 10 mm	Disposable	ADB-0610
Infinity™ Attachable Double Blade 6 mm x 11 mm	Disposable	ADB-0611
Infinity™ Reusable Double Blade Handle	Reusable	RDBH
Infinity™ Quad Tendon Stripper, Left	Reusable	QTSL
Infinity™ Quad Tendon Stripper, Right	Reusable	QTSR



Imagine what we could do together.



Imagine what we could do together.

CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com



Service & Support