



Hall® Lithium Battery System


Dependable Power and Convenient Sterilization

Our unique, UL-approved lithium batteries feature twice the power, twice the life¹ and the ability to be immediately autoclaved with other instruments in a single tray.



HALL®
POWERED INSTRUMENTS

M-CLASS®
BLADES



HALL 50™
OSCILLATING SAW



HALL 50™
RECIPROCATING SAW



DRILL-REAM
ATTACHMENTS



To learn more about these and other innovative products, call 800-448-6506 or visit ConMed.com.



POWERED
INSTRUMENTS
SYSTEM™

Dependable Power

Introducing Our Most Powerful Battery System

Hall® Powered Instruments is proud to introduce the next generation of battery innovation for powered instruments. In addition to being lighter, longer-lasting and faster-charging than standard batteries¹, our new lithium batteries are the **first on the market to be UL- approved – providing convenient autoclavability**. Additionally, the system's new sterilization tray is capable of accommodating handpieces, attachments and two batteries for streamlined processing and an *off-the-shelf, ready-to-go solution*. Backward and forward compatible, this system brings a new level of dependable power to every Hall battery-powered handpiece.

Hall® Lithium Battery

The Power To Improve Your Surgical Experience

Lithium-strength performance. Nearly twice the power of standard batteries¹ for uninterrupted use during high torque procedures.

Long-lasting power. Significantly longer battery and product life than standard autoclavable batteries¹ resulting in greater cost savings.

Compact size. 25% lighter than standard autoclavable batteries¹ for a more balanced instrument.

Environmentally responsible. Designed with environmentally acceptable materials, and recyclable.

¹ ConMed autoclavable NiCad PRO3011 battery

Hall® Charger

Battery Power You Can Count On
When It Matters Most



Fast charge time. Just 30 minutes to charge with no need to discharge first, for faster return to service.

Charge indicator. Displays charging progress to indicate battery readiness at a glance.

Advanced diagnostics. Shows temperature exposure and cycle count to help manage battery use.

Space-Saving Design. Wall mount capabilities optimize counter space. Wrap around cord reduces clutter.

POWERED
INSTRUMENTS
SYSTEM™



H Handpieces

B Batteries

A Attachments

Hall® Sterilization Tray

One-Step Sterilization

Unlike competitive batteries that need to be autoclaved separately at different sterilization cycles, **our batteries can actually be placed directly in the same tray as other instruments immediately after every procedure.** Our sterilization tray is designed to house two lithium batteries, eight attachments, two pin/wire drivers and three handpieces. Together, this provides OR staff members a new level of convenience that helps save time and improve efficiency.

HALL®
POWERED INSTRUMENTS



525 French Road, Utica, NY 13502

POWERED
INSTRUMENTS
SYSTEM™

HALL®
POWERED INSTRUMENTS

Ordering Information

Description	Catalog Number
Hall Lithium Battery	L3000LG
Hall Lithium Charger	L3000
Hall Sterilization Tray	PRO6001

Hall® Lithium Battery System

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

ConMed.com
info@mail.conmed.com