



NEW

For ablation... CONMED has the Edge™

Edge™
Bipolar Arthroscopic
Radiofrequency System



To learn more about this and other innovative products,
call 1-800-237-0169 or visit ConMed.com.

Edge™

Bipolar Arthroscopic Radiofrequency System


Faster Bulk Tissue Ablation with less Risk

The unique one piece electrode designed to prevent probe degradation allows the Edge bipolar Arthroscopic Radiofrequency system to safely use more power. This generates a stronger and more stable plasma field which data shows has a higher grams per minute ablation rate with a reduced risk of probe degradation.²

CONMED Edge 90S Probe




Smith and Nephew Super TurboVac 90 Probe



Better Visualization

The stronger and more stable plasma field allows for a greater aspiration rate without sacrificing ablation performance. Data shows¹ the large cross shape suction port of the AES-90S aspirates a higher mL per second than the Smith and Nephew Super TurboVac 90.

Ablation Rate in gm/min



¹Data on file SPR-13-157 (2014) ²Data on file TR-14-294 (2015)

To schedule an evaluation, call **800-433-6576**
or visit CONMED.com for more information.



525 French Road, Utica, NY 13502

Customer Service: 1-800-448-6506
International Sales: 1-315-797-8375
Fax: 1-800-438-3051

CONMED.com
info@mail.conmed.com