



Designed for Use with
the BioBrace® Implant

Radiolucent PEEK Material

Strong Pullout Resistance

Delivered Sterile and
Pre-assembled

Hercules® XE

Extended Eyelet Knotless Anchor for BioBrace®

PROCEDURE KIT NOW AVAILABLE




Extended eyelet specifically
designed for the 5mm x 250mm
BioBrace® Reinforced Implant

PRODUCT OVERVIEW





Hercules® XE

Extended Eyelet Knotless Anchor for BioBrace®




BioBrace® 5mm x 250mm Implant



Anterior Talofibular Ligament (ATFL)

The **Hercules® XE Knotless Anchor** features an extended eyelet intended for fixation of suture (soft tissue) to bone in the foot/ankle, and designed specifically for use with the 5mm x 250mm **BioBrace® Reinforced Implant** in medial and lateral ankle stabilizations and Achilles tendon repairs. The Hercules® XE Sterile Procedure Kit is a convenient, single-use package designed to promote efficiency and consistency in the OR.



Hercules® XE Sterile Procedure Kit

A40 SPB01

- A** 3.5mm Soft Tissue Protector
- B** 4.5mm Soft Tissue Protector
- C** Hercules®XE 3.5 x 15mm Anchor & Assembly
- D** Hercules®XE 4.5 x 18mm Anchor & Assembly
- E** 3.5mm Cannulated Drill Bit, 3.2mm
- F** 3.5mm Cannulated Drill Bit, 3.6mm
- G** 4.5mm Cannulated Drill Bit, 4.7mm
- H** 1.6mm Guide Wire, 8" (Qty. 2)
- J** Hercules®XE 3.5 x 15mm Anchor
- K** BioBrace® FA 5 x 250mm Implant

Sold Separately:

- L** BioBrace® FA 23 x 30mm Implant BB23X30-FA
- M** BioBrace® FA 40 x 60mm Implant BB40X60-FA



In2Bones USA
7750A Trinity Rd.
Cordova, TN 38018
844 602 6637
www.conmed.com/in2bones

In2Bones SAS
28, chemin du Petit Bois
69130 Ecully
France
+33 (0)4 72 29 26 26

CONMED Corporation
11311 Concept Blvd.
Largo, FL 33773
USA

Some products not available in all markets. This document provides information regarding medical devices and instruments in surgical procedures. It is not medical advice and healthcare professionals should use their own judgment before using to treat a particular patient. Healthcare professionals should be trained in the use of such devices before surgery and should always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any product.