



Designed for Use with  
the BioBrace® Implant

Radiolucent PEEK Material

Delivered Sterile and  
Pre-assembled

# Hercules® XE

Extended Eyelet Knotless Anchor  
for BioBrace®

**AUGMENTATION IMPLANT KIT NOW AVAILABLE**




Extended eyelet specifically  
designed for the 5mm x 250mm  
BioBrace® Reinforced Implant

PRODUCT OVERVIEW





# Hercules® XE

Extended Eyelet Knotless Anchor for BioBrace®




BioBrace® 5mm x 250mm Implant



Anterior Talofibular Ligament (ATFL)

The **Hercules® XE Knotless Anchor** features an extended eyelet intended for fixation of suture (soft tissue) to bone in the foot/ankle, and designed specifically for use with the 5mm x 250mm **BioBrace® Reinforced Implant** in medial and lateral ankle stabilizations and Achilles tendon repairs. The Hercules® XE Augmentation Implant Kit is a convenient, single-use package designed to promote efficiency and consistency in the OR.



## Hercules® XE Augmentation Implant Kit | A40 SPB01

- A** 3.5mm Soft Tissue Protector
- B** 4.5mm Soft Tissue Protector
- C** Hercules® XE 3.5mm x 15mm Anchor & Assembly (Loaded on Driver)
- D** Hercules® XE 4.5mm x 18mm Anchor & Assembly
- E** 3.5mm Cannulated Drill Bit, 3.2mm (Soft bone)
- F** 3.5mm Cannulated Drill Bit, 3.66mm (Standard bone)
- G** 4.5mm Cannulated Drill Bit, 4.67mm
- H** 1.6mm Guide Wire, 8" (Qty. 2)
- I** Hercules® XE 3.5mm x 15mm Anchor (Unloaded)
- J** BioBrace® FA 5mm x 250mm Implant (Tapered)



In2Bones USA  
7750A Trinity Rd.  
Cordova, TN 38018  
844 602 6637  
www.conmed.com/in2bones

In2Bones SAS  
28, chemin du Petit Bois  
69130 Ecully  
France  
+33 (0)4 72 29 26 26

CONMED Corporation  
11311 Concept Blvd.  
Largo, FL 33773  
USA

Some products not available in all markets. This document provides information regarding medical devices and instruments in surgical procedures. It is not medical advice and healthcare professionals should use their own judgment before using to treat a particular patient. Healthcare professionals should be trained in the use of such devices before surgery and should always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any product.