# 

ABC<sup>®</sup> Dissecting Electrodes<sup>™</sup>



HELIXAR<sup>™</sup> ELECTROSURGICAL GENERATOR WITH ABC<sup>®</sup>

ABC° HANDPIECES ARGON DISSECTING ELECTRODES BEND-A-BEAM° MALLEABLE ABC° HANDPIECES



To learn more about these and other innovative products, call **800-448-6506** or visit **CONMED.com** 



525 French Road, Utica, NY 13502

ABC® technology in contact mode allows surgeons to dissect while employing the benefits of Argon Beam Coagulation. Surgeons can also continue to use ABC® in non-contact mode simply by pulling the ABC® Dissecting Electrodes<sup>™</sup> away from tissue.

Simple to use, ABC° Dissecting Electrodes attach to a CONMED ABC° handpiece to convert from Argon Beam Coagulation to Argon Beam Dissect. By forcing argon gas over an electrode in contact or non-contact mode, ABC<sup>®</sup> Dissecting Electrodes<sup>™</sup> may provide the following benefits for surgical procedures:

### **Reduced Smoke from Surgical Field**

Argon gas used with ABC® Dissecting Electrodes carries surgical smoke away from the field-of-view.

## **Less Tissue Damage**

Argon gas flowing over the dissection may help to cool tissue at the perimeter of the intended cut. The latter is believed to lower the risk of post-operative infections.

#### **Clean Cut Hemostasis**

ABC<sup>®</sup> Dissecting Electrodes combine the benefits of argon gas with a tungsten electrode. This combination of electrosurgical modality provides a blend of cut and coagulation effects that may result in a cleaner cut with hemostasis.

**Ordering Information** 



ABC° Handpieces Compatible with Dissecting Electrodes

Quantity

Description	Per Box	Number
ABC® Dissecting Blade Electrode converts Argon Beam Coagulation to Argon Beam Dissect	20/Case	139330
ABC® Dissecting Blunt Needle Electrode converts Argon Beam Coagulation to Argon Beam Dissect	20/Case	139331
Laparoscopic 10 mm ABC® Probe For use in Laparoscopy	10/Case	134009
ABC® Single Use Nozzle (for use with ABC® Reusable Handpiece, #130500) For use in Open Procedures	20/Case	130343
ABC <sup>®</sup> Single Function Handpiece For use in Open Procedures	10/Case	130344
ABC® Angled 45° Handpiece For use in Open Procedures	10/Case	130345

Customer Service: 1-800-448-6506 International Sales: 1-315-797-8375 Fax: 1-800-438-3051

CONMED.com info@mail.CONMED.com

Catalog

#### ABC® DISSECTING **ELECTRODES**<sup>\*</sup>