n2bones

126 IS A GLOBAL EXTREMITY COMPANY

OSTEOCONDUCTIVE & OSTEOINDUCTIVE ALLOAID BP Demineralized Moldable Bone Putty

- Made from 100% Human Bone
- Shape and Conform to Bony Voids
- Retains Shape During Lavage
- Stores at Ambient Room Temperature
- Sterile Gamma Irradiated to SAL 10-6
- **AlloAid Crunch DBM** Putty **Sponge Cubes**



MADE FROM 100% HUMAN BONE

ALLOAID BP Demineralized Moldable Bone Putty

AlloAid Crunch, DBM Putty

- Available in flowable, formable putty and crunch in syringe
- Donor recovery and screening performed according to AATB and FDA guidelines
- Crunch consists of DBM putty & cortical cancellous bone chips

100% Cancellous Bone Sponge Cubes

- Moldable, sterile, demineralized/compressive bone matrix
- Cut and molded as needed, sponge-like fit and fill capabilities
- Natural porous scaffold for encourages osteoconductivity and promotes vascularization and stimulates new bone growth
- Add bone marrow aspirate or platelet rich plasma to create a composite graft

CAT NO.	STYLE	DELIVERY / VOLUME
CAA 001	. AlloAid DBM Crush Mix	Syringe 1.Occ
CAA 002	. AlloAid DBM Crush Mix	Syringe 2.5cc
CAA 005	. AlloAid DBM Crush Mix	Syringe 5cc
CAA 010	. AlloAid DBM Crush Mix	Syringe 10cc
PAA 001	. AlloAid DBM Putty	Syringe 1.Occ
PAA 025	. AlloAid DBM Putty	Syringe 2.5cc
PAA 050	. AlloAid DBM Putty	Syringe 5cc
PAA 010	. AlloAid DBM Putty	Syringe 10cc
AA014	. AlloAid Sponge Cube	14x14x14mm
SAA020	. AlloAid Sponge Cube	20x20x20mm

Representation of SAL (Sterility Assurance Level). Designates the occurrence a living microorganism surviving the sterilization process. SAL of 10⁻⁶ designates the possibility of finding an unsterile product will be 1 in a million. While other sterilization methods have also been shown to inactivate viruses, they can be detrimental to the biological and physical properties of allograft bone.



Sterile, Recyclable Instrument Set for Increased Efficiency; Saves Time, Saves Money*



- Sterile, Individually Packaged Implants
- Self-adhesive Label Traceability for Hospital & Patient Records
 - Sterile, Recyclable Instrument Set for Increased Efficiency; Saves Time, Saves Money*

*Projected savings of up to \$750 per case over standard hospital sterilization processes

844. 602. 6637 **OR EMAIL US NOW** i2b-USA.com Info@i2b-USA.com

All content contained herein is furnished for informational purposes only. i2b-USA does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgement and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and Precautions are listed in implant packaging insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure.
AlloAid name and logo are registered trademarks of i2b-USA or its subsidiaries. I.B.S., DuaFit, In2Bones and the In2Bones logo are registered trademarks of In2Bones or its subsidiaries. PEEK-OPTIMA® is a polymer produced by INVIBIO®.

i2b IS A GLOBAL EXTREMITY COMPANY

n

ones

CALL TO SCHEDULE A SURGERY TODAY

TRACEABLE

